



# Model 3550/3560 Series Integrating Nephelometer

Instruction Manual  
April 2001

<b>Product Overview</b>	<b>1</b>
<b>Unpacking and Setting Up the Hardware</b>	<b>2</b>
<b>Operations Overview</b>	<b>3</b>
<b>Using Nephelometer Software (DOS)</b>	<b>4</b>
<b>Using Nephelometer Software (Windows)</b>	<b>5</b>
<b>Using Serial Data Commands</b>	<b>6</b>
<b>Theory of Operation</b>	<b>7</b>
<b>Performing Maintenance</b>	<b>8</b>
<b>Contacting Customer Service</b>	<b>9</b>
<b>Appendixes and Index</b>	

# Manual History

The following is a history of the *Model 3550/3560 Series Integrating Nephelometer Instruction Manual*, part number 1933563.

<b>Revision</b>	<b>Date</b>
Preliminary	February 1994
Preliminary 1	August 1994
Preliminary 2	September 1994
Final	November 1994
A	May 1996
B	April 1997
B1	July 1998
C	December 2000
D	April 2001

This manual was first published, in preliminary form, February 1994.

This manual was published as a preliminary 1 in August 1994.

This manual was published as a preliminary 2 in September 1994, additions of text and figures were made throughout manual, Chapter 5 was added, and Chapter 8, "Troubleshooting Guide," was deleted.

The manual was published as a Final after Chapter 5, "Using Nephelometer Software (Windows)" was added in November 1994.

In revision A, TSI's "Limitation of Warranty and Liability" on page iii was updated.

In revision B, Tables A-1 and A-2 in Appendix A were revised.

In revision B1, TSI's area code changed from 612 to 651.

In revision C, TSI's Limitation of Warranty and Liability and TSI's Software License were updated.

In revision D, Analog Outputs were removed.

**Part Number**

1933563 / Revision D / April 2001

**Copyright**

©TSI Incorporated / 1994–2001 / All rights reserved.

**Address**

TSI Incorporated / 500 Cardigan Road / P.O. Box 64394 / St. Paul, MN 55164 / USA

**Fax No.**

(651) 490-3860

**Email Address**

particle@tsi.com

**Limitation of Warranty and Liability**

(effective July 2000)

Seller warrants the goods sold hereunder, under normal use and service as described in the operator's manual, shall be free from defects in workmanship and material for (12) months, or the length of time specified in the operator's manual, from the date of shipment to the customer. This warranty period is inclusive of any statutory warranty. This limited warranty is subject to the following exclusions:

- a. Hot-wire or hot-film sensors used with research anemometers, and certain other components when indicated in specifications, are warranted for 90 days from the date of shipment.
- b. Parts repaired or replaced as a result of repair services are warranted to be free from defects in workmanship and material, under normal use, for 90 days from the date of shipment.
- c. Seller does not provide any warranty on finished goods manufactured by others or on any fuses, batteries or other consumable materials. Only the original manufacturer's warranty applies.
- d. Unless specifically authorized in a separate writing by Seller, Seller makes no warranty with respect to, and shall have no liability in connection with, goods which are incorporated into other products or equipment, or which are modified by any person other than Seller.

The foregoing is IN LIEU OF all other warranties and is subject to the LIMITATIONS stated herein. **NO OTHER EXPRESS OR IMPLIED WARRANTY OF FITNESS FOR PARTICULAR PURPOSE OR MERCHANTABILITY IS MADE.**

TO THE EXTENT PERMITTED BY LAW, THE EXCLUSIVE REMEDY OF THE USER OR BUYER, AND THE LIMIT OF SELLER'S LIABILITY FOR ANY AND ALL LOSSES, INJURIES, OR DAMAGES CONCERNING THE GOODS (INCLUDING CLAIMS BASED ON CONTRACT, NEGLIGENCE, TORT, STRICT LIABILITY OR OTHERWISE) SHALL BE THE RETURN OF GOODS TO SELLER AND THE REFUND OF THE PURCHASE PRICE, OR, AT THE OPTION OF SELLER, THE REPAIR OR REPLACEMENT OF THE GOODS. IN NO EVENT SHALL SELLER BE LIABLE FOR ANY SPECIAL, CONSEQUENTIAL OR INCIDENTAL DAMAGES. SELLER SHALL NOT BE RESPONSIBLE FOR INSTALLATION, DISMANTLING OR REINSTALLATION COSTS OR CHARGES. No Action, regardless of form, may be brought against Seller more than 12 months after a cause of action has accrued. The goods returned under warranty to Seller's factory shall be at Buyer's risk of loss, and will be returned, if at all, at Seller's risk of loss.

Buyer and all users are deemed to have accepted this LIMITATION OF WARRANTY AND LIABILITY, which contains the complete and exclusive limited warranty of Seller. This LIMITATION OF WARRANTY AND LIABILITY may not be amended, modified or its terms waived, except by writing signed by an Officer of Seller.

**Software License**  
(effective March 1999)

1. GRANT OF LICENSE. TSI grants to you the right to use one copy of the enclosed TSI software program (the "SOFTWARE"), on a single computer. You may not network the SOFTWARE or otherwise use it on more than one computer or computer terminal at the same time.
2. COPYRIGHT. The SOFTWARE is owned by TSI and is protected by United States copyright laws and international treaty provisions. Therefore, you must treat the SOFTWARE like any other copyrighted material (e.g., a book or musical recording) except that you may either (a) make one copy of the SOFTWARE solely for backup or archival purposes, or (b) transfer the SOFTWARE to a single hard disk provided you keep the original solely for backup or archival purposes.
3. OTHER RESTRICTIONS. You may not rent or lease the SOFTWARE, but you may transfer the SOFTWARE and accompanying written material on a permanent basis, provided you retain no copies and the recipient agrees to the terms of this Agreement. You may not reverse-engineer, decompile, or disassemble the SOFTWARE.
4. DUAL MEDIA SOFTWARE. If the SOFTWARE package contains multiple types of media, then you may use only the media appropriate for your single-user computer. You may not use the other media on another computer or loan, rent, lease, or transfer them to another user except as part of the permanent transfer (as provided above) of all SOFTWARE and written material.
5. U.S. GOVERNMENT RESTRICTED RIGHTS. The SOFTWARE and documentation are provided with RESTRICTED RIGHTS. Use, duplication, or disclosure by the Government is subject to the restrictions set forth in the "Rights in Technical Data and Computer Software" Clause at 252.227-7013 and the "Commercial Computer Software - Restricted Rights" clause at 52.227-19.
6. LIMITED WARRANTY. TSI warrants that the SOFTWARE will perform substantially in accordance with the accompanying written materials for a period of ninety (90) days from the date of receipt.
7. CUSTOMER REMEDIES. TSI's entire liability and your exclusive remedy shall be, at TSI's option, either (a) return of the price paid or (b) repair or replacement of the SOFTWARE that does not meet this Limited Warranty and which is returned to TSI with proof of payment. This Limited Warranty is void if failure of the SOFTWARE has resulted from accident, abuse, or misapplication. Any replacement SOFTWARE will be warranted for the remainder of the original warranty period or thirty (30) days, whichever is longer.
8. NO OTHER WARRANTIES. TSI disclaims all other warranties, either express or implied, including, but not limited to implied warranties of merchantability and fitness for a particular purpose, with regard to the SOFTWARE and the accompanying written materials.
9. NO LIABILITY FOR CONSEQUENTIAL DAMAGES. In no event shall TSI be liable for any damages whatsoever (including, without limitation, special, incidental, consequential or indirect damages for personal injury, loss of business profits, business interruption, loss of information or any other pecuniary loss) arising out of the use of, or inability to use, this SOFTWARE.

**Service Policy**

Knowing that inoperative or defective instruments are as detrimental to TSI as they are to our customers, our service policy is designed to give prompt attention to any problems. If any malfunction is discovered, please contact your nearest sales office or representative, or call TSI's Particle Instrument Division at 1-800-874-3893 (USA) or (651) 490-3893.

# Safety


This section gives instructions to promote safe and proper handling of the Model 3550/3560 Series Integrating Nephelometers.



## **C a u t i o n**

High temperatures that can cause burns.

To avoid personal injury, disconnect power to the Nephelometer and allow the halogen lamp and the lamp housing to cool before handling.



## **W A R N I N G**

High voltages that can shock or burn.

The box that houses the photomultiplier tubes (PMT) contains a high voltage source. Components within this box should only be serviced by a qualified technician. Disconnect all power to the Nephelometer and contact TSI personnel before disassembling the PMT box.



## **C a u t i o n**

Light sensitive components.

To avoid damage to the photomultiplier tubes (PMT), remove power from the Nephelometer and open the PMT box in subdued lighting. Exposing the photomultiplier tubes to direct light may affect tube performance.




## **W A R N I N G**

Moving parts that can cause serious personal injury.

Disconnect power to Nephelometer before inserting anything into the aerosol inlet. The ball valve and position sensor are motorized and can move if power is applied.

# Labels

Figure 1 shows the special warning and caution labels and their approximate location inside the Nephelometer.



**Figure 1**  
Nephelometer Warning and Caution Labels

# Contents

<b>Manual History</b> .....	ii
<b>Warranty</b> .....	iii
<b>Software License</b> .....	iv
<b>Safety</b> .....	v
Labels .....	vi
<b>About This Manual</b> .....	xvii
Purpose.....	xvii
Reusing and Recycling .....	xvii
Submitting Comments .....	xvii
Notational Conventions .....	xviii

---

## Chapters/Appendixes

<b>1 Product Overview</b> .....	1-1
Product Description .....	1-1
Applications .....	1-1
Beer-Lambert Law .....	1-2
How the Nephelometer Operates .....	1-2
<b>2 Unpacking and Setting Up the Hardware</b> .....	2-1
Packing List .....	2-1
Unpacking Instructions.....	2-2
Connecting Power .....	2-2
Connecting Line Voltage to the Power Supply Module.....	2-2
Connecting the Power Supply Module to the Nephelometer .....	2-3
Using Another Power Supply .....	2-4
Connecting a Computer .....	2-4
Mounting the Nephelometer .....	2-5
Space Requirements.....	2-7
Mounting Holes.....	2-8
Inlet and Outlet Ports.....	2-8
Environmental Concerns .....	2-8
Connecting an External Blower .....	2-9
Calibrating the Nephelometer.....	2-11
<b>3 Operations Overview</b> .....	3-1
Serial Communications.....	3-1
Hardware Components.....	3-1
Power Failures .....	3-1
Lamp Power .....	3-2

Photomultiplier Tubes (PMTs).....	3-2
Heater .....	3-2
Blower Settings .....	3-3
Humidity and Temperature Sensors .....	3-3
Setting Up Parameters .....	3-4
Operating Parameters.....	3-5
Sensitivity .....	3-5
Backscatter Measurements .....	3-5
Zero Baseline Measurements.....	3-6
Autozero Modes.....	3-7
Time Parameters .....	3-7
Averaging Time.....	3-8
Zero Time .....	3-8
Auto Zero Period.....	3-8
Valve Blanking Time.....	3-8
Polled and Unpolled Communications .....	3-9
<b>4 Using Nephelometer Software (DOS).....</b>	<b>4-1</b>
Hardware and Software Requirements .....	4-1
Installing and Starting the Software .....	4-1
Keys and Conventions.....	4-3
Menus and Command Options: A Reference Guide .....	4-4
Main Menu.....	4-5
Main/Set Commands Menu .....	4-6
Set Command/K Calibration Constants Option.....	4-6
Set Command/Calibration Label String Option .....	4-7
Set Command/Mode for BackScatter Shutter Option ....	4-8
Set Command/Lamp Power Option .....	4-9
Set Command/Time for Averaging Option.....	4-9
Set Command/Time and Date Option.....	4-10
Set Command/Voltage on Photomultiplier Tube Option .....	4-11
Main/Calibration Option .....	4-12
Main/Data Collection (Polled) Option .....	4-12
Main/Configuration Menu.....	4-14
Configuration/Display Current Config. Data Option ....	4-14
Configuration/Read Config. Data from Nephelometer Option .....	4-16
Configuration/Read Config. Data from Disk File Option .....	4-17
Configuration/Save Config. Data to Nephelometer Option .....	4-17
Configuration/Save Config. Data to Disk File Option...	4-17
Configuration/Print Config. Data Option .....	4-17



Main/Terminal Mode Option .....	4-17
Main/Com Port and Data Path Option .....	4-18
Performing Calibration .....	4-19
Hardware Setup .....	4-19
Using the Main/Calibration Option .....	4-19
<b>5 Using Nephelometer Software (Windows) .....</b>	<b>5-1</b>
Hardware and Software Requirements.....	5-1
Installing and Starting the Software .....	5-2
Quick Start .....	5-4
Step 1. Verify Com Port .....	5-5
Step 2. Start Data Collection .....	5-6
Step 3. View Data and Nephelometer Status.....	5-7
Step 4. Print a Graph .....	5-10
Step 5. Perform a Background (Zero Baseline)	
Measurement .....	5-10
Basic Nephelometer Operations.....	5-12
Clean Air Test .....	5-14
Menus and Commands: A Reference Guide .....	5-15
Using the File Menu .....	5-16
Open Log File Command .....	5-16
Log Data.....	5-16
Print.....	5-18
Print Screen .....	5-19
Exit .....	5-20
Using the Config Menu .....	5-20
Nephelometer .....	5-20
Com Port .....	5-24
Using the View Menu.....	5-24
Data .....	5-24
Graph.....	5-26
Photon Frequency .....	5-29
Neph State .....	5-31
Neph Config Data .....	5-31
Log Data Table .....	5-33
Using the Run Menu .....	5-39
Calibration .....	5-39
Data Collection.....	5-47
Terminal Mode .....	5-49
Using the Window Menu.....	5-49
Using the Pause/Resume Menu .....	5-50
Using the Help Menu.....	5-50
Contents .....	5-50
About .....	5-51

<b>6</b>	<b>Using Serial Data Commands</b> .....	6-1
	Pin Connectors.....	6-1
	Baud Rate.....	6-2
	Parity (8-Bits Even) .....	6-2
	Commands.....	6-3
	Set Commands.....	6-5
	Action Commands.....	6-19
	Read Commands (Polled).....	6-22
	Unpolled Commands .....	6-27
	Unpolled Record Formats .....	6-28
	Time Record (T) .....	6-28
	Photon Count Records (B, G, and R).....	6-29
	Data Record (D).....	6-30
	Auxiliary Status Record (Y).....	6-30
	Zero Background Data Record (Z).....	6-31
	T, B, G, R, D, Y and Z Records.....	6-31
	How to Input Commands and Troubleshoot the Results .	6-32
	Input Guidelines .....	6-32
	Troubleshooting Input .....	6-33
<b>7</b>	<b>Theory of Operation</b> .....	7-1
	History.....	7-1
	Construction.....	7-2
	Theory of Integration.....	7-4
	Theory of Wavelength Dependence .....	7-7
	Signal Processing .....	7-10
	Description.....	7-10
	Equations/Calculations .....	7-10
	Firmware Calculations .....	7-11
	Calibration Calculations.....	7-15
	Detection Optics.....	7-18
	Calibration.....	7-18
<b>8</b>	<b>Performing Maintenance</b> .....	8-1
	Removing Nephelometer Covers: Top, Bottom, and PMT.....	8-2
	Removing the Top Cover.....	8-2
	Removing the Bottom Cover .....	8-3
	Removing the PMT Cover.....	8-5
	Calibrating the Nephelometer.....	8-6
	Cleaning Fan Filters.....	8-7
	Replacing the Main Microprocessor EPROM.....	8-8
	Replacing the Motor Control Microprocessor .....	8-11
	Replacing the Lamp.....	8-13
	Replacing Aerosol Filters .....	8-15
	Cleaning the Light Pipe Lens .....	8-19

Checking for Leaks.....	8-22
Cleaning the Reference Chopper.....	8-23
Cleaning the Light Pipe Outlet and the Backscatter Shutter.....	8-31
Cleaning the Flocked Paper.....	8-41
Replacing a Temperature or Humidity Sensor.....	8-47
Replacing a Photomultiplier Tube.....	8-49
Checking, Cleaning, and Replacing Bandpass Filters.....	8-52
Cleaning the Lens.....	8-56
<b>9 Contacting Customer Service.....</b>	<b>9-1</b>
Technical Contacts at TSI.....	9-1
Returning the Nephelometer for Service.....	9-1
<b>A Specifications.....</b>	<b>A-1</b>
<b>B Bibliography.....</b>	<b>B-1</b>
Instrument Principles, Design and Calibration.....	B-2
Applications of the Integrating Nephelometer to Studies of Atmospheric Visibility.....	B-4
Light Scattering Measurements Related to Aerosol Mass Concentration.....	B-4
Vertical Column and Horizontal Profile Studies with the Integrating Nephelometer.....	B-5
Apportioning of Light Scattering Coefficient to Chemical Species and to Sources.....	B-6
Use of the Integrating Nephelometer in Studies of Physical and Chemical Properties of Aerosols.....	B-7
Observation of the Effects of Humidity on Light Scattering Coefficient; Measurements on Fogs and Clouds.....	B-8
Observation of Aerosols in Remote Locations with the Integrating Nephelometer.....	B-8
The Use of Integrating Nephelometer Data in Radiation and Climate Studies.....	B-9
<b>C Schematics.....</b>	<b>C-1</b>
<b>D Windows Program Files.....</b>	<b>D-1</b>

---

## Figures

1 Nephelometer Warning and Caution Labels.....	vi
1-1 TSI Incorporated's Nephelometer.....	1-1
1-2 A Cross-Section View of the Nephelometer.....	1-3

2-1	Power Supply Module .....	2-3
2-2	Ventilation Requirements .....	2-4
2-3	COM PORT Connector .....	2-5
2-4	Nephelometer Resting on its Feet.....	2-6
2-5	Nephelometer in the Optimal Orientation .....	2-6
2-6	As-Built Dimensions (Space Requirements) .....	2-7
2-7	Blower Bracket .....	2-10
2-8	Mounting the External Blower .....	2-10
3-1	Humidity and Temperature Sensor Locations .....	3-4
4-1	Main Menu Screen.....	4-2
4-2	Software Command Conventions .....	4-3
4-3	Software Menu Structure.....	4-4
4-4	Main Menu Screen.....	4-5
4-5	Set Commands Menu .....	4-6
4-6	K Calibration Constants Option.....	4-7
4-7	Calibration Label Option.....	4-8
4-8	BackScatter Shutter Mode Option .....	4-8
4-9	Lamp Power Option .....	4-9
4-10	Time for Averaging Option .....	4-10
4-11	Time and Date Option.....	4-10
4-12	Voltage on Photomultiplier Tube Option .....	4-11
4-13	Data Collection Polled Option .....	4-13
4-14	Configuration Menu.....	4-14
4-15	Display Current Config. Data Option .....	4-15
4-16	Display Current Config. Data Option .....	4-16
4-17	Com Port and Data Path Option .....	4-18
4-18	Calibration Option.....	4-20
4-19	Plugs for Inlet and Outlet and Gas Line With the Blue DQ Filter .....	4-22
5-1	Select the Destination Directory for Program Installation..	5-3
5-2	Setup Creates a Program Group and Icon.....	5-4
5-3	The Main Menu of the Windows-Based Nephelometer Program .....	5-5
5-4	The COM PORT Setup Display.....	5-5
5-5	Choosing Data Collection from the Run Menu Starts Data Collection.....	5-7
5-6	Select Data Display to View Data as it is Collected by the Nephelometer .....	5-8
5-7	The Nephelometer State Display Shows the Status of the Nephelometer .....	5-8
5-8	The Photon Frequency Display Shows the Current Photon Frequency Calculations .....	5-9

5-9	The Graph Display Shows Current Data .....	5-9
5-10	The Terminal Mode Display .....	5-11
5-11	The Set Nephelometer Parameters Display Shows the Current Parameters of the Nephelometer .....	5-14
5-12	Select a Filename After Choosing the Open Log File Command.....	5-17
5-13	Select a Filename After Choosing the Log Data Command.....	5-17
5-14	Select the Print Command from the File Menu.....	5-18
5-15	Types of Graphs. Scatter Data Graph (top) and Status Graph (bottom).....	5-19
5-16	Choose Nephelometer from the Config Menu.....	5-21
5-17	Use the COM Port Setup Display to Select the Appropriate Communications Port.....	5-24
5-18	The Data Display Shows Current (raw) Data Values Collected by the Nephelometer.....	5-25
5-19	Graph Data of the Current Log Data File .....	5-27
5-20	The Graph Options Window Lets You Customize Your Graphs .....	5-28
5-21	Photon Frequency Display .....	5-29
5-22	The Photon Frequency Graph Displays Current Measurements.....	5-30
5-23	The Status of the Nephelometer is Displayed by Selecting Neph State from the View Menu.....	5-31
5-24	Nephelometer Configuration Data is Displayed by Selecting Neph Config Data from the View Menu .....	5-32
5-25	Display Log File Records in Tabular Format by Selecting Log Data Table from the View Menu .....	5-34
5-26	Select the Type of Graph from the Log Data Graph Options .....	5-36
5-27a	Operating Status Graph.....	5-36
5-27b	Photon Counts (Calibrator) Graph .....	5-36
5-27c	Photon Counts (Dark) Graph.....	5-37
5-27d	Photon Counts (Measure) Graph.....	5-37
5-27e	Scattering Coefficients Graph.....	5-38
5-28	The Destination Dialog Box .....	5-38
5-29	The Nephelometer Calibration Data Display.....	5-39
5-30	Calibration Parameters Setup Screen .....	5-41
5-31	Selecting Data Collection from the Run Menu.....	5-48
5-32	Select Terminal Mode to Issue Serial Data Commands to the Nephelometer .....	5-49
5-33	The Help Contents Display Provides a List of Help Topics .....	5-51

5-34	About Shows the Current Version of Your Software and the Copyright Notice .....	5-51
6-1	COM PORT Pin Designations .....	6-1
6-2	Timelines for Data Measurement .....	6-12
6-3	Comparison of Manual, Normal, and Air Chop Modes.....	6-13
7-1	Nephelometer Schematic .....	7-2
7-2	Reference Chopper Zones .....	7-3
7-3	Nephelometer Geometry.....	7-6
8-1	Nephelometer Showing Top, Bottom, and PMT Covers .....	8-2
8-2	Top Cover Screws, Bottom Cover Screws, Fan Cover Screws, and I/O Panel Screws.....	8-3
8-3	PMT Cover.....	8-5
8-4	Fan Covers and Foam Inserts .....	8-7
8-5	Top Cover Screws .....	8-9
8-6	Microprocessor and EPROM Locations on the Digital Circuit Board.....	8-12
8-7	Top Cover Screws .....	8-13
8-8	Lamp Shield .....	8-14
8-9	Top Cover Screws .....	8-16
8-10	Bottom Cover Screws Securing End Plate .....	8-16
8-11	White HEPA Filter.....	8-17
8-12	Blue DQ Filter .....	8-18
8-13	Top Cover Screws .....	8-19
8-14	Lamp Shield .....	8-20
8-15	Light Pipe Lens .....	8-21
8-16	Top Cover Screws, Bottom Cover Screws, Fan Cover Screws, and I/O Panel Screws.....	8-24
8-17	Circuit Board Mounting Screws .....	8-25
8-17a	PMT End Plate Screws .....	8-26
8-18	Aperture Assembly .....	8-27
8-19	Aperture Plates and Circuit Board .....	8-28
8-20	Reference Chopper.....	8-28
8-21	Top Cover Screws, Bottom Cover Screws, Fan Cover Screws, and I/O Panel Screws.....	8-32
8-22	Lamp Shield .....	8-32
8-23	Disconnect Connectors and Tubing .....	8-34
8-24	Nephelometer Sections Together .....	8-35
8-25	Nephelometer Sections Apart.....	8-35
8-26	Lamp Base .....	8-36
8-27	Removing the Lamp Base Assembly From the Sensor .....	8-37
8-28	Light Pipe Outlet and Backscatter Shutter.....	8-38
8-29	O-ring on the Sensor Backbone .....	8-39

8-30 Top Cover Screws, Bottom Cover Screws, Fan Cover Screws, and I/O Panel Screws .....	8-41
8-31 Disconnect Connectors and Tubing .....	8-43
8-32 Nephelometer Sections Together .....	8-44
8-33 Nephelometer Sections Apart .....	8-44
8-34 Flocked Paper Outside the Nephelometer .....	8-45
8-35 Top Cover Screws, Bottom Cover Screws, Fan Cover Screws, and I/O Panel Screws .....	8-47
8-36 Humidity and Temperature Sensor Locations .....	8-48
8-37 PMT Cover .....	8-50
8-38 Red, Blue and Green PMTs with Setscrews .....	8-51
8-39 PMT Cover .....	8-53
8-40 Red, Blue and Green PMTs with Setscrews .....	8-54
8-41 Bandpass Filter .....	8-55
8-42 PMT Cover .....	8-57
8-43 Red, Blue and Green PMTs with Setscrews .....	8-58

---

## Tables

2-1 Packing List with Accessories .....	2-1
3-1 Typical Time Settings .....	3-10
4-1 Keyboard Functions .....	4-3
4-2 Calibration Procedure .....	4-20
5-1 Basic Windows Nephelometer Operations .....	5-12
5-2 Nephelometer Setup Parameters .....	5-21
5-3 Data Display Parameters .....	5-25
5-4 Graph Options for the Graph Display .....	5-28
5-5 Nephelometer Configuration Data .....	5-33
5-6 Calibration Parameters .....	5-42
5-7 Calibration Procedure .....	5-43
6-1 Signal Connections for RS-232 Configurations .....	6-2
6-2 Serial Commands .....	6-4
7-1 Examples of Raw Count Data .....	7-11
7-2 Comparisons Between Gate Width and Physical Dimensions .....	7-12
7-3 Counts Normalize to Photon Frequency (Hz) .....	7-12
7-4 Boxcar Averages .....	7-13
8-1 Maintenance Overview .....	8-1

**A-1 Specifications for the Nephelometer Sensor ..... A-1**



# About This Manual

---

## Purpose

This is an instruction manual for the operation and maintenance of the Model 3550/3560 Series Integrating Nephelometers.

---

## Reusing and Recycling



As part of TSI Incorporated's effort to have a minimal negative impact on the communities in which its products are manufactured and used:

- This manual uses recyclable paper.
- This manual has been shipped, along with the instrument, in a reusable carton.

---

## Submitting Comments

TSI values your comments and suggestions on this manual. Please use the comment sheet, located on the last page of this manual, to send us your opinion on the manual's usability, to suggest specific improvements, or to report any technical errors.

If the comment sheet has already been used, mail or fax your comments on another sheet of paper to:

Particle Instrument Division  
TSI Incorporated  
P.O. Box 64394  
St. Paul, MN 55164  
Fax: (651) 490-3860

---

## Notational Conventions

This manual uses the following conventions when describing software:

- <Enter>** Denotes the “Return” or “Enter” key on the keyboard.
- <xx>** Denotes an alphanumeric key on the keyboard. For instance, **<F4>** means press the function “F4” on the keyboard.
- examples* Examples of what you see on the screen and the text you type are shown in monospace type resembling computer output.
- italic* Variable information in computer responses, commands, and options you must supply and type are shown in italics in a font that resembles computer font.
- UPPERCASE Uppercase letters indicate names of programs, files or commands.
- Bold** Bold type indicates names of fittings as they appear on the instrument and commands appearing on the display.

## CHAPTER 1

# Product Overview

This chapter contains a product description and a list of features for the Model 3550/3560 Series Integrating Nephelometer, as well as a brief description of how the instrument works.

---

## Product Description

The Model 3550/3560 Series Integrating Nephelometer, shown in Figure 1-1, is designed for long-term monitoring of visual range and air quality in ground-based and airborne studies. It continuously monitors the light scattering coefficient of airborne particles.



**Figure 1-1**  
TSI Incorporated's Nephelometer

---

## Applications

TSI Integrating Nephelometers are designed specifically for studies of direct radiative forcing of the Earth's climate by aerosol particles, or studies of ground-based or airborne atmospheric visual air quality in clean areas. They may also be used as an analytical detector for aerosol particles whenever the parameter of interest is the light-scattering coefficient of the particles after a pretreatment step, such as heating, humidification, or segregation by size.

The light-scattering coefficient is a highly variable aerosol property. Integrating Nephelometers measure the angular integral of light scattering that yields the quantity called the *scattering coefficient*, used in the Beer-Lambert Law to calculate total light extinction.

### **Beer-Lambert Law**

$$I/I_o = e^{(-\sigma x)}$$

where:

$I_o$  = intensity of light source

$I$  = intensity of light after passing through atmospheric path

$x$  = thickness of medium through which light passes


$\sigma$  = total extinction coefficient

(= scattering coefficient + absorption coefficient)

---

## **How the Nephelometer Operates**

A small, turbine blower draws an aerosol sample through the large diameter inlet into the measurement volume (Figure 1-2). There, the sample is illuminated over an angle of 7 to 170° by a halogen lamp directed through an optical light pipe and opal glass diffuser.



**Figure 1-2**  
A Cross-Section View of the Nephelometer

The sample volume is viewed by three photomultiplier tubes (PMT) through a series of apertures set along the axis of the main instrument body. Aerosol scattering is viewed against the backdrop of a very efficient light trap. The light trap, apertures, and a highly light-absorbing coating on all internal surfaces of the instrument combine to give a very low scatter signal from the walls of the instrument.

Dichroic filters, in front of the PMT tubes, split and direct the light, which has been scattered by aerosol. The light is directed into three bandpass filters, blue, green and red. A constantly rotating reference chopper has separate areas to provide three types of signal detection. The first area gives a measure of the aerosol light-scattering signal allowed by an opening in the rotating chopper. The second area blocks all light from detection and gives a measurement of the PMT dark current that which is subtracted from the measurement signal. The third area is a translucent portion of the chopper, illuminated by the halogen lamp, which provides a measure of the light-source signal. In this way, over time, any change in the light source or in detector efficiency is compensated.

In backscatter mode, the backscatter shutter rotates under the lamp to block light in the 7 to 90° range. When light is blocked, only light scattered in the backward direction is transmitted to the PMT detectors. The backscatter signal can be subtracted from the total signal to calculate forward-scattering signal data. When this measurement is not of interest, the backscatter shutter can be “parked” in the total scatter position.

Periodically, an automated valve built into the inlet can be activated to divert all of the aerosol sample through a high-efficiency (HEPA) filter. This gives a measure of the clean-air signal for the local environment. This signal is subtracted, along with the PMT dark-current signal, from the aerosol-scatter signal to give only that portion of the scatter signal provided by the sample aerosol. Particle-scattering parameters for all three wavelengths of total and backscatter signal are continuously averaged and passed to a computer or data logger for permanent storage.

## CHAPTER 2

# Unpacking and Setting Up the Hardware

Use the information in this chapter to unpack and set up the hardware components of a Model 3550/3560 Series Integrating Nephelometer.

This chapter contains these main sections:

- Packing List
- Unpacking instructions
- Connecting power
- Connecting a computer
- Mounting the Nephelometer
- Connecting an external blower

---

## Packing List

Table 2-1 gives a packing list for the Nephelometer and the power supply module. The packing list for the Data Analysis Center is included in one of the two or three Data Analysis Center shipping cartons.

**Table 2-1**  
Packing List with Accessories

<b>Qty</b>	<b>Description</b>	<b>Part Number</b>
1	Nephelometer	3551, 3553, 3561 or 3563
1	Power supply module	3590
1	Power supply cable	1035564
1	Power cord (24 volts dc)	1035551
1	Blower bypass	1035545
1	Serial data cable (9-pin, 12 ft.)	962002
1	Serial cable adapter (25F-9M)	962003
1	Program disk (3.5 in. format)	1906118
1	Instruction manual	1933563
1	White HEPA filter	1602051
1	Blue DQ filter	1602080
1	Set (2) fan filters (foam inserts)	1602071
1	Lamp (halogen)	2201111

---

## Unpacking Instructions

The power supply module and power cord are in one box; the Nephelometer sensor, power cable, manual, software disk, and other accessories are in another box. Keep the packing material in case you have to send the Nephelometer back to TSI.

The Nephelometer sensor comes fully assembled. If anything is missing or appears to be damaged, contact your TSI representative, or contact Customer Service at 1-800-874-3893 (USA) or (651) 490-3893. Chapter 9, "Contacting Customer Service," gives instructions for returning the Nephelometer to TSI Incorporated.

---


## Connecting Power

This section describes connecting power to the power supply module and connecting the power supply module to the Nephelometer sensor, as well as using another power source for the Nephelometer sensor.

### Connecting Line Voltage to the Power Supply Module

Use the line cord (supplied) to connect the TSI power supply module to any line voltage from 85–260 volts AC at 50–60 Hz (Figure 2-1). The auto-switching power supply automatically adjusts to the AC voltage provided.






**Figure 2-1**  
Power Supply Module

## Connecting the Power Supply Module to the Nephelometer

The power supply module includes a 4-meter (12-ft.) cord with four-conductor, quarter-turn quick connectors. Before connecting the cord, make sure the power switch is in the Off position. Connect the pin end to the power supply and the socket end to the POWER AC INPUT connector (Figure 2-1) on the Nephelometer sensor.

The TSI power supply module can be oriented in any direction, but the cooling fan intake and exhaust vents should be free from obstructions at all times (Figure 2-2).

**Note:** *The power supply module contains no user-serviceable parts. If the module is not operating properly, use the information in Chapter 9, “Contacting Customer Service,” to contact TSI.*



**Figure 2-2**  
Ventilation Requirements

## Using Another Power Supply


As an option, you can connect the Nephelometer sensor to a 24–28 volts DC source. The lamp power is maximum at the higher end of the DC range. You may order an extra cable (1035564) and wire one end by connecting the red wire to +, the white wire to –, and the green wire to chassis ground.

---

## Connecting a Computer

Connect the serial port of an IBM-compatible computer to the COM PORT connector on the Nephelometer sensor (Figure 2-3). Use the 4-meter cable provided, and if you need additional length, use a standard IBM 9-pin, serial extension cable.

**Note:** Refer to Chapter 4 for instructions on using TSI software and refer to Chapter 6 for information on using serial data commands.




**Figure 2-3**  
COM PORT Connector

---


## Mounting the Nephelometer

You can mount the Nephelometer in a variety of positions depending on the application. Normally, the Nephelometer rests on its feet (Figure 2-4), but the best position is a vertical mount with the PMT box at the top and the light trap end at the bottom (Figure 2-5).

**Note:** Other positions are possible but consider dust collection opposite the lamp or on the lens. Dust raises the background noise level and necessitates cleaning more often.




**Figure 2-4**  
Nephelometer Resting on its Feet



**Figure 2-5**  
Nephelometer in the Optimal Orientation

## Space Requirements

The Nephelometer sensor requires a rectangular space (Figure 2-6). You also need to allow room for the power supply module, air ventilation, and a cable to reach the sensor. If you are using a computer to set up the sensor, allow space for that as well.



**Figure 2-6**  
As-Built Dimensions (Space Requirements)

## Mounting Holes

Use the four mounting holes, two on each leg, to mount the Nephelometer sensor (Figure 2-6c). Use four 6 mm (¼ in.) diameter bolts or screws to mount the instrument to a secure bracket, if needed.

## Inlet and Outlet Ports

The inlet and outlet ports are 1 in. diameter female NPT (pipe threaded) fittings (Figure 2-6a). These fittings can be used with adapters to mate to smaller or larger rigid tubing or can be fastened directly to 1 in. pipe with a 1 in. NPT end. Rigid PVC tubing is often used to bring a sample into the instrument and to exhaust the sample from the instrument.

Keep in mind that the inlet and outlet tubing should have at least 8 inches of straight length to allow the cover to be easily removed for maintenance.

## Environmental Concerns

The Nephelometer is a highly sensitive instrument that detects light-scattering from air (gas) molecules and particles present in the sample chamber. The Nephelometer also detects undesirable scatters from insects, birds, and rodents that may come into the sampling chamber. Pests can severely limit the Nephelometer's ability to measure scattering from particles.

One way to minimize pests is to use insect screens, fly, or insect traps in the upstream sampling line to prevent insects or small pests from coming into the sampling chamber. Use metal screens to prevent larger pests from coming into the sampling chamber.

It is good to keep in mind that humidities above 50–70% will enhance scattering extinction (normally assumed to be small) by particles. You may want to consider a heater or desiccant on the inlet flow to keep the humidity below this level. Although the Nephelometer incorporates a heater in the body, this is intended to compensate for any cooling effects of the instrument body. It is *not* intended to reduce the humidity through the instrument.

If the instrument will be used to sample from a warm, moist environment, while located in a cooled enclosure, you should consider insulating the inlet tubing and using the Nephelometer heater feature (see H command in Chapter 6).

In designing the inlet and plumbing systems, be aware that changes in building pressure or a venturi effect caused by high winds, can cause flow changes and reversals through the Nephelometer.

**Note:** *Minimize flow restrictions and particle loss when taking these measures.*

---


## Connecting an External Blower

You can use the Nephelometer with an external blower to move the sample through the instrument. If you use an external blower, you can remove the onboard blower and replace it with the blower bypass fitting to reduce pressure drop through the sampling volume.

**Note:** *You need custom software to control the external blower. As of this date, TSI software does not include an interface to control an external blower.*


To remove the onboard blower and install the external blower:

1. Remove the top cover of the Nephelometer using the instructions in Chapter 8, "Performing Maintenance."
2. Remove the two screws attaching the blower bracket to the Nephelometer (Figure 2-7).
3. O-ring seals on the ports hold the blower in place. *Firmly* pull the blower away from the mounting plate.
4. Insert the blower bypass fitting into the pair of holes previously occupied by the blower.



**Figure 2-7**  
Blower Bracket

5. Install the external blower bypass fitting by pressing it into the hole in the mounting plate (Figure 2-8) and reattaching the blower bracket.



**Figure 2-8**  
Mounting the External Blower



---

## Calibrating the Nephelometer

After you set up the Nephelometer, take a sample reading and make a printout of the reading. To check the Nephelometer calibration, use the sample reading to compare with a reading taken during calibration. See the "Performing Calibration" section in Chapter 4 for the calibration procedure.



## CHAPTER 3

# Operations Overview

Before you set up the Nephelometer, it is important to consider how you will use the data you collect. This chapter provides an operations overview with these main sections:

- ❑ Serial Communications
- ❑ Hardware Components
- ❑ Setting Up Parameters

For more detailed information on operating and timing parameters, refer to Chapter 6, "Using Serial Data Commands."

---

## Serial Communications

Nephelometer operation is controlled through the serial communications interface (COM PORT). After the Nephelometer is configured, data can be logged by a computer connected to the serial interface (COM PORT).

Data gathered by a computer connected to the serial interface provides the greatest flexibility in the post processing of data.

---

## Hardware Components

This section describes some of the main hardware components of the Nephelometer from an operations standpoint. Refer to Chapter 8, "Performing Maintenance," for information on taking the Nephelometer apart and cleaning or replacing components.

### Power Failures

The Nephelometer is designed to be a long term monitoring device with watch-dog circuitry. This circuitry allows the Nephelometer to restart itself automatically in the event of a power failure or disruption. Operating parameters are preserved in battery backed-up RAM so that you do not have to reconfigure the Nephelometer after a power failure.

## **Lamp Power**

Lamp power is typically set to 75 watts power. This setting allows for a lamp life of approximately 2000 hours.

In some cases, like airborne measurements of short duration, when fast response time (short averaging time) and high sensitivity are desired, you can increase the lamp power. Increasing the lamp power to 90 watts will improve sensitivity but will reduce the lamp life.

## **Photomultiplier Tubes (PMTs)**

The Nephelometer PMTs convert scattered photons of light into electronic pulses that can be counted by internal processing electronics. The gain of each PMT is controlled by an applied voltage between 0 and 1200 volts DC. The voltages are set through serial data commands sent to the Nephelometer through its communications interface. PMT voltages are set at TSI for optimal sensitivity and usually do not require adjustment.

Generally, increasing the voltage increases the gain, which results in higher sensitivity. However, increasing the voltage also increases the “dark current” of the PMT until a point is reached where the dark current is increasing faster than signal gain. As this occurs, sensitivity decreases.

The PMT dark current is also affected by temperature. If the temperature of the environment is much above 70° F (21° C), as the dark current increases, sensitivity decreases and the PMT voltage may not be optimal. In this case, voltages may need to be lowered.

## **Heater**

If the Nephelometer is installed in an air conditioned area and is sampling hot, humid air, use the Nephelometer heater to prevent condensation. Generally the Nephelometer produces enough heat to maintain a sample temperature that is equal to or above the inlet temperature. However, if the sample temperature is below that of the inlet temperature, condensation may occur, which affects particle scattering inside the instrument.

## **Blower Settings**


Unless the blower is turned off, the Nephelometer provides full power to the blower during an auto zero cycle. This allows the Nephelometer to have the blower set to a lower power during measurement (reducing contamination build-up) and a fast purge during the auto zero cycle.

Response time is the time required for the Nephelometer to respond to a step change in aerosol scattering. Response time is a combination of the time for the actual aerosol to change within the instrument (as low as a few seconds) and the averaging time. The speed setting of the blower and restrictions in the flow path determine how long it actually takes the aerosol to be exchanged.

## **Humidity and Temperature Sensors**

The Nephelometer has one humidity and two temperature sensors. The humidity sensor and one temperature sensor are located near the sample outlet and the other temperature sensor is located at the sample inlet.

The outlet temperature sensor is positioned as close to the sample volume as possible to provide an accurate reading of the sample temperature. The outlet temperature sensor is used in the calculations as the sample temperature.



**Figure 3-1**  
Humidity and Temperature Sensor Locations

---

## Setting Up Parameters

This section discusses some considerations you should make before collecting data with the Nephelometer. This section contains these main subsections:

- Operating Parameters
- Timing Parameters
- Polled and Unpolled Communications

## Operating Parameters

All of the Nephelometer's operating parameters are set with the serial interface through the COM PORT. While many of the operating parameters can be adjusted using a computer and the DOS-based software described in Chapter 4, all of the parameters can be set using the firmware commands described in Chapter 6 or by using the terminal emulation program on a computer.

**Note:** *Although many of the operating parameters are mentioned in this section, refer to Chapters 4 and 6 for instructions on actually setting the parameters.*

## Sensitivity

The aerosol scattering coefficient is the primary measurement of the Nephelometer. The coefficient is the light scattering by particles in units of inverse meters. Particle scattering is measured by taking the total scattering value of the sample and then subtracting contributions of scattering from air molecules (Rayleigh scattering) and the instrument background. Ultimately, the maximum sensitivity that can be achieved is relative to the amount of light detected or captured from the sample.

Many factors relating to sensitivity have been addressed during the design of the Nephelometer. Sample averaging time and lamp power are the two primary user-controlled parameters affecting sensitivity.

## Backscatter Measurements

Nephelometer Models 3561 and 3563 are equipped with a backscatter shutter. The motor-operated shutter either runs synchronously with the reference chopper (in backscatter mode) or is parked (in total scatter mode).

When operating in backscatter mode, the Nephelometer is measuring total scatter half of the time and measuring the back scatter fraction the other half. In backscatter mode, Nephelometer's total scatter sensitivity drops by 0.7 [ $1/\sqrt{2}$ ] because the sample time is effectively cut in half. For this reason the backscatter mode should be used only if backscatter data is desired.

Vibration can also be a consideration in choosing backscatter mode. During backscatter operation, vibration or rapid shifts in the position of the Nephelometer (when airborne, for example) can cause a loss of synchronization between the reference chopper and the backscatter shutter. No data is taken until synchronization is reestablished. Although loss of motor synchronization can be monitored via the status flags accessible through the serial interface (see RF command in Chapter 6), use the backscatter mode only when backscatter data is desired.

### **Zero Baseline Measurements**

At times, the light scattering from particles can be very small compared to that of the air and the instrument background. The Nephelometer has built-in pressure and temperature sensors that provide compensation for the changes in the scattering of air as its density changes. This compensation, however, cannot correct for the value of the instrument's background scattering, which will also change over time as dirt and contamination from the sample deposit on the walls of the Nephelometer.

The Nephelometer has a HEPA particle filter and a motorized valve that provide the ability to accurately measure the scattering from air and the instrument background. This zero baseline measurement (sometimes referred to as zero measurement) occurs when the Nephelometer switches the motorized valve at the inlet, causing all of the sampled air to flow through the high efficiency filter. Particles are eliminated from the air flowing through the Nephelometer and the scattering value measured by the Nephelometer is from air and instrument scatter.

Particle-free air is measured for several minutes to obtain the zero baseline scattering value. The valve is again switched to remove the filter from the sample flow. Normal sampling resumes. The zero baseline value is used to obtain the actual instrument scatter. This value along with the scatter of air (determined from its temperature and pressure) are subtracted from the sample's total scatter signal to determine the scattering from particles.



## **Autozero Modes**

Zero baseline measurements can be initiated in three ways as determined by the autozero mode setting of the Nephelometer: manual, normal, and air-chop. Refer to the SMZ command in Chapter 6.

### **Manual Mode**

In manual mode, an automatic zero measurement occurs only when initiated by a direct command to the communications port (see the Z command in Chapter 6).

Once the command is received, the Nephelometer changes the position of the filtered air valve, measures the zero baseline value, and then repositions the valve for normal sampling.

### **Normal Mode**

In normal mode, the Nephelometer makes an auto zero measurement at preset intervals (see "Time Parameters" in this chapter).

### **Air-Chop**

Air chop mode is similar to normal mode, but the zero baseline calculation is based on two or more auto zero cycles. This mode is generally used in very "clean" (low particle scatter) environments.

## **Time Parameters**

This section discusses the several user-controlled time parameters for Nephelometer operation. Also refer to Figures 6-2 and 6-3 in Chapter 6. These figures detail the Nephelometer timing sequences and the relationships between the time parameters described in this section.

## **Averaging Time**

The most important parameter is the “sample averaging time” or just “averaging time.” The averaging time value can be set between 1 and 9999 seconds by using the software or by using the STA command (Chapter 6).

Longer averaging times increase the Nephelometer's sensitivity at the sacrifice of response time. Sensitivity increases by the square root of the averaging time, so quadrupling of the sample time results in a doubling of the sensitivity. Typically, averaging times of less than 30 seconds are used only for testing purposes, or if the data is to be post-processed.

## **Zero Time**

The “zero time” setting is the time the Nephelometer takes to make a zero baseline measurement. This time, which can be set between 1 and 9999 seconds, is typically 300 seconds (see the STZ command in Chapter 6). Generally the zero time should *not* be less than the sample averaging time setting.

## **Auto Zero Period**

The auto zero period occurs between subsequent zero measurements during normal or air chop modes. The auto zero period must be small enough to follow changes in instrument background due to the build up of contaminants. Generally a value of 3600 seconds is sufficient for ground-based measurements. For extremely clean environments where air chop mode is used, the auto zero period should be shorter to provide equal zero and measurement periods. See Table 3-1, “Typical Operating Parameter” in this chapter, and the STP command in Chapter 6.

## **Valve Blanking Time**

When a zero measurement cycle is initiated, there is a period of time, called blanking time, during which the filtered air valve turns and the Nephelometer switches between clean and sample air (or visa-versa).

During blanking time, the data is in transition and should not be used. Blanking time represents the number of seconds in which data is disregarded and not used for either the zero or measurement data. Typically, blanking time is set to 30 seconds, but use a longer blanking time if the flowrate through the Nephelometer is reduced causing a longer purge or transition time.

See the STB command in Chapter 6.

## **Polled and Unpolled Communications**

The Nephelometer communicates in two ways: polled mode and unpolled mode. There are advantages and disadvantages with either mode. In polled mode, the Nephelometer waits for a command from an external computer or terminal and responds only when a command is sent. (TSI Nephelometer software uses polled mode.) In unpolled mode, the Nephelometer automatically sends data to the communications port once every averaging period.

Although the Nephelometer updates data records as fast as once per second, this data is averaged in a running or boxcar average. The averaged data is completely updated once every averaging time period. Therefore, to avoid double averaging, the computer should poll the Nephelometer at a rate greater than the averaging time.

To reduce the burden on the external computer to provide accurate timing, the Nephelometer has an unpolled communications mode. In unpolled mode the Nephelometer automatically reports selected data at intervals equal to the averaging time.

When data is recorded in unpolled mode, you can use a shorter averaging time to provide more detail in the data (faster time response) at the expense of having more data to process. This also provides more flexibility in post processing of data as short averaged data can be combined into longer averages. This is useful in airborne work where several hours of fast response data may be taken. Ground-based measurements may opt for longer averaging times to reduce the amount of data gathered over long-term operation.

The following chart represents some typical settings for the previously mentioned time parameters.

**Table 3-1**  
Typical Time Settings

<b>Parameter</b>	<b>Ground Based</b>	<b>Ground Based (clean)</b>	<b>Airborne</b>
Avg Time	60 to 300	1500	30
Zero	300	300	300
Blanking	30	30	30
AZ Period	3600	720	Manual
AZ Mode	Normal	Air Chop (SMZ = 5)	Manual

When using long averaging times (greater than 300 seconds), all timing parameters should be set to integer multiples of the boxcar size (see the STA command under "Set Commands" in Chapter 6). This minimizes data loss that may occur if autozero transitions occur during the middle of a measurement boxcar.

## CHAPTER 4

# Using Nephelometer Software (DOS)

This chapter contains installation and operating instructions for the DOS-based software, designed for exclusive use with the Model 3550/3560 Series Integrating Nephelometer.

This chapter contains these main sections:

- Hardware and software requirements
- Installing and starting the software
- Keys and conventions
- Menus and command options
- Performing Calibration

---

## Hardware and Software Requirements

To operate the Nephelometer software you need an IBM personal computer or 100% IBM-compatible computer with the following:

- 80386 or higher processor
- One hard drive with at least 100 Mb storage capacity (recommended for data logging)
- One 3.5-in. floppy drive
- 4 Mb or more of random-access memory (RAM)
- One RS-232 serial port for connection to the Nephelometer
- DOS version 3.3 or higher.

**Note:** *The DOS-based Nephelometer software may be unstable if used with multiple applications under Microsoft® Windows.*

---

## Installing and Starting the Software

To load the Nephelometer software on the computer hard drive:

1. Apply power to the computer.

**2.** Insert the 3.5 inch software disk in the floppy drive of the computer.

**3.** Go to the root directory on the C drive (hard drive):

```
C:\>
```

**4.** Make a subdirectory for the Nephelometer, like "nephdos":

```
C:\nephdos
```

**5.** Copy the contents to the software disk into the subdirectory. These files are copied:

```
NEPH.DAT
NEPHCALB.EXE
NEPHCFG.EXE
NEPHCNFG.DAT
NEPHMAIN.EXE
NEPHPOLL.EXE
NEPHSET.EXE
NEPHTERM.EXE
RAYSCAT.DAT
```

**6.** Start the software by typing:

```
nephmain.exe
```

The main menu screen is displayed (Figure 4-1).

```
+-----+
| TSI Nephelometer Software      |
+-----+
+-----Main Menu-----+
| Set Commands Menu             |
| Calibration                    |
| Data Collection (Polled)       |
| Configuration Menu            |
| Terminal Mode                 |
| Com Port and Data Path        |
| Exit                           |
+-----+
```

**Figure 4-1**  
Main Menu Screen

**Note:** The DOS-based Nephelometer software may be unstable if used along with multiple applications under Microsoft® Windows.

## Keys and Conventions

This section gives key functions, and command conventions for the Nephelometer software. Table 4-1 gives the keyboard functions; software command conventions are shown in Figure 4-1.

**Table 4-1**  
Keyboard Functions

Key	Function
Hot keys	Type the first letter of a menu option to cursor to that option.
<b>F10</b>	Press <b>F10</b> to quit either the calibration, terminal, or polled mode.
<b>Esc</b>	Press <b>Esc</b> to return at any time to the previous screen, or to return to the DOS prompt from the main menu screen.
<b>Arrow</b> keys (←↑→↓)	Press an <b>Arrow</b> key to move within a screen field, or to move up and down a list of menu options.
<b>Tab</b>	Press <b>Tab</b> to move from field to field within a screen.
<b>Enter</b>	Use <b>Enter</b> in combination with other keys (like Arrow keys) or conventions (like <Read>) to select an option, to enter a parameter, or to read a parameter.

```

+-----+
|           < Set >           < Read >           < Cancel >           |
+-----+


```

**Figure 4-2**  
Software Command Conventions

- <Set>      Sends a parameter, set by the user, to the Nephelometer.
- <Cancel>   Returns to the previous screen without entering any selections or new parameters.
- <Read>     Reads the parameter that is currently being used by the Nephelometer.

# Menus and Command Options: A Reference Guide

This section describes the software menus and command options in the order they appear on menu screens. Figure 4-3 shows a structural overview with the three menus (Main, Set Commands, and Configuration).



**Figure 4-3**  
Software Menu Structure



## Main Menu

The Main Menu screen is displayed when you start the Nephelometer software (Figure 4-4). When you select any menu option, a message describing the option is displayed at the bottom of the screen.

**Note:** Start the software by typing `nephmain.exe` under the subdirectory created for the Nephelometer program.

```
+-----+
| TSI Nephelometer Software v1.x |
+-----+

+-----Main Menu-----+
| Set Commands Menu      |
| Calibration            |
| Data Collection (Polled)|
| Configuration Menu     |
| Terminal Mode          |
| Com Port and Data Path |
| Exit                   |
+-----+
```

**Figure 4-4**  
Main Menu Screen

## Main/ Set Commands Menu

The Set Commands Menu screen (Figure 4-5) contains options that allow you to set up the operating parameters for the Nephelometer. These options correspond to some of the set commands described in Chapter 6.

**Note:** You can use the Terminal option under the Main Menu to enter any serial data commands for parameters not included in the Set Commands Menu.

Go to the Set Commands Menu by selecting the Set Up Nephelometer option from the Main Menu screen. When you select an option from the Set Commands Menu, a message describing the option is displayed at the bottom of the screen.

```
+-----Set Commands Menu-----+
| K Calibration Constants          |
| Calibration Label String        |
| Mode for BackScatter Shutter    |
| Lamp Power                      |
| Time for Averaging              |
| Time and Date                   |
| Voltage on Photomultiplier Tube |
| Exit                            |
+-----+
```

**Figure 4-5**  
Set Commands Menu

## Set Command/ K Calibration Constants Option

The K Calibration Constants option (Figure 4-6) allows you to choose calibration constants for each of the three detection colors (blue, green, and red). The K calibration constants, 1 through 4, are further defined in the SK command under "Set Commands" in Chapter 6. In order to read the current K Calibration Constants for a particular channel (blue, green, or red), you have to select the channel and then <Read> followed by <CR>.

**Note:** Models 3551/3561 only require calibration constants for the green channel.

```

+----- K Constants -----+
| K1 (1-65535 picoseconds):  |
|+-----+                   |
| 0                           |
|+-----+                   |
| K2 (Total Scatter):        |
|+-----+                   |
| 0.000E+0                   |
|+-----+                   |
| K3 (Rayleigh Scatter):     |
|+-----+                   |
| 0.000E+0                   |
|+-----+                   |
| K4 (Backscatter):          |
|+-----+                   |
| 0.000E+0                   |
|+-----+                   |
|           -- Channel --    |
| ( ) Blue                   |
| ( ) Green                   |
| ( ) Red                     |
|+-----+                   |
| < Read >   < Set >   < Cancel > |
|+-----+                   |

```

**Figure 4-6**  
K Calibration Constants Option

### **Set Command/ Calibration Label String Option**

The Calibration Label String option allows you to send an identifying label for the calibration constants to the Nephelometer (Figure 4-7). As an example, "Last Calibration 1/30/93". The Calibration Label String option corresponds to the SL command under "Set Commands" in Chapter 6.

```

+----- Calibration Label -----+
|
| Label String (up to 50 char.):
|+-----+
|SL
|
+-----+
|
+-----+
|           < Set >           < Read >           < Cancel >
|
+-----+

```

**Figure 4-7**  
Calibration Label Option

### **Set Command/ Mode for BackScatter Shutter Option**

The Mode for BackScatter Shutter option (Figure 4-8) allows you to disable the backscatter shutter for total scatter measurements, or to enable the backscatter shutter so that you can take total scatter and backscatter measurements. This option corresponds to the SMB command under "Set Commands" in Chapter 6.

Refer to the Backscatter Measurements section in Chapter 3, "Operations Overview," to determine how to set this option.

```

+-- BackScatter Shutter Mode --+
|
|           -- Mode --
|  ( ) Total Scatter Only
|  ( ) Backscatter
|
+-----+
| < Set > < Read > < Cancel >
|
+-----+

```

**Figure 4-8**  
BackScatter Shutter Mode Option

### Set Command/ Lamp Power Option

The Lamp Power option (Figure 4-9) allows you to set the intensity of the halogen lamp, which is used as a light source inside the Nephelometer. The setting for the Lamp Power option directly affects the life of the lamp. The Lamp Power option corresponds to the SP command under "Set Commands" in Chapter 6.

Refer to the Lamp Power section in Chapter 3, "Operations Overview," to determine how to set this option.

```
+----- Lamp Power -----+
|                               |
| Power in watts (0 - 150):   |
|                               |
| +-----+                   |
| |0|                               |
| +-----+                   |
|                               |
+-----+
| < Set > < Read > < Cancel > |
+-----+
```

**Figure 4-9**  
Lamp Power Option

### Set Command/ Time for Averaging Option

The Time for Averaging option allows you set the averaging time, in seconds, over which the sample is measured (Figure 4-10). The Time for Averaging option corresponds to the STA command under "Set Commands" in Chapter 6.

Refer to the Averaging Time section in Chapter 3, "Operations Overview," to determine how to set this option.

**Note:** *Averaging times in the valid range that do not fall on a boxcar boundary, will be rounded up to the next valid boxcar boundary averaging time. That is, if set to 604, because the boxcar size is 3 seconds, the next valid averaging time is 606. The Nephelometer uses this number.*

```

+----- Averaging Time -----+
|                               |
|   Time in seconds (1 - 9960): |
|                               |
|   +-----+                 |
|   | 0      |                 |
|   +-----+                 |
|                               |
+-----+
| < Set > < Read > < Cancel > |
+-----+

```

**Figure 4-10**  
Time for Averaging Option

### **Set Command/ Time and Date Option**

The Time and Date option allows you to set the time and date for the Nephelometer's internal clock (Figure 4-11). The Time and Date option can be accessed with the STT command, the RT command, or in unpolled mode with the UT command. Refer to the "Commands" section in Chapter 6.

```

+----- Date/Time -----+
|                               |
|   Date:                       |
|   +-----+                 |
|   |05-10-1994   |                 |
|   +-----+                 |
|                               |
|   Time:                       |
|   +-----+                 |
|   |13:20:47     |                 |
|   +-----+                 |
|                               |
+-----+
| < Set > < Read > < Cancel > |
+-----+

```

**Figure 4-11**  
Time and Date Option

## Set Command/ Voltage on Photomultiplier Tube Option

The Voltage on Photomultiplier Tube option (Figure 4-12) allows you to set the voltage that drives each PMT tube, blue, green, and red. The range for this option is 1-1200 volts. This option corresponds to the SV command under "Set Commands" in Chapter 6. In order to read the current Voltage on Photomultiplier Tube option for a particular channel (blue, green, or red), you have to select the channel and then <Read> followed by <CR>.

**Note:** Models 3551/3561 only require voltages for the green tube.

Refer to the Photomultiplier Tubes section in Chapter 3, "Operations Overview," to determine how to set this option.

```
+----- PMT Voltage -----+
|           -- Color --           |
| ( ) Blue                         |
| ( ) Green                        |
| ( ) Red                          |
|                                   |
|           PMT Voltage            |
| +-----+                        |
| |0                                     |
| +-----+                        |
|                                   |
+-----+
| < Read > < Set > < Cancel > |
+-----+
```

**Figure 4-12**  
Voltage on Photomultiplier Tube Option

## **Main/ Calibration Option**

Refer to the "Performing Calibration" section at the end of this Chapter.

## **Main/ Data Collection (Polled) Option**

The Data Collection (Polled) option allows you to read data records from the Nephelometer and display the records in a format that is easy to interpret (Figure 4-13). This option also enables you to store the data, in an unpolled record format, to a specified file as it is received from the Nephelometer. This format is displayed in the "Unpolled Record Formats" section in Chapter 6. The specified file is saved to a directory set by the Com Port and Data Path option (refer to that option in this chapter).

**Note:** *The display screen is updated approximately once per second but data is written to the data files at intervals set by the averaging time (see the "Time for Averaging Option" in this chapter).*



```

+-----Data Record-----+-----Status Record-----+
|Current Mode:   Normal Measurement   ||Green Cal (in Hz)           3.196E+05 | |
|               BackScatter           ||                             ||
|Time Remaining: 129                 ||Pressure (MB)               971.3 |
|(current state)                    ||Sample Temp. (K)            300.4 |
|                                       ||Inlet Temp. (K)             297.6 |
|               Scattering Coefficient ||Rel. Humidity(%)            33.8 |
|   Total Scatter   Back Scatter     ||Lamp Voltage                 12.9 |
|                                       ||Lamp Current                  5.9 |
|B  -4.98700E-07    1.08700E-07     ||BNC Input Voltage            16 |
|G   4.71600E-08   -1.45200E-07     ||Status Flags                  OK |
|R  -9.66400E-08   -1.83400E-08     ||                             ||
|                                       ||                             ||
+-----+-----+
+-----Photon Frequency in Hz (Blue, Green, Red)-----+
|               Total Scatter                BackScatter                Sample|
|               Cal   Meas   Dark            Cal   Meas   Dark            Press Temp.|
|               |               |               |               |               |               |
|B  187332      1205      13            115844      601       10            971.3 300.4|
|G  317605      1031      26            197815      493        20            971.3 300.4|
|R  195066       534      219           120578      404        209           971.3 300.4|
|               |               |               |               |               |               |
+-----+-----+
(L)og Data   Log (F)ilename   (A)vg Time   Mode           F10-Quit
  OFF        LOGDATA.LOG      60           POLLED        9/20/1994 10:14:40

```

**Figure 4-13**  
Data Collection Polled Option

## Main/ Configuration Menu

The Configuration Menu (Figure 4-14) contains options that allow you to display current configuration parameters, and to read the configuration file from either a disk or from the Nephelometer. The options under the Configuration Menu are useful if the Nephelometer configuration data is ever lost or corrupted. The options allow you to quickly analyze a situation because you are able to view all the parameters as a group. You can also store, then easily retrieve different configurations.

**Note:** *You cannot change the contents of Configuration screens using Configuration Menu options, but you can make changes using either the Terminal Mode option or using options under the Set Commands Menu*

```
+-----Configuration Menu-----+
| Display Current Config. Data      |
| Read Config. Data from Nephelometer |
| Read Config. Data from Disk File  |
| Save Config. Data to Nephelometer  |
| Save Config. Data to Disk File    |
| Print Config. Data                |
| Exit                              |
+-----+
```

**Figure 4-14**  
Configuration Menu

## Configuration/ Display Current Config. Data Option

The Display Current Config. Data option displays the configuration data in the software buffer, either data read from the Nephelometer or data from the disk file (Figures 4-15 and 4-16).

```

+-----Nephelometer Configuration Data-----+
|Analog Output Range:           Zero V   Full Scale   |
|      Relative Humidity  ->      0 -     0           |
|      Barometric Pressure ->      0 -     0           |
|      Sample Temperature ->      0 -     0           |
|      Scattering         ->      0 -     0           (Offset=0.000E+00) |
|                                                                           |
|Channel  Type                               Channel  Type   |
|   0     0 (External PC Control)           5     0 (External PC Control) |
|   1     0 (External PC Control)           6     0 (External PC Control) |
|   2     0 (External PC Control)           7     0 (External PC Control) |
|   3     0 (External PC Control)           8     0 (External PC Control) |
|   4     0 (External PC Control)           |
|                                                                           |
|Calibration Points:           Low Bits  Low Value   High Bits  High Value |
|      Relative Humidity  ->           0         0         0         0   |
|      Barometric Pressure->           0         0         0         0   |
|      Sample Temperature ->           0         0         0         0   |
|      Inlet Temperature  ->           0         0         0         0   |
|                                                                           |
|Calibration Constants:   K1         K2         K3         K4   |
|      Blue ->            0  0.00000E+00  0.00000E+00  0.00000E+00 |
|      Green->            0  0.00000E+00  0.00000E+00  0.00000E+00 |
|      Red  ->            0  0.00000E+00  0.00000E+00  0.00000E+00 |
+Screen 1 of 2-----Press any key to Next Screen-----+

```

**Figure 4-15**  
Display Current Config. Data Option

```

+-----Nephelometer Configuration Data-----+
|Data Delimiter: Comma
|
|Calibration Label:
|
|Backscatter Shutter Mode:          Total Scatter Mode ONLY
|Autozero Baseline Measurement:    Manual Mode
|
|Lamp Power (Watts)                 :      0
|Auxiliary BNC Output (millivolts):  0
|
|Time for Averaging (sec.)          :      0
|Time to Blank Valve (sec.)         :      0
|Time between Autozeros (sec.)     :      0
|Time for AutoZero Measurement (sec.):  0
|
|PMT Voltage: Blue =      0
|                Green =    0
|                Red  =    0
|
|Zero Baseline Values:  Total Scatter      BackScatter      Rayleigh Scatter
|                Blue ->  0.000E+00      0.000E+00      0.000E+00
|                Green -> 0.000E+00      0.000E+00      0.000E+00
|                Red  ->  0.000E+00      0.000E+00      0.000E+00
+Screen 2 of 2-----Press any key to return to Menu-----+

```

**Figure 4-16**  
Display Current Config. Data Option

### **Configuration/ Read Config. Data from Nephelometer Option**

The Read Config. Data from Nephelometer option allows you to read the configuration data residing in the Nephelometer to the software buffer. This offers a tool to help you quickly analyze any operating problems.

### **Configuration/ Read Config. Data from Disk File Option**

The Read Config. Data from Disk File option allows you to read the contents of the disk file into the software buffer. You can then select “Display Current Config. Data” to review this data before you send it to the Nephelometer.

### **Configuration/ Save Config. Data to Nephelometer Option**

The Save Config. Data to Nephelometer allows you to save the configuration data residing in the software buffer to the Nephelometer. The operating parameters in this file become the operating parameters for the Nephelometer.

### **Configuration/ Save Config. Data to Disk File Option**

The Save Config. Data to Disk File option allows you to save the configuration data residing in the software buffer to a disk file for easy access should the Nephelometer’s configuration become corrupt.

### **Configuration/ Print Config. Data Option**

Sends a copy of the Config. Data to the printer attached to your parallel port.

### **Main/ Terminal Mode Option**

The Terminal Mode option provides a basic serial interface to the Nephelometer, where you can type in serial data commands and send them directly to the Nephelometer. This option is useful for troubleshooting Nephelometer operating parameters, the calibration process, or configuration data.

Examples of the commands you might enter in terminal mode are commands that are normally set once, and not changed during the course of taking sample data. These include any of the set commands for mapping analog outputs (like SA, SB, and SX), action commands that control the mechanical components, read commands and unpolled commands (refer for Chapter 6, "Using Serial Data Commands").

## Main/ Com Port and Data Path Option

The Com Port and Data Path option allow you select the serial communications port of the computer and to designate a data path for storage of data files (Figure 4-17). Use DOS conventions in naming a directory path. Any files stored will go to the directory set by this option.

```
+----- Data Path -----+
|                          |
|              COM Settings |
|                          |
|  ( ) Com1                |
|  ( ) Com2                |
|                          |
| Directory Path For Data Files (*.dat): |
|                          |
| +-----+ |
| |C:\SOFT | |
| +-----+ |
|                          |
|-----+-----+
|              < OK >       < Cancel > |
|-----+-----+
```

**Figure 4-17**  
Com Port and Data Path Option

---

## Performing Calibration

You should calibrate the Nephelometer before an intensive experiment, calibrate periodically to verify no drift has occurred, calibrate if the reference chopper is dirty or scratched, or if you clean the chopper as part of periodic maintenance.

The Calibration option under the Main menu allows you to easily calibrate the Nephelometer using two span gases. The Calibration option allows you to compare the results between air (low span) and CO<sub>2</sub> (high span).

### Hardware Setup

Set up the gas before begin you begin calibration. Gas should be approximately room temperature as it reaches the Nephelometer. One way to ensure this is to use a 6- to 9-meter (20- to 30-foot) length of 6-mm (¼-in.) diameter tubing from the CO<sub>2</sub> tank to the Nephelometer. (The length of the tubing affects the temperature of the CO<sub>2</sub>.) Gas of 99.9% purity or better is typically used with ambient air for calibration.



#### Caution

Use calibration gases only in a well-ventilated area or exhaust the gases outside. Many gases used for calibration can cause asphyxiation if used in a confined area.

### Using the Main/ Calibration Option

Calibrate the Nephelometer by selecting the Calibration command from the Main Menu . The interactive calibration screen is displayed (Figure 4-18).

- The top section of the calibration screen displays the current calibration settings stored in the Nephelometer and any new calibration settings not yet saved to the Nephelometer.
- The middle section of the screen displays calibration data in frequency (Hz) for low span and high span gas in each wavelength (blue, green, and red). Calibration data is displayed for both total scatter and backscatter modes.

- The bottom section of the screen lists the steps you will use to calibrate the Nephelometer. Perform each step in the order given. As you complete each step, press <Enter>. A check (✓) is displayed next to each command as the command is executed.

```

+-----Calibration Parameters Table-----+
|          Current(B)  New(B)          Current(G)  New(G)          Current(R)  New(R)          |
|K2      4.325E-03     0.000E+00     4.484E-03     0.000E+00     4.614E-03     0.000E+00     |
|K4      9.533E-03     0.000E+00     9.717E-03     0.000E+00     1.035E-02     0.000E+00     |
+-----+
+-----Calibration Data in Hz-----+
|          Total Scatter          BackScatter          |
|          Cal      Meas      Dark      Cal      Meas      Dark      |
|LOW SPAN GAS          |
|Blue  0.00000E+00  0.00000E+00  0.00000E+00  0.00000E+00  0.00000E+00  0.00000E+00  |
|Green  0.00000E+00  0.00000E+00  0.00000E+00  0.00000E+00  0.00000E+00  0.00000E+00  |
|Red    0.00000E+00  0.00000E+00  0.00000E+00  0.00000E+00  0.00000E+00  0.00000E+00  |
|HIGH SPAN GAS          |
|Blue  0.00000E+00  0.00000E+00  0.00000E+00  0.00000E+00  0.00000E+00  0.00000E+00  |
|Green  0.00000E+00  0.00000E+00  0.00000E+00  0.00000E+00  0.00000E+00  0.00000E+00  |
|Red    0.00000E+00  0.00000E+00  0.00000E+00  0.00000E+00  0.00000E+00  0.00000E+00  |
+-----+
+-----Calibration Steps-----Elapsed Time(s):---+
1. set low span gas      (Air      )          9. feed CO2 into line   (manual)
2. set high span gas    (CO2     )          A. start data collection-high gas
3. switch valve filter                                     B. clear buffer
4. start data collection-low gas                          C. store high span gas data
5. clear buffer                                           D. write calibration label
6. store low span gas data                                E. save cal data (K2 & K4) to Neph
7. toggle blower to off (OFF)                            F. zero w/ clean air
8. plug inlet & outlet (manual)                          F8-Avg Time F9-Restart F10-Quit

```

**Figure 4-18**  
Calibration Option

**Note:** If you are using air as the low span gas and CO<sub>2</sub> as the high span gas, perform the calibration steps in the order given. If you are using a gas other than air as the low span gas (like helium), refer to calibration process given at the end of the Calibration Procedure.

**Table 4-2**  
Calibration Procedure

Step	Explanation
1. Set low span gas	Select to choose the low span gas. Air is the default. When you select step 1, a menu is displayed with different gases.
2. Set high span gas	Select to choose a high span gas, or a gas with a higher scattering coefficient than the low span gas, such as CO <sub>2</sub> . When you select step 2, a menu is displayed with different gases.


(continued)



**Table 4-2**  
Calibration Procedure (continued)

<b>Step</b>	<b>Explanation</b>
<b>3.</b> Switch valve filter	Select to switch the filtered air valve to the zero (closed) position and to switch the blue DQ filter in-line with the air inlet. When you select step 3, a blanking time begins. The blanking time is the time it takes the filtered air valve to move to the zero position plus the time it takes to purge the Nephelometer with filtered air. The blanking time is set using the STB command (see Chapter 6).
<b>4.</b> Start data collection-low gas	Select to begin low span gas data collection. After you select step 4, data collection begins after the blanking time described in step 3. The low span gas area of the calibration screen is updated, but not the high span area.  If you want to observe the filtered air purge, press <F8> to select a smaller averaging time of 5 seconds. The number in the Meas field should fluctuate, not rise and fall. When the Meas field is stable, press <F8> to select a 300-second averaging time and proceed to step 5.
<b>5.</b> Clear buffer	Select to clear the Nephelometer internal buffers of the data collected during step 4. Wait at least 300 seconds (see the Elapsed Time(s) field), then proceed to step 6.
<b>6.</b> Store low span gas data	Select to store in the software, the 300-seconds of low span gas data collected in step 5. Selecting step 6 causes the low span gas area of the calibration screen to freeze and low span data is no longer updated.
<b>7.</b> Toggle blower to off (OFF)	Select to turn off the blower. High span gas should be regulated at the source.  <b>Note:</b> Step 7 is a toggle function and can also be used to toggle on the blower.
<b>8.</b> Plug inlet & outlet (manual)	Remove the top cover of the Nephelometer using the procedure in Chapter 8. Plug the aerosol inlet and outlet of the Nephelometer (Figure 4-19).
<b>9.</b> Feed CO2 into line (manual)	Remove the tubing from the inlet of the blue DQ filter and connect the tubing from the high span gas to either the inlet tubing or to the inlet block (Figure 4-19). Turn on the gas until the flowrate is approximately 5 lpm. You should feel gas exiting through the DQ filter in the opposite direction of the arrow on the filter. Proceed to step A.
<b>A.</b> Start data collection-high gas	Select to begin high span gas data collection. The high span gas area of the calibration screen is updated, but not the low span area.  If you want to observe the high span gas purge, press <F8> to select smaller averaging time of 5 seconds. The number in the Meas field should fluctuate, not rise and fall. When the Meas field is stable with minimal fluctuation, press <F8> to select a 300-second averaging time and proceed to step B.

(continued)



**Figure 4-19**  
Plugs for Inlet and Outlet and Gas Line With the Blue DQ Filter

**Table 4-2**  
Calibration Procedure (*continued*)

<b>Step</b>	<b>Explanation</b>
<b>B.</b> Clear buffer	Select to clear the Nephelometer internal buffers of the data collected during step A. Remain at this step for at least 300 seconds (see the Elapsed Time(s) field), then proceed to step C.
<b>C.</b> Store high span gas data	Select to store to the software the 300 seconds of high span gas data collected in step B. Selecting step C causes the high span gas area of the calibration screen to freeze and high span data is no longer updated. <b>Note:</b> <i>At this time, new calibration constants (K2 and K4) are calculated and the top section of the screen is updated. If the new constants are not satisfactory, you can recalibrate without changing the constants, as long as you do not perform step E.</i>
<b>D.</b> Write calibration label	Select to write an identifying label for the calibration performed. The label can include the serial number of the Nephelometer, the calibration date, the low span and high span gases used, and the initials of the person performing the calibration.

(*continued*)

**Table 4-2**  
Calibration Procedure (*continued*)

<b>Step</b>	<b>Explanation</b>
<b>E.</b> Save cal data (K2 & K4) to Neph	<p>If the new K2 and K4 calibration constants are satisfactory, select E to write the new K2 and K4 constants to the Nephelometer.</p> <p>After selecting step E,</p> <ol style="list-style-type: none"> <li><b>1.</b> Remove the high span gas and the inlet and outlet plugs from the Nephelometer.</li> <li><b>2.</b> Reconnect the DQ filter to the inlet block, with the arrow of the filter pointing away from the inlet block.</li> <li><b>3.</b> Replace the top cover of the Nephelometer, attaching the cover with the four screws.</li> <li><b>4.</b> Select 7 to toggle on the blower and to purge the Nephelometer of high span gas.</li> </ol> <p><b>Note:</b> <i>To allow the Nephelometer to be completely purged of high span gas, wait at least 600 seconds with the blower on before proceeding to step F.</i></p>
<b>F.</b> Zero w/clean air	<p><i>Optional</i>—Select to purge the Nephelometer with clean air and to perform a zero background measurement.</p> <p>&lt;<b>F9</b>&gt; resets the check marks. &lt;<b>F10</b>&gt; exits the calibration.</p> <p>Go to the Main Menu and choose the Configuration option for options to store the calibration file to disk.</p>
<p>If you are using a low span gas other than air, perform the calibration steps in this order:</p>	<ol style="list-style-type: none"> <li><b>1.</b> See step 1.</li> <li><b>2.</b> See step 2.</li> <li><b>3.</b> See step 3.</li> <li><b>7.</b> See step 7.</li> <li><b>8.</b> See step 8.</li> <li><b>9.</b> Remove the tubing from the inlet of the blue DQ filter and connect the tubing from the low span gas to either the inlet tubing or to the inlet block (Figure 4-19). Turn on the gas until the flowrate is approximately 5 lpm. You should feel gas exiting through the DQ filter in the opposite direction of the arrow on the filter.</li> <li><b>4.</b> See step 4.</li> <li><b>5.</b> See step 5.</li> <li><b>6.</b> See step 6.</li> <li><b>9.</b> See step 9.</li> <li><b>A.</b> See step A.</li> <li><b>B.</b> See step B.</li> <li><b>C.</b> See step C.</li> <li><b>D.</b> See step D.</li> <li><b>E.</b> See step E.</li> <li><b>F.</b> See step F.</li> </ol>



## CHAPTER 5

# Using Nephelometer Software (Windows)

This chapter describes the Windows™ version of the Nephelometer software including:

- ❑ The setup procedure that loads the program onto your computer.
- ❑ A “Quick Start” section to walk you through the steps necessary to collect data and perform other basic operations.
- ❑ A reference section to describe each program menu and its commands.

This chapter assumes you have Microsoft® Windows™ on your computer and that you are familiar with how Windows works. If you are not familiar with Windows, please refer to the documentation and other information that came with Windows before you load and use this program.

The Nephelometer software is provided in a Windows-based format because Windows offers a popular and easy-to-use interface to work from.

---

## Hardware and Software Requirements

To use this software program you will need an IBM personal computer or 100% IBM-compatible computer with the following minimum features, components, and software:

- ❑ A 386 microprocessor (80386) or better.
- ❑ A hard drive large enough to accommodate Windows, the Nephelometer program (approximately 2.5 Mb), and Nephelometer data files.
- ❑ One 3.5-in. floppy drive.
- ❑ 4 Mb or more of random access memory (RAM).
- ❑ A mouse or similar pointing device.
- ❑ An RS-232 serial interface port (in addition to the one that may be required for the mouse).
- ❑ MS-DOS, version 5.0 or higher.
- ❑ Microsoft Windows version 3.1 or Microsoft Windows for Workgroups version 3.11.

---

## Installing and Starting the Software

Set up the Windows-based Nephelometer software as follows:

1. Verify that Windows version 3.1 or Windows for Workgroups 3.11 is installed and running and that your MS-DOS version is at least 5.0. To verify your version of Windows, select *About Windows* from the Windows Help menu or type:

WINVER

at the DOS prompt and press <Enter>.

To verify your MS-DOS version, type:

VER

at the DOS prompt and press <Enter>.

2. Insert the Nephelometer setup disk in your 3.5-in. floppy drive (generally this is your A: or B: drive).
3. Start the setup program by selecting **RUN** from the File menu and typing the following in the command box:

A:\SETUP <Enter> (or B:\SETUP <Enter>).

**Note:** *The Nephelometer program is compressed on the setup disk, you cannot load it without running the Setup program.*

4. A setup screen appears and indicates the drive and directory in which the program will be installed (the default directory), Figure 5-1. If you want the program installed in another directory or on another drive, you must enter the drive and/or directory name.

**Note:** *You must select a hard drive for program installation. You cannot install the Nephelometer program onto floppy diskettes or run it from a floppy drive.*



**Figure 5-1**  
Select the Destination Directory for Program Installation

5. Select **Continue**. The program files are automatically decompressed and stored on your hard drive.


**Note:** *If you are re-installing the Nephelometer software and the NEPHCNFG.DAT file on your computer is different than the one on the setup disk, you will be asked if you want to overwrite the existing file. Select **Yes** if you want to overwrite the NEPHCNFG.DAT file. See “Neph Config Data” section later in this chapter for more information on NEPHCNFG.DAT.*

When setup is complete, the README.TXT file opens automatically. The README.TXT file contains information that could not be included in this manual. Read the file and then close it. You can access the file anytime using a text editor such as NOTEPAD.

The setup program creates a directory called NEPHWIN on your hard disk (assuming you accepted the default directory name) that contains the required program files. Refer to Appendix E for a list of program files and their contents.

The Setup program also creates a new group called TSI APPLICATIONS and an icon for the program called NEPH for Windows (Figure 5-2).

**Note:** *Before creating a TSI APPLICATIONS group, the Setup program checks for an existing one. If one is present, it adds the Nephelometer icon only.*



**Figure 5-2**  
Setup Creates a Program Group and Icon

Once the program icon has been installed, you can start the Nephelometer program by double-clicking on its icon.

---

## Quick Start

This section is designed to give you a quick introduction to the basic operation of the Nephelometer Windows program. Using these instructions, you will:

- Verify the Com Port setting of the Nephelometer.
- Start collecting sample data
- View the data as it is being collected and check the status of the Nephelometer.
- Pause and Restart data collection from the Nephelometer.
- Print a graph of the collected data.
- Perform a manual background (zero baseline) measurement.

Before you begin, verify the Nephelometer is powered up and connected to the computer.

Start the Nephelometer program by double-clicking on its icon. The main menu appears (Figure 5-3).






**Figure 5-3**  
The Main Menu of the Windows-Based Nephelometer Program

## Step 1. Verify Com Port

Choose the **Com Port** command from the Config menu. The COM Port Setup display is shown (Figure 5-4).



**Figure 5-4**  
The COM PORT Setup Display

Verify that the communications port indicated on this display is the one connected to your Nephelometer (refer to Chapter 2 if necessary). The default is COM2. Select another communications port if necessary.

**Note:** *Most computers have only two active communications ports, COM1 and COM2. Typically, COM1 is used to connect a serial mouse, so COM2 is most likely available for the Nephelometer connection.*

Select **OK** to exit this display.

## Step 2. Start Data Collection

To start collecting data you should first open a log file and then start collecting data.

1. Select the **Log Data** command from the File menu. When the Log Data Setup window appears, type in:


ONETEST.DAT

and press **OK**. Watch the lower right corner of the display, it will indicate Data Logging is to ONETEST.DAT.

2. Select the **Data Collection** command from the Run menu. Immediately, four windows are automatically opened and minimized on your screen (Figure 5-5) and data collection begins.

The icons that appear on the screen represent the Data, Graph, Photon Frequency, and Neph State commands available under the View menu.

**Note:** *If you receive a message that the program is unable to communicate with the Nephelometer, the wrong COM PORT may have been selected in Step 1. Verify cabling between the computer and the Nephelometer, then verify the correct COM PORT is selected and that the Nephelometer is on.*



**Figure 5-5**  
Choosing Data Collection from the Run Menu Starts Data Collection


### Step 3. View Data and Nephelometer Status

Follow these steps to view data and the status of the Nephelometer:

1. Double-click on the **Data Display** icon to display current data (Figure 5-6).


This display lets you view Nephelometer data as it is being collected, and the information displayed should resemble that shown in Figure 5-6.

Data is updated on this display at approximately one second intervals but data is written to the data files at intervals set by the averaging time. If no data is displayed after 60 seconds, check that the Nephelometer is running.



**Figure 5-6**  
Select Data Display to View Data as it is Collected by the Nephelometer


2. Double-click on the **Neph State** icon (you don't need to close the Data Display window) to display the Neph State window, Figure 5-7.



**Figure 5-7**  
The Nephelometer State Display Shows the Status of the Nephelometer

This window displays the status of various Nephelometer components. They should all be marked with checkmark (✓) to indicate they are functional. If any are marked with an ✗, then there is a problem with the component. Refer to Chapters 2, 6, 8 and 9 for information to help you troubleshoot the problem.


3. Double-click on the **Photon Frequency** icon (you don't need to close either the Data Display window or the Neph State window) to display the Photon Frequency window (Figure 5-8).



**Figure 5-8**  
The Photon Frequency Display Shows the Current Photon Frequency Calculations

This window displays the raw photon frequency calculated by the software from the raw photon counts collected by the Nephelometer.

4. Double-click on the **Graph Display** icon (you don't need to close any of the other windows) to display the information that has been collected by the Nephelometer so far (Figure 5-9).



**Figure 5-9**  
The Graph Display Shows Current Data

The graph displays the data that has been collected by the Nephelometer up to this point.

## Step 4. Print a Graph

There are two ways to print this graph. You can select either **Print** or **Print Screen** from the File menu.

Before doing either, however, select **Pause** from the main menu. You should always Pause data collection before printing.

Before using the Print Screen function, the Print Screen key must be pressed to capture the current active screen to the clipboard. The Print Screen function will print the entire screen, including menu bars, elevators and the like. It is being used to print log data graphs and may be suitable for draft prints or quick snaps of data.


When you choose **Print**, you are asked to select either a scatter data graph or a status graph and to enter a header message.

Choose **Scatter Data Graph** and enter the header message "This is test data collected by (enter your initials)." Then press **OK**. The data will be printed in landscape mode.

## Step 5. Perform a Background (Zero Baseline) Measurement

Zero baseline measurements are important so that the Nephelometer's background scattering can be deducted from sample data measurements. Normally, zero baseline measurements are performed automatically; however, to introduce you to the Terminal function of the Windows program, you should perform a manual zero baseline measurement as follows.

1. Select the **Terminal Mode** command from the Run menu. The Terminal Display appears, Figure 5-10. This display allows you to communicate with the Nephelometer using the Serial Data Commands, described in Chapter 6.)



```
RU
T,1994,10,19,14,16,37
B,2129112,49109,274,2769,1334311,24621,256,2769,963.6,300.0
G,3369857,39024,357,2769,2127474,19613,357,2769,963.6,300.0
R,2582601,27932,5172,2769,1611354,22169,5051,2769,963.6,300.0
D,ZBXX,59,2.555e-5,1.264e-5,6.282e-6,1.265e-5,6.215e-6,4.020e-6
Y,253120,963.6,300.0,296.8,29.3,12.8,5.9,18,0000
STA
300
SP
75
SK
ERROR
SKG
20000,3.923e-3,1.226e-5,.479
SKR
20000,3.557e-3,4.605e-6,.438
SKB
20000,3.940e-3,2.789e-5,.475
SL
SEPT 30 1994 5/N 1004 AIR CO2.
SMB
1
STT
1994,10,19,14,19,12
```

**Figure 5-10**  
The Terminal Mode Display

2. Enter a **Z** command (type the letter z) and press **<Enter>**. This begins a zero measurement.
3. Select the **Terminal Mode** command from the Run menu to close Terminal Display.
4. Look at the Data Display, refer to Figure 5-6. It should show the Time Remaining counting down toward 0 (zero). The current mode should first be changed to Blank Period, then to the Zero Measurement, and eventually to Normal Measurement when the time remaining reaches 0. This indicates that zero measurement has been taken and stored to the log data file.

## Basic Nephelometer Operations

Although it does not list all of the functions available in the Windows Nephelometer program, Table 5-1 lists many of the most basic operations you might want to perform. Use the instructions provided to get better acquainted with how the program works, or as a quick reference.

Although you may be comfortable using the program from these instructions, take the time to read the detailed information provided in the reference portion of this chapter so you are familiar with it when you have specific questions about program operation.

**Table 5-1**  
Basic Windows Nephelometer Operations

To...	Perform the Following Steps
Log data to a file	<ol style="list-style-type: none"> <li>1. Select <b>Log Data</b> from the File menu.</li> <li>2. Enter the filename and select <b>OK</b>.</li> <li>3. Select <b>Data Collection</b> from the Run menu.</li> </ol> <p><i>or</i></p> <p>Perform Step 3 then 1 then 2.</p>
Start collecting data (not logged to a file)	Select <b>Data Collection</b> from the Run menu.
Pause/resume data collection	Select <b>Pause</b> or <b>Resume</b> from the main menu.
View raw data currently being collected by the Nephelometer	<ol style="list-style-type: none"> <li>1. Verify Data Collection is proceeding.</li> <li>2. Double-click on the Data Display icon.</li> </ol>
View the Nephelometer configuration data	Select <b>Neph Config Data</b> from the View menu.
View a graph of current Nephelometer activity	<ol style="list-style-type: none"> <li>1. Verify Data Collection is proceeding.</li> <li>2. Double-click on the <b>Graph Display</b> icon.</li> </ol>
View data already saved to a file	<ol style="list-style-type: none"> <li>1. Verify Data Collection is OFF.</li> <li>2. Select <b>Open Log File</b> from the Log menu.</li> <li>3. Enter the name of the file and select <b>OK</b>.</li> </ol>
View a graph of data already saved (closed file)	<ol style="list-style-type: none"> <li>1. Verify Data Collection is OFF.</li> <li>2. Select <b>Open Log File</b> from the Log menu.</li> <li>3. Enter the name of the file and select <b>OK</b>.</li> <li>4. Select <b>Graph</b> from the display.</li> </ol>

(continued)



**Table 5-1**Basic Windows Nephelometer Operations (*continued*)

<b>To...</b>	<b>Perform the Following Steps</b>
Send serial data commands to the Nephelometer	<ol style="list-style-type: none"> <li>1. Select <b>Terminal Mode</b> from the Run menu.</li> <li>2. Type the Serial Data Command (see Chapter 6) and press &lt;<b>Enter</b>&gt;.</li> </ol>
Print a graph	<ol style="list-style-type: none"> <li>1. View a graph.</li> <li>2. Select <b>Print</b> from the File menu.</li> <li>3. Select <b>Info to Print</b> and <b>Header Message</b> (if desired)</li> <li>4. Select <b>Print</b>.</li> </ol>
Print whatever is displayed on the screen	<ol style="list-style-type: none"> <li>1. Press &lt;<b>Print Screen</b>&gt; key.</li> <li>2. Select <b>Print Screen</b> from the File menu (or press &lt;<b>F9</b>&gt;).</li> </ol>
Run a background (zero baseline) measurement manually	<ol style="list-style-type: none"> <li>1. Verify Data Collection is OFF or paused.</li> <li>2. Select <b>Terminal Mode</b> from the Run menu.</li> <li>3. Type z and press &lt;<b>Enter</b>&gt;.</li> <li>4. Close the Terminal Display.</li> <li>5. Resume Data Collection.</li> </ol>
Get online help	Select the <b>Contents</b> command from the Help menu <i>or</i> Press < <b>F1</b> >.
Exit the program	<ol style="list-style-type: none"> <li>1. Verify Data Collection is stopped.</li> <li>2. Close all windows.</li> <li>3. Select <b>Exit</b> from the File menu.</li> </ol>
Arrange windows or icons on the screen	Select <b>Cascade</b> , <b>Tile</b> , or <b>Arrange Icons</b> command from the Window menu.
Select a communications port for the Nephelometer	<ol style="list-style-type: none"> <li>1. Verify Data Collection is OFF.</li> <li>2. Select <b>Com Port</b> from the Config menu.</li> <li>3. Select a communications port and press <b>OK</b>.</li> </ol>
Set up the Nephelometer (or make changes to its set up)	<ol style="list-style-type: none"> <li>1. Select <b>Nephelometer</b> from the Config menu.</li> <li>2. Make changes to the setup and Exit the screen.</li> </ol>

## Clean Air Test

Before you begin using the Nephelometer, you should conduct a clean air test to make certain the instrument is calibrated properly. Proceed with this test step only after you have verified the program is communicating with the Nephelometer and it is able to collect data.

1. Select **Pause** from the main menu to temporarily stop data collection.
2. Select the **Nephelometer** option from the Config menu. The Set Nephelometer Parameters display appears (Figure 5-11).

The screenshot shows the 'Set Nephelometer Parameters' interface. It includes several configuration sections: 'Analog Output Range' with fields for Type (Bar. Pressure [mbar]), Full Scale [10V] Output (1,200), Zero Volts Output (0), and Offset Value; 'Analog Output Channel' with Channel No (0) and Analog Data Type (Relative Humidity); 'BackScatter Mode' with radio buttons for Total Scatter Only and Total/Backscatter; 'AutoZero Mode' with radio buttons for Manual Mode, Normal Mode, and Air-Chop Mode; 'Valve Position' with radio buttons for Normal and Zero; 'Data Delimiter' with radio buttons for Comma, Space, and Tab; 'Date/Time' with Date (11-03-1994) and Time (11:32:02) fields; and 'Time' with Averaging (5), Blank Valve (30), Autozero Period (3,600), and Autozero Measurement (300) fields. At the bottom, there are status indicators for Power, Lamp, Fan, Heater, Blower (0%), and External Host Onv, along with an 'EXIT' button.

**Figure 5-11**  
The Set Nephelometer Parameters Display Shows the Current Parameters of the Nephelometer

3. Change the Valve Position parameter to Zero (in the middle of the display) and note the Averaging Time (just to the right). The Averaging Time is the interval between data collection points.

4. Exit from the screen (use the Exit door in the lower right).
5. Resume data collection.

Clean air is now entering the Nephelometer and the Scattering Coefficient numbers, refer to Figure 5-6, should read in the E-7 or E-8 range. This is the normal background noise level of the Nephelometer. The numbers should fluctuate about zero (from negative to positive).

If the noise level is within this range, the Nephelometer is properly zeroed and calibrated and ready for use.

If the numbers are greater than E-7 or E-8 (an occasional E-6 is okay), the Nephelometer is not calibrated properly. (The Nephelometer is calibrated before shipment, but the software settings may have been corrupted or lost.)

Refer to the calibration procedure described later in this chapter to use the Windows-based calibration procedure or refer to Chapter 4 to use the DOS-based calibration procedure. The procedures themselves are identical.

---

## Menus and Commands: A Reference Guide

This section describes the menus and commands of the Windows-based Nephelometer software. The section is organized like the main menu—commands available under each menu follow the menu name. The seven menu items are:

- File
- Config
- View
- Run
- Window
- Pause/Resume
- Help

## Using the File Menu

The File menu commands let you perform functions related to software file maintenance.

### Open Log File Command

Choose the **Open Log File** command from the File menu to select the filename of an existing log file to open and display its data (Figure 5-12). You can type the filename (including drive and directory) in the filename box or select a file using the list boxes. Log files have a .DAT extension.

If you want the AutoZero Data stored in the file (if available) to be displayed along with the sample data, select (mark) the Include AutoZero Data box.

Select the **OK** button when done and watch the bottom of the display for confirmation that the file is being processed. Large files can take a long time to open. After the file is opened, the file's data appears in the Log Data Table window.


Refer to the description of the Log Data Table command under the View menu, see below, for the options that are available for viewing, graphing, copying and printing this data.

### Log Data


Choose the Log Data command from the File menu to identify the file you want the Nephelometer to log data to (Figure 5-13). The filename defaults to the last name you specified (file extension .DAT).

You can type the filename (including drive and directory) in the filename box or select a file using the list boxes. The file can be a new file or an existing file.

**Note:** *Before data is logged to a file, you MUST have selected a filename using the Log Data command. You can begin data collection (select Data Collection from the Run menu), but unless you specify a filename here, data is not saved to any file.*



**Figure 5-12**  
Select a Filename After Choosing the Open Log File Command



**Figure 5-13**  
Select a Filename After Choosing the Log Data Command

After you select the file and press **OK**, observe the bottom right corner of the display. It should indicate the name of the currently selected log file.


To stop logging data to a file, select the **Log Data** command again and verify that the logging status is OFF.

You can begin data collection with or without selecting log data file.

## Print

Choose the **Print** command from the File menu to print graphs of the data currently being logged (the file selected with the Log Data command).

You can print a scatter data graph or a status graph and, if desired, add a message header to describing the graph (Figure 5-14). The Date/Time Stamp box lets you automatically imprint the time and date the graph was printed. Graphs are printed in landscape mode.




**Figure 5-14**  
Select the Print Command from the File Menu

Before you print a graph, be certain you Pause data collection. (Select **Pause** from the main menu.) Once you send the graph to the printer, select **Resume** to continue collecting data.

Select the **Setup** button from the Print display to set up the printer.

**Note:** Normally, you will not need to use Setup or you will only need to perform Setup once, since any settings you make will be stored with the program. Refer to your Windows documentation to change printer options and setup parameters.

Figure 5-15 shows a typical scatter data graph (top) and status graph (bottom).



**Figure 5-15**  
Types of Graphs. Scatter Data Graph (top) and Status Graph (bottom)

## Print Screen

Press the **<Print Screen>** key to capture the information currently on your display and then select the **Print Screen** command from the File menu or press **<F9>** to print it.

When you select the command, the screen information is captured and sent to the printer. Only the information visible on your display monitor is printed.

## Exit

Choose the **Exit** command to exit (close) the program and return to Windows.

## Using the Config Menu

The Config menu commands let you configure the Nephelometer and select the communications port that the Nephelometer is connected to.

## Nephelometer

Choose the **Nephelometer** command from the Config menu to read the current parameters of the Nephelometer or make changes to the Nephelometer settings (Figure 5-16).

**Note:** *Setup is not the same as calibration. Select the Calibration command from the Run menu to calibrate or recalibrate the Nephelometer.*

Table 5-2 lists the parameters that can be set or modified from this display. For more detailed information about these settings refer to Chapter 6, "Serial Data Commands."

Changes take affect immediately. Use the Exit button when changes are completed.






Figure 5-16  
Choose Nephelometer from the Config Menu

Table 5-2  
Nephelometer Setup Parameters

Parameter	Equivalent Serial Data Command	Description
<b>Analog Output Range</b>	<i>SACm,v,b.bbbe-b</i>	This group of parameters is used to set the ranges for PROGRAMMABLE I/O (BNC). Ranges can be set for relative humidity, barometric pressure, sample and inlet temperatures, and scattering values.  Select the <b>Set</b> button after making changes to a specific type of port.
Type		Select <b>Barometric Pressure, Relative Humidity, Sample &amp; Inlet Temperature</b> or <b>Scattering</b> .
Full Scale (10V) Output		Enter a value within these ranges: Barometric Pressure (1-1200 Mbar) Relative Humidity (1-100%) Sample & Inlet Temperature (1-400K) Scattering (1-10)

(continued)

**Table 5-2**  
Nephelometer Setup Parameters *(continued)*

<b>Parameter</b>	<b>Equivalent Serial Data Command</b>	<b>Description</b>
Zero Volts Output		Enter a value within these ranges: Barometric Pressure (0-1199 Mbar) Relative Humidity (0-99%) Sample & Inlet Temperature (0-399K) Scattering (4-10)
Offset Value		Relates only to Scattering. The offset value added to the scattering data allows negative values (0.000e-11 to 9.999e-3).
<b>Analog Output Channel</b>	SBC,t	This group of parameters is used to map (assign) an analog data type to a hardware output on the PROGRAMMABLE I/O (BNC) connector.  Select the <b>Set</b> button after assigning a data type to each channel.
Channel No.		0 = PROGRAMMABLE I/O port
Analog Data Type		Select one of the following data types for each channel (the default channel number is shown in parenthesis): External Host Value (not assigned) Barometric Pressure (1) Sample Temperature (2) Blue $\sigma_{sp}$ (3) Green $\sigma_{sp}$ (4) Red $\sigma_{sp}$ (5) Blue $\sigma_{bsp}$ (6) Green $\sigma_{bsp}$ (7) Red $\sigma_{bsp}$ (8) Relative Humidity (0) Inlet Temperature (not assigned)
<b>Back Scatter Mode</b>	SMBz	Select <b>Total Scatter Only</b> to measure only the total scatter.  Select Total/Backscatter to enable the backscatter shutter, measuring total scatter <b>and</b> backscatter.


*(continued)*

**Table 5-2**  
Nephelometer Setup Parameters (*continued*)

<b>Parameter</b>	<b>Equivalent Serial Data Command</b>	<b>Description</b>
<b>AutoZero Mode</b>	SMZ <i>v</i>	Sets the mode used to zero baseline drift using filtered air.  In Manual Mode, zeroing of the baseline only occurs when a Z command is given.  In Normal Mode, an autozero is performed at the Autozero Period [0 to 9999]  In Air-Chop Mode, [2 to 24]. Similar to normal mode except the last number of autozero baselines are averaged together to provide a reading.
<b>Valve Position</b>	V <i>c</i>	Normal is used for normal measurements.  Zero is used for zeroing of the baseline.
<b>Data Delimiter</b>	SD <i>x</i>	Sets the character used to delimit data in the output file. Select comma, space, or tab
<b>Date/Time</b>	STT	Sets the date and time of the Nephelometer's internal clock.  Enter date as mm-dd-yyyy  Enter time as hh:mm:ss
<b>Time</b>		Set the time intervals for various parameters/operations.
Averaging	STA <i>ttt</i>	Sets the averaging time, in seconds, over which the sample is measured.
Blank Valve	STB <i>ttt</i>	Sets the blanking time, in seconds, when the valve is switching.
Autozero Period	STP <i>ttt</i>	Sets the time, in seconds, between autozeros.
Autozero Measurement	STZ <i>ttt</i>	Sets the time, in seconds, the Nephelometer spends measuring filtered air during a zero baseline measurement.
<b>Power</b>	L	Turns the Nephelometer on or off.
<b>Lamp</b>	SP	Turns the Nephelometer lamp on or off and allows you to adjust lamp power.
<b>Fan</b>	F	Turns the fan on or off.
<b>Heater</b>	H	Turns the heater on or off.
<b>Blower (xx%)</b>	B	Allows you to adjust blower power from 0 (off) to 255 (full power)
<b>External Host x<i>mv</i></b>	SX	Sets external host analog value to a specified voltage (0 to 5000).

## Com Port

Choose the **Com Port** command from the Config menu to select the communications port that the Nephelometer is attached to (Figure 5-17). The default is Com Port 2.



**Figure 5-17**

Use the COM Port Setup Display to Select the Appropriate Communications Port

If you need to change the communications port, select another communications port and press **OK**.

## Using the View Menu

The View menu commands let you display information collected by the Nephelometer and check Nephelometer configuration and status.

**Note:** *The Data, Graph, Photon Frequency, and Neph State commands can be selected only when data collection is ongoing.*

## Data

Select the **Data** command from the View menu to display data as it is being collected by the Nephelometer (Figure 5-18).

**Note:** *When you start data collection, the Data Display window is automatically opened and minimized. You can double-click on the Data Display icon to open its window instead of selecting it from the View menu.*



**Table 5-3**  
Data Display Parameters (continued)


<b>Parameter</b>	<b>Description</b>
<b>Status Record</b>	The information displayed here are current parameters calculated by the Nephelometer.
Pressure (MB)	The current barometric pressure in millibars (mb)
Sample Temp. (K)	The current sample temperature in °K
Inlet Temp (K)	The current temperature at the Nephelometer inlet in °K
Rel. Humidity	The relative humidity of the sample in %
Lamp Voltage	The current voltage applied to the lamp in volts DC
Lamp Current	The current amperage applied to the lamp in amperes
BNC Input Voltage	The current input voltage to the BNC

## Graph

Select the **Graph** command from the View menu to display Nephelometer data graphically as it is being collected (Figure 5-19).


**Note:** *When you start data collection, the Graph window is automatically opened and minimized. You can double-click on the Graph Display icon to open its window instead of selecting it from the View menu.*

The default graph options show the scattering coefficients over time in the top graph and status parameters in the lower graph.



**Figure 5-19**  
Graph Data of the Current Log Data File

Select **Graph Options** from the Graph Display to customize the graph (Figure 5-20). The parameters you can select are described in Table 5-4. To print the graph, select **Print** from the File menu.



**Figure 5-20**  
The Graph Options Window Lets You Customize Your Graphs

**Table 5-4**  
Graph Options for the Graph Display

Parameter	Description
<b>Type</b> (scattering coefficients only)	Select linear or log(rithmic) to set the x-axis scale of the graph. The maximum value is automatically determined by the maximum value of the data displayed.
<b>Scale</b>	Select AutoScale or User-defined to set the y-axis of the graph. If you select User-defined, you must also set the following parameters: Y-Axis Min.—The minimum value for the y-axis. Y-Axis Max.—The maximum value for the y-axis. Y-Axis Ticks—The value between each mark (tick) on the y-axis.
<b>Scattering Coefficients</b>	Indicates the type of data to be graphed. The default is all data. Select (deselect) each checkbox.

(continued)




**Table 5-4**  
Graph Options for the Graph Display (*continued*)

Parameter	Description
<b>Maximum No. of Data Points</b>	The maximum number of points to include on the graph. If the Log Data file includes more than this number of data points, only the most current points are displayed.
<b>Offset</b>	Offset the y-axis from 0 (zero) by this amount.
<b>Refresh Rate</b>	How often (in percentage) the graph display is updated (refreshed). This is independent of the averaging data time.
<b>Title</b>	Enter a title for the graph. Maximum 60 characters.
<b>X-Axis Text</b>	Enter a label for the x-axis.
<b>Y-Axis Text</b>	Enter a label for the y-axis.
<b>Legend Text</b>	Check this box to include the legend on the graph. The legend identifies the type of data displayed by each line.

## Photon Frequency

Choose the **Photon Frequency** command from the View menu to display the photon frequency data (in Hz) as it is collected from the Nephelometer (Figure 5-21) and calculated by the software.

**Note:** When you start data collection, the Photon Frequency window is automatically opened and minimized. You can double-click on the Photon Frequency icon to open its window instead of selecting it from the View menu.



Photon Frequency in Hz						
	Total Scatter			BackScatter		
	Cal	Meas	Dark	Cal	Meas	Dark
Blue	1.58419E+5	1.05855E+3	1.15270E+1	9.93232E+4	5.29145E+2	1.16071E+1
Green	2.50964E+5	8.39585E+2	1.90516E+1	1.58559E+5	4.23738E+2	1.70104E+1
Red	1.92582E+5	6.04036E+2	2.57637E+2	1.20173E+5	4.80481E+2	2.58597E+2


**Figure 5-21**  
Photon Frequency Display

The Photon Frequency display shows the raw photon frequency calculated by the software from the raw photon counts collected by

the Nephelometer. The raw photon frequency is displayed for each section of the calibrate shutter (calibrate, measure, dark), at each wavelength (blue, green, red). The data is updated approximately once a second.

**Note:** For more details about photon frequency and calibrate shutter, see “Signal Processing” in Chapter 7.

Select **Graph** from this display to graphically view the data that has been collected in the current Log Data file (Figure 5-22).




**Figure 5-22**  
The Photon Frequency Graph Displays Current Measurements

Select **Graph Options** from the Photon Frequency Graph display to customize this graph. The options are nearly identical to those described in Table 5-4 except you must select a color (wavelength) to indicate the data you want shown on the graph.

## Neph State

Choose the **Neph State** command from the View menu to display the operational status of selected Nephelometer components (Figure 5-23). Use this display to verify general Nephelometer status and to troubleshoot problems with the Nephelometer.

**Note:** When you start data collection, the Neph State window is automatically opened and minimized. You can double-click on the Neph State icon to open its window instead of selecting it from the View menu.



**Figure 5-23**

The Status of the Nephelometer is Displayed by Selecting Neph State from the View Menu

All items on the display should be marked with a check mark (✓). If they are marked with an X, the component/parameter is reporting a hardware fault. Review Chapters 2 and 3 and refer to Chapter 8 for troubleshooting procedures to correct the problem.


## Neph Config Data

Choose the **Neph Config Data** command from the View menu to display Nephelometer configuration parameters (Figure 5-24). Current Nephelometer configuration parameters are stored in a text file named NEPHCNFG.DAT, which is the default config file used by the software.

You can display current configuration parameters or read a configuration file from disk when the **Use Default Config File** option is selected.

The information in this display is useful for analyzing system problems because you can view all of the parameters as a group. It is also useful because you can store, then easily retrieve, different configurations if the current Nephelometer configuration should become lost or corrupted.

**Note:** You cannot make changes to the Nephelometer Configuration Data from this display.



**Figure 5-24**  
Nephelometer Configuration Data is Displayed by Selecting Neph Config Data from the View Menu

The buttons across the bottom of the window let you select the source and destination of the displayed configuration data. If the **Use Default Config. File** is marked (X), the information will be written to or read from the default configuration file named NEPHCNFG.DAT. **Do not remove this file from your hard drive or change its name.**

**Table 5-5**  
Nephelometer Configuration Data

Read Data from Disk	Select this button to read Nephelometer Configuration Data from a file other than the default configuration file. You can then review this data before downloading it to the Nephelometer.  Enter the complete path and filename of the file to be read.
Write Data to Disk	Select this button to write the currently displayed Nephelometer Configuration Data to a file to save it.  <b>Caution:</b> You can write over the default configuration file! If the "Use Default Config. File" box is checked and you select Write Data to Disk, the currently displayed configuration information will overwrite the default configuration file. The old default configuration file will then become a backup configuration file (NEPHCNFG.BAK).
Read Data from Neph	Select this button to read the configuration data currently being used by the Nephelometer. Use this button to verify the operating configuration parameters.
Write Data to Neph	Select this button to download the displayed configuration data to the Nephelometer, replacing the configuration parameters that are stored in the Nephelometer.
Send Data to Printer	Select this button to print the currently displayed Nephelometer Configuration Data.

## Log Data Table

Choose the Log Data Table command from the View menu to display data from a previously saved log file or the data currently being collected (Figure 5-25).

**Note:** To view a previously saved data file, you must open the log file first by choosing the Open Log File command from the File menu. Large files can take a long time to process and load. Watch the bottom of the screen to verify the file is being processed.

Scat Coeff. Data Table					
Graph Data - tso0929.log					
Row	Date/Time	TotScat-Blue	TotScat-Green	TotScat-Red	BkScat-Blue
1	09/29/94 9:18:40	-4.41400E-7	6.71700E-8	1.12600E-7	02700E-7
2	09/29/94 9:19:43	-3.23400E-7	-1.67100E-7	1.80000E-7	82600E-7
3	09/29/94 9:20:44	-2.97700E-7	2.41400E-7	9.61200E-8	47800E-7
4	09/29/94 9:21:45	6.48200E-8	-4.27000E-9	-2.54100E-8	98800E-8
5	09/29/94 9:22:46	-1.77500E-7	-5.28900E-8	6.37800E-8	76700E-7
6	09/29/94 9:23:47	4.35600E-8	1.17400E-7	2.34000E-7	59900E-7
7	09/29/94 9:24:48	9.38400E-8	-1.74900E-7	2.19700E-7	11000E-8
8	09/29/94 9:25:49	-9.15400E-8	-1.04000E-7	1.90100E-7	96200E-7
9	09/29/94 9:26:52	-1.58300E-7	-4.09000E-9	1.43500E-7	17400E-8
10	09/29/94 9:27:53	1.12000E-7	-4.87900E-8	2.23000E-7	21700E-8
11	09/29/94 9:28:54	-1.89500E-7	-2.90200E-7	-2.95100E-8	03700E-8
12	09/29/94 9:29:56	-2.91200E-7	-9.40900E-8	3.47800E-7	26000E-7
13	09/29/94 9:30:58	-6.59800E-7	1.88400E-7	2.68200E-7	30200E-8
14	09/29/94 9:31:59	2.20000E-7	-6.22200E-8	-1.25400E-7	28000E-9
15	09/29/94 9:33:01	-1.69900E-7	-2.20900E-7	-8.03000E-9	18200E-7
16	09/29/94 9:34:01	-3.97100E-7	-3.72800E-7	2.10800E-7	11300E-7
17	09/29/94 9:35:03	-1.16300E-7	-1.01200E-7	2.60400E-7	63000E-7
18	09/29/94 9:36:04	-2.26300E-7	-2.65700E-7	1.14600E-8	19900E-7

Graph Copy/Print Close Start Row: 1 End Row: 809 Total Rows: 809

**Figure 5-25**

Display Log File Records in Tabular Format by Selecting Log Data Table from the View Menu

Up to eighteen (18) data records are displayed in the window at one time. Use the <↑> or <↓> and <→> or <←> arrow keys to view all the data and all rows of data in the file. The bottom right corner of the window indicates the total number of rows of data.

The information displayed in the table depends on whether you are in the analysis operating mode or polled operating mode.

In the analysis mode you are viewing a file that has been saved (not the current log file), and the following information is displayed:

- Total Backscattering Coefficients for each wavelength
- Backscattering Coefficients for each wavelength.
- Green Sensitivity
- Barometric Pressure
- Sample Temperature
- Inlet Temperature
- Relative Humidity

- Lamp Voltage
- Lamp Current
- BNC Input Voltage
- Status Flags

In the polled data mode you are viewing the current log file and the following information is displayed:

- Total Backscattering Coefficients for each wavelength
- Backscattering Coefficients for each wavelength.
- Sample Temperature
- Barometric Pressure
- BNC Input Voltage
- Relative Humidity

The buttons and indicators at the bottom of the table allow you to graph the displayed data, copy the data to the clipboard or another file, or print the data.


You can graph, copy, or print all rows of data in the table or select a range of rows. The Start Row and End Row boxes indicate the row of data that will be graphed, copied, or printed.

Use the following procedures to identify the rows to include:


- Point the cursor to the first row of data and highlight it by clicking the mouse button. Continue holding the mouse button and drag the cursor to the last record you want included.
- Use the <↑> and <↓> arrows on the Start Row box and End Row box to display the first and last record you want graphed.

Once you have defined the rows of data, select **Graph** or **Copy/Print**.


If you select **Graph**, the Log Data Graph Options display appears (Figure 5-26). Select the appropriate type of graph and press **OK**. Figure 5-27 shows the type of information shown in each graph.



**Figure 5-26**  
Select the Type of Graph from the Log Data Graph Options




**Figure 5-27a**  
Operating Status Graph




**Figure 5-27b**  
Photon Counts (Calibrator) Graph






**Figure 5-27c**  
Photon Counts (Dark) Graph




**Figure 5-27d**  
Photon Counts (Measure) Graph



**Figure 5-27e**  
Scattering Coefficients Graph

If you select **Copy/Print**, the Destination dialog box appears (Figure 5-28).

The Destination dialog box lets you select the type of data to be copied (only the record types marked will be copied), the delimiter that will be used between fields, and the destination of the data (clipboard or disk file). If you mark the “Include Header Message,” the header information—Row, Date/Time, TotScat...etc.—is copied along with the data.



**Figure 5-28**  
The Destination Dialog Box

If you select **Disk File** from the Destination dialog box, a Save Disk File dialog box appears.

If you select **Clipboard**, the data is copied to the Clipboard and can be pasted into another Windows application such as Microsoft Excel.®

If you select **Printer**, the data is sent to the printer.

## Using the Run Menu

The Run menu commands let you calibrate the Nephelometer, start data collection, and control the Nephelometer through a terminal (using the serial data commands described in Chapter 6).

### Calibration

Choose the Calibration command from the Run menu to calibrate (or recalibrate) the Nephelometer (Figure 5-29).

		Current (B)	New (B)	Current (G)	New (G)	Current (R)	New (R)
K2		3.49900E-3		3.26700E-3		3.13200E-3	
K4		5.06000E-1		4.62000E-1		5.87000E-1	

		TotS-Cal	TotS-Meas	TotS-Dark	BkS-Cal	BkS-Meas	BkS-Dark
Low Span Gas	Blue						
	Green						
	Red						
High Span Gas	Blue						
	Green						
	Red						

1 <input type="checkbox"/> Set Low Span Gas	Air	±	8 <input type="checkbox"/> Plug inlet and outlet (manual)	<input type="button" value="Graph Meas Data"/> <input type="button" value="Pause"/> <input type="button" value="New Avg Time"/> <input type="button" value="Setup"/> <input type="button" value="Exit"/>
2 <input type="checkbox"/> Set High Span Gas	CO2	±	9 <input type="checkbox"/> Feed High Span Gas into line	
3 <input type="checkbox"/> Switch Valve Filter	(NORMAL)		10 <input type="checkbox"/> Start Data Collection-High Gas	
4 <input type="checkbox"/> Start Data Collection-Low Gas			11 <input type="checkbox"/> Clear Buffer	
5 <input type="checkbox"/> Clear Buffer			12 <input type="checkbox"/> Store High Span Gas Data	
6 <input type="checkbox"/> Store Low Span Gas Data			13 <input type="checkbox"/> Write Calibration Label	
7 <input type="checkbox"/> Toggle Blower Off	(ON)		14 <input type="checkbox"/> Save Cal Data to Neph	
			15 <input type="checkbox"/> Zero w/ Clean Air	

Data Logging OFF	Calibration
---------------------	-------------

Figure 5-29  
The Nephelometer Calibration Data Display

The information displayed on the Calibration Data screen is the same as that described in Chapter 4, “Performing Calibration,” and you use the same procedure for calibration.

- The top section of the window shows the current calibration settings stored in the Nephelometer and any new calibration settings not yet saved to the Nephelometer.
- The middle section of the window shows calibration data in frequency (Hz) for low span and high span gas in each wavelength (blue, green, and red). Calibration data is displayed for both total scatter and backscatter modes.
- The bottom section of the window lists the steps you will use to calibrate the Nephelometer. Perform each step in the order given. As you complete each step, select the button for the appropriate step. A check (✓) is displayed next to the button as it is selected or executed.

Calibration allows you to easily calibrate the Nephelometer using two span gases and compare the results between air (low span) and CO<sub>2</sub> (high span).

You should calibrate the Nephelometer before an intensive experiment, calibrate periodically to verify no drift has occurred, calibrate if the reference chopper is dirty or scratched, or if you clean the chopper as part of periodic maintenance.


The Graph Meas Data button lets you observe (graphically) the low span or high span data as calibration is under way. This allows you to see graphically the low span gas purged and replaced by the high span gas.

The Pause button lets you stop data collection during calibration. The button then displays “Continue.” Select **Continue** to restart the data collection.

The New Avg Time button is available on this display to let you set/reset a new averaging time without going to the Nephelometer Setup display. This allows you to select a lower averaging time to have a better response time in viewing the low span gas being purged and replaced by the high span gas.

**Note:** If the averaging time is changed to a low number (such as 5 sec) to view purging, be sure to change the averaging time back to a higher number (default is 300 sec). Otherwise, a low averaging time may generate very noisy data and result in bad calibration.

Before calibration, use the Setup button to set up selected calibration parameters (Figure 5-30). From the Calibration Parameter display you can set up Calibration Points, Calibration Constants and PMT Voltages you want to use and enter a Calibration Label String.



**Figure 5-30**  
Calibration Parameters Setup Screen

Table 5-6 describes each of the parameters you can set. After entering values for each device/parameter type, you must press the Set button for the new values to be stored.

**Table 5-6**  
Calibration Parameters

<b>Parameter</b>	<b>Description</b>
<b>Calibration Points</b>	These parameters let you calibrate the four internal sensors by entering low and high calibration points.
Sensor Type	Select <b>Barometric Pressure, Inlet Temperature, Relative Humidity, or Sample Temperature.</b>
r1	Sets the low value of raw A/D converter bits.
v1	Sets the data value (actual value × 10) for low scale calibration point.
r2	Sets the high value of raw A/D converter bits.
v2	Sets the data value (actual value × 10) for high scale calibration point.
<b>Calibration Constants</b>	These parameters let you enter calibration constants for each of the detector wavelengths.
Color Type	Select <b>Blue, Green, or Red</b>
K1	Sets the photon count input pulse width (dead time).
K2	Sets the total scatter calibration of the reference chopper in units of inverse meters.
K3	Sets the Rayleigh scatter of air at 273 °K and 1013.25 mb.
k4	Sets the backscatter calibration of the reference chopper in units of inverse meters.
<b>PMT Voltage</b>	These parameters set the voltage that drives each photomultiplier tube.
Color Type	Select <b>Blue, Green, or Red.</b>
Voltage	Sets the voltage level applied to the photomultiplier (0-1200).
<b>Calibration Label String</b>	Enter text to describe the calibration parameters, for example: Last Calibration 12/25/94. Up to 80 characters.

Select the **Exit** button when you are finished with the Calibration Parameters display.

## Hardware Setup

Set up the gas before you begin calibration. Gas should be approximately room temperature as it reaches the Nephelometer. One way to ensure this is to use a 6- to 9-meter (20- to 30-foot) length of 6-mm (¼-in.) diameter tubing from the CO<sub>2</sub> tank to the Nephelometer. (The length of the tubing affects the temperature of the CO<sub>2</sub>.)



### Caution

Use calibration gases only in a well-ventilated area or exhaust the gases outside. Many gases used for calibration can cause asphyxiation if used in a confined area.

**Note:** *If you are using air as the low span gas and CO<sub>2</sub> as the high span gas, perform the calibration steps in the order given. If you are using a gas other than air as the low span gas (like helium), refer to calibration process given at the end of this procedure.*

**Table 5-7**  
Calibration Procedure

Step	Explanation
1. Set Low Span Gas	Select the low span gas. Air is the default. The following types of gas can be selected from the list menu: air, CO <sub>2</sub> , CCl <sub>2</sub> F <sub>2</sub> , SF <sub>6</sub> , R22, He and Ar.
2. Set High Span Gas	Select a high span gas, or a gas with a higher scattering coefficient than the low span gas, CO <sub>2</sub> is the default. The following types of gas can be selected from the list menu: air, CO <sub>2</sub> , CCl <sub>2</sub> F <sub>2</sub> , SF <sub>6</sub> , R22, He and Ar.
3. Switch Valve Filter	Select to switch the filtered air valve to the zero (closed) position and to switch the large white HEPA filter in-line with the air inlet. When you select step 3, a blanking time begins. The blanking time is the time it takes the filtered air valve to move to the zero position plus the time it takes to purge the Nephelometer with filtered air. The blanking time is set using the Configure Nephelometer Parameters window or the STB command (see Chapter 6).

(continued)

**Table 5-7**  
Calibration Procedure (*continued*)

<b>Step</b>	<b>Explanation</b>
<b>4.</b> Start Data Collection-Low Gas	<p>Select to begin low span gas data collection. After you select step 4, data collection begins after the blanking time described in step 3.</p> <p>The low span gas area of the calibration screen is updated, but not the high span area.</p> <p>If you want to observe the filtered air purge, select New Avg Time to select a smaller averaging time of 5 seconds. The number in the Meas field should fluctuate, not rise and fall. When the Meas field is stable, select New Ave Time to select a 300-second averaging time and proceed to step 5.</p>
<b>5.</b> Clear Buffer	<p>Select to clear the Nephelometer internal buffers of the data collected during step 4. Wait at least 300 seconds, then proceed to step 6.</p>
<b>6.</b> Store Low Span Gas Data	<p>Select to store in the software, the 300-seconds of low span gas data collected in step 5. Selecting step 6 causes the low span gas area of the calibration screen to freeze and low span data is no longer updated.</p>
<b>7.</b> Toggle Blower To Off (OFF)	<p>Select to turn off the blower. High span gas should be regulated at the source.</p> <p><b>Note:</b> <i>Step 7 can also be used to toggle on the blower.</i></p>
<b>8.</b> Plug inlet & outlet (manual)	<p>Remove the top cover of the Nephelometer using the procedure in Chapter 8. Plug the aerosol inlet and outlet of the Nephelometer (Figure 4-19).</p>
<b>9.</b> Feed High Span Gas into line (manual)	<p>Remove the tubing from the inlet of the blue DQ filter and connect the tubing from the high span gas to either the inlet tubing or to the inlet block (Figure 4-19). Turn on the gas until the flowrate is approximately 5 lpm. You should feel gas exiting through the DQ filter in the opposite direction of the arrow on the filter. Proceed to step 10.</p>

*continued*



**Table 5-7**  
Calibration Procedure (*continued*)

<b>Step</b>	<b>Explanation</b>
<b>10. Start Data Collection-High Gas</b>	<p>Select to begin high span gas data collection. The high span gas area of the calibration screen is updated, but not the low span area.</p> <p>If you want to observe the high span gas purge, select New Avg Time to select smaller averaging time of 5 seconds. The number in the Meas field should fluctuate, not rise and fall. When the Meas field is stable with minimal fluctuation, select New Avg Time to select a 300-second averaging time and proceed to step 11.</p>
<b>11. Clear Buffer</b>	<p>Select to clear the Nephelometer internal buffers of the data collected during step 10. Remain at this step for at least 300 seconds, then proceed to step 12.</p>
<b>12. Store High Span Gas Data</b>	<p>Select to store to the software the 300 seconds of high span gas data collected in step 11. Selecting step 12 causes the high span gas area of the calibration screen to freeze and high span data is no longer updated.</p> <p><b>Note:</b> <i>At this time, new calibration constants (K2 and K4) are calculated and the top section of the screen is updated. If the new constants are not satisfactory, you can recalibrate without changing the constants, as long as you do not perform step 14.</i></p>
<b>13. Write Calibration Label</b>	<p>Select to write an identifying label for the calibration performed. The label can include the serial number of the Nephelometer, the calibration date, the low span and high span gases used, and the initials of the person performing the calibration.</p>

(*continued*)

**Table 5-7**  
Calibration Procedure (*continued*)

<b>Step</b>	<b>Explanation</b>
<p><b>14.</b> Save Cal Data (K2 &amp; K4) to Neph</p>	<p>If the new K2 and K4 calibration constants are satisfactory, select 14 to write the new K2 and K4 constants to the Nephelometer.</p> <p>After selecting step 14,</p> <ol style="list-style-type: none"> <li><b>1.</b> Remove the high span gas and the inlet and outlet plugs from the Nephelometer.</li> <li><b>2.</b> Reconnect the DQ filter to the inlet block, with the arrow of the filter pointing away from the inlet block.</li> <li><b>3.</b> Replace the top cover of the Nephelometer, attaching the cover with the four screws.</li> <li><b>4.</b> Select 7 to toggle on the blower and to purge the Nephelometer of high span gas.</li> </ol> <p><b>Note:</b> <i>To allow the Nephelometer to be completely purged of high span gas, wait at least 600 seconds with the blower on before proceeding to step 15.</i></p>
<p><b>15.</b> Zero w/Clean Air</p>	<p><i>Optional</i>—Select to purge the Nephelometer with clean air and to perform a zero background measurement.</p>

(*continued*)

**Table 5-7**  
Calibration Procedure (*continued*)

<b>Step</b>	<b>Explanation</b>
<p>If you are using a low span gas other than air, perform the calibration steps in this order:</p>	<ol style="list-style-type: none"> <li><b>1.</b> See step 1.</li> <li><b>2.</b> See step 2.</li> <li><b>3.</b> See step 3.</li> <li><b>7.</b> See step 7.</li> <li><b>8.</b> See step 8.</li> <li><b>9.</b> Remove the tubing from the inlet of the blue DQ filter and connect the tubing from the low span gas to either the inlet tubing or to the inlet block (Figure 4-19). Turn on the gas until the flowrate is approximately 5 lpm. You should feel gas exiting through the DQ filter in the opposite direction of the arrow on the filter.</li> <li><b>4.</b> See step 4.</li> <li><b>5.</b> See step 5.</li> <li><b>6.</b> See step 6.</li> <li><b>9.</b> See step 9.</li> <li><b>A.</b> See step 10.</li> <li><b>B.</b> See step 11.</li> <li><b>C.</b> See step 12.</li> <li><b>D.</b> See step 13.</li> <li><b>E.</b> See step 14.</li> <li><b>F.</b> See step 15.</li> </ol>

Select **Exit** when calibration is complete.


### **Data Collection**

Choose the **Data Collection** command from the Run menu to begin collecting data from the Nephelometer.

**Note:** *You can collect data without logging it to a file. If you start Data Collection without specifying a file in which to log data, no data is logged. To log data, choose Log Data from the File menu and enter the name of the file in which you want to collect data.*

When you select **Data Collection**, four windows are opened and minimized (Figure 5-31). Each window represents one of the following commands of the View menu: Data, Graph, Photon Frequency and Neph State. Refer to the description of each of these commands, above. You can maximize (restore) these icons to view the data being collected or check Nephelometer status.

To verify that data collection is underway, look in the lower right corner of the display. It should indicate Polled.



**Figure 5-31**  
Selecting Data Collection from the Run Menu

To stop data collection, select **Data Collection** from the Run menu again.


Data files can become quite large if you have a short averaging time and allow data to accumulate over a long period of time. As a rule or thumb, do not allow files to get much larger than 450 kbytes.

## Terminal Mode

Choose the **Terminal Mode** command from the Run menu to issue serial data commands directly to the Nephelometer (Figure 5-32).

You may want to use this method of controlling the Nephelometer during troubleshooting.

The serial data commands are described in Chapter 6. Refer to it for details.



```
Terminal Display
RU
T,1994,10,19,14,16,37
B,2129112,49109,274,2769,1334311,24621,256,2769,963,6,300,0
G,3369857,39024,357,2769,2127474,19613,357,2769,963,6,300,0
R,2582601,27932,5172,2769,1611354,22169,5051,2769,963,6,300,0
D,2800,59,2.555e-5,1.264e-5,6.282e-6,1.265e-5,6.215e-6,4.020e-6
Y,253120,963,6,300,0,296,8,29,3,12,8,5,9,18,0000
STA
300
SP
75
SK
ERROR
SKG
20000,3.923e-3,1.226e-5,.479
SKR
20000,3.557e-3,4.605e-6,.438
SKB
20000,3.940e-3,2.789e-5,.475
SL
SEPT 30 1994 5/N 1004 AIR CO2.
SMB
1
STT
1994,10,19,14,19,12
```

**Figure 5-32**  
Select Terminal Mode to Issue Serial Data Commands to the Nephelometer

To issue a command, type the command and press **<Enter>**. The Nephelometer responds with the information requested, an **OK** acknowledgment, or an **ERROR** acknowledgment.

## Using the Window Menu

The Window menu commands help you organize the open windows or arrange icons in open windows.

The commands available under this menu (Cascade, Tile, and Arrange Icons) are standard Windows commands. Please refer to your Windows documentation for complete information.

## Using the Pause/Resume Menu

The Pause/Resume menu toggles data collection on and off. After you begin data collection (select the **Data Collection** command from the Run menu), data is collected from the Nephelometer. To stop data collection temporarily without quitting data collection, select **Data Collection** from the Pause menu.

When you select to pause, the menu name changes to “Resume,” so when you want to start collecting data again, select **Data Collection** from the Resume menu. Data collection then continues and data is logged to the same log file as before.

## Using the Help Menu


The Help menu commands provide access to the program help database. This help function is patterned after Windows Help, so if you are familiar with Help in Windows, you should find this function easy to use. If you are unfamiliar with Windows Help, review your Windows documentation.

Help is context sensitive. If you select Help anywhere in the program, you will find help for the area you are currently using. For example, if you press <F1>, while you are in the Terminal screen, the Help screen that details the serial data commands will appear.

## Contents

Choose the Contents command from the Help menu to view the contents of online help (Figure 5-33). To select the relevant help topic, click on the keyword. You can select any word or group of words that are underlined (in green color).

**Note:** You can open help contents anytime by pressing the <F1> key.



**Figure 5-33**  
The Help Contents Display Provides a List of Help Topics

## About

Choose the **About** command from the Help menu to display the version level and copyright notice for this Windows-based Nephelometer program (Figure 5-34).



**Figure 5-34**  
About Shows the Current Version of Your Software and the Copyright Notice





## CHAPTER 6

# Using Serial Data Commands

This chapter contains information you need if you are writing your own software for a computer or data acquisition system. The main sections are:


- Pin connectors
- Baud rate
- Parity
- Command definitions, syntax and examples, as well as input and troubleshooting directions.

**Note:** *If you are using the software that TSI provides, refer to Chapter 4, "Using Nephelometer Software."*

---

## Pin Connectors

The Nephelometer has a single 9-pin, D-subminiature connector port on the back panel labeled COM PORT (Figure 6-1). This communication port is configured at the factory to work with RS-232 type devices. Table 6-1 gives the signal connections.



**Figure 6-1**  
COM PORT Pin Designations

**Table 6-1**  
Signal Connections for RS-232 Configurations

<b>Pin Number</b>	<b>RS-232 Signal</b>
1	—
2	Transmit Output
3	Receive Input
4	—
5	GND
6	—
7	—
8	—
9	—

---

## Baud Rate

The baud-rate setting is the rate of communication in terms of bits per second (baud). The Nephelometer uses a baud rate setting of 9600. For proper communications, make sure that all software used with the instrument is also set at this rate.

---

## Parity (8-Bits Even)

Parity is the additional bit that accompanies the seven data bits to confirm that they are transmitted correctly. It is set so that the number of “1” bits (high) in a transmitted character is always an even number. The Nephelometer uses even parity as the only setting.

---

# Commands

The Nephelometer uses an ASCII-based communications protocol that utilizes the RS-232 port of a computer to transmit commands in the form of strings.

The four types of commands are:

- Set commands, which set all the operating parameters for the Nephelometer
- Action commands, which control mechanical components of the Nephelometer
- Read (polled) commands, in which the Nephelometer sends data in response to a specific request from the computer
- Unpolled commands, in which the Nephelometer automatically outputs data records at specific intervals.

No line feed characters are transmitted. Either the requested data or an "OK" is returned if the command is understood. The word "ERROR" is returned if the command is not understood or if the command has an invalid parameter.

Table 6-2 gives a quick reference of all the serial commands. Command definitions, syntax, and examples begin after Table 6-2. Directions for inputting commands and troubleshooting commands are given at the end of this section.

**Table 6-2**  
Serial Commands

**Set Commands**

<b>SA</b>	Set <u>A</u> nalog output range <b>SA</b> <i>cm,v,b.bb</i> <i>be-b</i>
<b>SB</b>	Set the analog output channel <b>SB</b> <i>c,t</i>
<b>SC</b>	Set <u>C</u> alibration points <b>SC</b> <i>xr1,v1,r2,v2</i>
<b>SD</b>	Set <u>D</u> ata <u>D</u> elimiter <b>SD</b> <i>x</i>
<b>SK</b>	Set <u>K</u> calibration constants <b>SK</b> <i>caaaaa,b.bb</i> <i>be-b,c.cc</i> <i>ce-c,dd</i> <i>de-d</i>
<b>SL</b>	Set calibration <u>L</u> abel string
<b>SMB</b>	Set <u>M</u> ode for <u>B</u> ackscatter shutter <b>SMB</b> <i>z</i>
<b>SMZ</b>	Set <u>M</u> ode for auto <u>Z</u> ero baseline measurement <b>SMZ</b> <i>v</i>
<b>SP</b>	Set lamp <u>P</u> ower <b>SP</b> <i>www</i>
<b>STA</b>	Set <u>T</u> ime for <u>A</u> veraging <b>STA</b> <i>ttt</i>
<b>STB</b>	Set <u>T</u> ime to <u>B</u> lank valve <b>STB</b> <i>ttt</i>
<b>STP</b>	Set autozero <u>P</u> eriod <b>STP</b> <i>tttt</i>
<b>STT</b>	Set <u>T</u> ime and date <b>STT</b> <i>yyyy,mm,dd,hh,nn,ss</i>
<b>STZ</b>	Set <u>T</u> ime for auto <u>Z</u> ero measurement <b>STZ</b> <i>tttt</i>
<b>SV</b>	Set <u>V</u> oltage on the photomultiplier tube <b>SV</b> <i>cb</i> <i>bb</i>
<b>SX</b>	Set e <u>X</u> ternal host analog value <b>SX</b> <i>v</i> <i>v</i> <i>v</i>
<b>SZ</b>	Set <u>Z</u> ero baseline (background) <b>SZ</b> <i>cx.xx</i> <i>xe-x,y.yyye-y,r.rr</i> <i>re-r</i>

**Read (Polled) Commands**

<b>RA</b>	Read raw <u>A</u> nalog bit values
<b>RB</b>	Read <u>B</u> arometric pressure
<b>RD</b>	Read scatter <u>D</u> ata
<b>RF</b>	Read status <u>F</u> lags
<b>RI</b>	Read <u>I</u> nlet temperature
<b>RL</b>	Read <u>L</u> amp voltage and current
<b>RN</b>	Read <u>N</u> oise (sensitivity) level
<b>RO</b>	Read accumulated <u>O</u> n time of Nephelometer
<b>RP</b>	Read all <u>P</u> hoton counts (for blue, green, and red)
<b>RPG</b>	Read <u>P</u> hoton counts for <u>G</u> reen
<b>RR</b>	Read <u>R</u> elative Humidity
<b>RS</b>	Read <u>S</u> ample temperature
<b>RT</b>	Read <u>T</u> ime and date
<b>RU</b>	Read all enabled <u>U</u> npolled records
<b>RV</b>	Read firmware <u>V</u> ersion
<b>RX</b>	Read au <u>X</u> iliary PROGRAMMABLE I/O input port
<b>RY</b>	Read Auxiliar <u>Y</u> data record
<b>RZ</b>	Read <u>Z</u> ero background data record

**Unpolled Commands**

<b>UB</b>	<u>U</u> npolled operation <u>B</u> egins
<b>UD</b>	<u>U</u> npolled <u>D</u> ata record
<b>UE</b>	<u>U</u> npolled mode <u>E</u> nds
<b>UP</b>	<u>U</u> npolled <u>P</u> hoton counts <b>UP</b> <i>n</i>
<b>UT</b>	<u>U</u> npolled <u>T</u> ime and date
<b>UY</b>	<u>U</u> npolled Auxiliar <u>Y</u> status data record
<b>UZ</b>	<u>U</u> npolled <u>Z</u> ero background data record

**Action Commands**

<b>B</b>	<u>B</u> lower control <b>B</b> <i>xxx</i>
<b>F</b>	<u>F</u> an control <b>F</b> <i>c</i>
<b>H</b>	<u>H</u> eater control <b>H</b> <i>c</i>
<b>L</b>	<u>L</u> amp control <b>L</b> <i>c</i>
<b>PD</b>	<u>P</u> ower <u>D</u> own
<b>PU</b>	<u>P</u> ower <u>U</u> p
<b>V</b>	<u>V</u> alve position <b>V</b> <i>c</i>
<b>Z</b>	<u>Z</u> ero command

## Set Commands

Set commands allow you to set up operating parameters for the Nephelometer. If a set command is sent with no parameter, the current parameter is echoed.

**Note:** *Some of the commands directly affect or are affected by other commands. Refer to other commands where indicated.*

<p><b>SA</b></p> <p><b>see also</b></p> <p><b>SB</b></p> <p><b>SX</b></p> <p><b>RX</b></p>	<p><u>Set Analog output range</u></p> <p>SA sets the range scales of the PROGRAMMABLE I/O (BNC) port. Ranges may be set for relative humidity, barometric pressure, sample and inlet temperatures, and scattering values.</p> <p><b>S</b><i>Acm,v,b.bbbe-b</i></p> <p>where:</p> <ul style="list-style-type: none"> <li><i>c</i> = analog output <ul style="list-style-type: none"> <li>R for relative humidity</li> <li>P for barometric pressure</li> <li>T for sample and inlet temperatures</li> <li>S for scattering</li> </ul> </li> <li><i>m</i> = data value represented by zero (0) volts output <ul style="list-style-type: none"> <li>For R, the integer relative humidity in percent (0 - 99)</li> <li>For P, the integer pressure in mbar (0-1199)</li> <li>For T, the integer temperature in K (0-399)</li> <li>For S, the absolute value of the scattering coefficient exponent (integer 4-10). All scattering channels share the same range.</li> </ul> </li> <li><i>v</i> = data value represented by full scale output (BNC = 5 volts) <ul style="list-style-type: none"> <li>For R, the integer relative humidity in percent (1 - 100)</li> <li>For P, the integer pressure in mbar (1-1200)</li> <li>For T, the integer temperature in K (1-400)</li> <li>For S, the integer number of volts per decade of scattering (1-10 acceptable)</li> </ul> </li> </ul> <p><b>Note:</b> For SAR, SAP, and SAT commands, <i>v</i> must be greater than <i>m</i>.</p> <p><i>b.bbbe-b</i> = offset value added to the scattering data so that negative scattering values are allowed (0.000e-11 to 9.999e-3). This parameter is only used with the SAS command.</p>
--	---

*(continued on next page)*

(continued)

<b>SA</b>	<p><b>Examples</b></p> <p>To set the barometric pressure output to 0.0 volts at zero (0) pressure, and 5.0 volts (BNC) at 1200 mb pressure: SAP0,2400</p> <p>To set the sample temperature output to 0.0 volts at 220 K, and 5.0 volts BNC at 320 K: SAT220,320</p> <p>To set the scattering outputs range to 0.0 volts at <math>10^{-7}m^{-1}</math>, and to 5.0 volts analog at <math>10^{-2}m^{-1}</math> (offset by <math>1.0e-7m^{-1}</math>): SAS7,2,1.00e-7</p>
-----------	--

<b>SB</b>  <b>see also</b> <b>SA</b> <b>SX</b> <b>RX</b>	<p>Set the analog output channel</p> <p>SB maps an analog data type to a hardware output on the PROGRAMMABLE I/O (BNC) connector.</p> <p><b>SB<sub>c,t</sub></b> where: <i>c</i> = analog channel number 0 = PROGRAMMABLE I/O port <i>t</i> = analog data type 0 = external host analog value (set by the SX command) 1 = Barometric pressure 2 = Sample temperature 3 = Blue <math>\sigma</math>sp 4 = Green <math>\sigma</math>sp 5 = Red <math>\sigma</math>sp 6 = Blue <math>\sigma</math>bsp 7 = Green <math>\sigma</math>bsp 8 = Red <math>\sigma</math>bsp 9 = Relative humidity 10 = Inlet temperature</p>
---	--

(continued on next page)

*(continued)*

<b>SB</b>	The following is the factory default configuration:			
	<table><thead><tr><th><u>Channel</u></th><th><u>Type</u></th></tr></thead><tbody><tr><td>0</td><td>9 (Relative humidity)</td></tr></tbody></table>	<u>Channel</u>	<u>Type</u>	0
<u>Channel</u>	<u>Type</u>			
0	9 (Relative humidity)			
	<b>Examples</b>			
	To make channel 0 computer controlled: SB0,0			
	<b>Note:</b> The SX command automatically achieves the same results as SB0,0.			
	To assign Barometric Pressure to channel 0: SB0,1			
	To echo the channel 1 assignment: SB1			

*continued on next page*

<p><b>SC</b></p> <p><b>see also</b></p> <p><b>RA</b></p> <p><b>RB</b></p> <p><b>RI</b></p> <p><b>RR</b></p> <p><b>RS</b></p>	<p><u>S</u>et <u>C</u>alibration points</p> <p>SC allows you to calibrate the four internal sensors by entering low and high calibration points.</p> <p><b>SCxr1,v1,r2,v2,</b></p>
--	--

*(continued on next page)*



(continued)

<b>SC</b>	<p><b>Examples</b></p> <p>To set the inlet temperature calibration at 0 bits to 0 K and at 38991 bits to 297.9 K: SCI0,0,38991,2979</p> <p>To set the barometric pressure calibration at 15359 bits to 126.6 mb, and at 55291 bits to 994.4 mb: SCP15359,1266,55291,9944</p> <p>To set the relative humidity calibration at 250 bits to 12.0% and at 645 bits to 96.8%: SCR250,120,645,968</p> <p>To set the sample temperature calibration at 0 bits to 0 °, and at 39175 bits to 297.9 °: SCS0,0,39175,2979</p> <p>To echo the current inlet temperature calibration: SCI</p>
-----------	---


<b>SD</b>	<p><u>Set Data Delimiter</u></p> <p>SD sets the character used to delimit data in the output responses.</p> <p><b>SDx</b></p>
-----------	---

<b>SK</b>	<p><u>S</u>et <u>K</u> calibration constants</p> <p>SK allows you to enter calibration constants for each of the detection colors (blue, green, and red).</p> <p><b>SK</b><i>caaaaa,b.bbbe-b,c.ccce-c,ddde-d</i></p>
-----------	--


<b>SL</b>	<p><u>S</u>et calibration <u>L</u>abel string</p> <p>SL enters an ASCII calibration string into the Nephelometer, where the string is up to 80 characters terminated with a carriage return.</p> <p><b>Examples</b></p> <p>To enter the string "Last calibration 12/25/93" :</p> <pre>SLLast Calibration 12/25/93</pre> <p>To echo the calibration string:</p> <pre>SL</pre>
-----------	--

<b>SMB</b>	<p><u>S</u>et <u>M</u>ode for <u>B</u>ackscatter shutter</p> <p><b>SMBz</b>  where:  z = backscatter mode (0 or 1)  0 sets the Nephelometer to measure only in total scatter mode, where the backscatter shutter is parked in total scatter position.  1 enables operation of the backscatter shutter, where both total scatter and backscatter are measured.</p> <p><b>Examples</b>  To measure only the total scatter:  SMB0</p> <p>To enable the backscatter shutter, measuring total scatter and backscatter:  SMB1</p> <p>To echo the backscatter mode:  SMB</p>
------------	---

<b>see  also  STB  STP  STZ  SZ  V  Z</b>	<p><u>S</u>et <u>M</u>ode for auto<u>Z</u>ero baseline measurement</p> <p>SMZ sets the method used to zero baseline drift using filtered air. See Figures 6-2 and 6-3.  <b>Note:</b> <i>The duration of zero is set by the STZ command.</i></p> <p><b>SMZv</b>  where:  v = autozero mode (0 to 24)  0 = manual mode where a zeroing of the baseline only occurs when a Z command is given.  1 = normal mode where an autozero is performed at intervals set by the STP command.  2 - 24 = air-chop mode, which is similar to normal mode (1), except that the last "v" measured autozero baselines are averaged together.</p> <p><b>Examples</b>  To set the autozero mode to air-chop, with a baseline zero value based on the average of the last twelve readings:  SMZ12</p> <p>To echo the autozero mode:  SMZ</p>
---	---



**Figure 6-2**  
Timelines for Data Measurement



**Figure 6-3**  
Comparison of Manual, Normal, and Air Chop Modes

<b>SP</b>  <b>see also</b> <b>F</b> <b>L</b>	<p><u>S</u>et <u>l</u>amp <u>P</u>ower</p> <p>SP sets the input power of the lamp.</p> <p><b>SP</b><i>www</i></p> <p>where:  <i>www</i> = watts (0 to 150)</p> <p><b>Examples</b></p> <p>To set the input lamp power to 75 watts:      SP75</p> <p>To echo the lamp power setting:      SP</p>
--	--

<b>STA</b>	<p><u>S</u>et <u>T</u>ime for <u>A</u>veraging</p> <p>STA sets the running average length, in seconds, over which the sample is measured. STA also affects the time between unpolled reports. See Figure 6-2.</p> <p><b>STA</b><i>ttt</i></p> <p>where:</p> <p><i>ttt</i> = averaging time (1 - 9960 seconds)</p> <p>Internally, the Nephelometer uses a sliding average called "boxcar," where the latest value bumps the earliest value from the average. Up to 300 values can be stored in each boxcar. The following shows averaging time versus boxcar size:</p> <table style="margin-left: auto; margin-right: auto;"> <thead> <tr> <th style="text-align: left;"><u>Averaging Time</u></th> <th style="text-align: left;"><u>Boxcar Size</u></th> </tr> </thead> <tbody> <tr> <td>1 to 300 seconds</td> <td>1 second</td> </tr> <tr> <td>302 to 600 seconds</td> <td>2 seconds</td> </tr> <tr> <td>603 to 900 seconds</td> <td>3 seconds</td> </tr> <tr> <td>906 to 1800 seconds</td> <td>6 seconds</td> </tr> <tr> <td>1810 to 3000 seconds</td> <td>10 seconds</td> </tr> <tr> <td>3020 to 6000 seconds</td> <td>20 seconds</td> </tr> <tr> <td>6030 to 9000 seconds</td> <td>30 seconds</td> </tr> <tr> <td>9060 to 9960 seconds</td> <td>60 seconds</td> </tr> </tbody> </table> <p><b>Note:</b> <i>Averaging time entered must be an integer multiple of the boxcar size for the range in which the averaging time falls. If not, the Nephelometer will use the next smaller valid averaging time.</i></p> <p>Read commands RD and RP report new data only as it changes in the running average. Using the information above, if the averaging time is set to 3600 seconds, the RD and RP will reflect a change every 20 seconds. The analog outputs are also updated once each boxcar.</p> <p><b>Examples</b></p> <p>To set the averaging time to 60 seconds:</p> <pre>STA60</pre> <p>To echo the averaging time:</p> <pre>STA</pre>	<u>Averaging Time</u>	<u>Boxcar Size</u>	1 to 300 seconds	1 second	302 to 600 seconds	2 seconds	603 to 900 seconds	3 seconds	906 to 1800 seconds	6 seconds	1810 to 3000 seconds	10 seconds	3020 to 6000 seconds	20 seconds	6030 to 9000 seconds	30 seconds	9060 to 9960 seconds	60 seconds
<u>Averaging Time</u>	<u>Boxcar Size</u>																		
1 to 300 seconds	1 second																		
302 to 600 seconds	2 seconds																		
603 to 900 seconds	3 seconds																		
906 to 1800 seconds	6 seconds																		
1810 to 3000 seconds	10 seconds																		
3020 to 6000 seconds	20 seconds																		
6030 to 9000 seconds	30 seconds																		
9060 to 9960 seconds	60 seconds																		

<b>STB</b>  <b>see also SMZ STP STZ SZ V Z</b>	<p><u>S</u>et <u>T</u>ime to <u>B</u>lank valve</p> <p>STB sets the blanking time, in seconds, when the valve is switching. This includes the turning time of the valve as well as the purge time of the Nephelometer. During the STB interval, no data is taken. See Figures 6-2 and 6-3.</p> <p><b>STB</b><i>ttt</i></p> <p>where:  <i>ttt</i> = blanking time (5 - 999 seconds)</p> <p><b>Note:</b> <i>The STB value must be less than (STP - STZ)/2</i></p> <p><b>Examples</b></p> <p>To set the blanking time to 30 seconds:  STB30</p> <p>To echo the blanking time:  STB</p>
--	---

<b>STP</b>  <b>see also SMZ STB STZ SZ V Z</b>	<p><u>S</u>et autozero <u>P</u>eriod</p> <p>STP sets the time, in seconds, between autozeros. The value set with this command is only used if the Nephelometer autozero mode is normal or air-chop (SMZ &gt; 0). See Figures 6-2 and 6-3.</p> <p><b>STP</b><i>tttt</i></p> <p>where:  <i>tttt</i> = the range (10 - 9999 seconds)</p> <p><b>Note:</b> <i>The STP value must be greater than 2(STB) + STZ.</i></p> <p><b>Examples</b></p> <p>To set the time between autozeros to 3600 seconds:  STP3600</p> <p>To echo the autozero period:  STP</p>
--	--

<b>STT</b>	<p><u>S</u>et <u>T</u>ime and date</p> <p>STT sets the time and date of the Nephelometer's internal clock.</p> <p><b>STT</b><i>yyyy,mm,dd,hh,nn,ss</i>  where:  <i>yyyy</i> is the year (1994 - 2999)  <i>mm</i> is the month (1 - 12)  <i>dd</i> is the date (1 - 31)  <i>hh</i> is the hour (0 - 23)  <i>nn</i> is the minute (0 - 59)  <i>ss</i> is the second (0 - 59)</p> <p><b>Examples</b></p> <p>To set the internal clock to December 7, 1994, 10 pm.:</p> <pre>STT1994,12,7,22,0,0</pre> <p>To echo the internal clock time and date:</p> <pre>STT</pre>
------------	--

<b>STZ</b>  <b>see</b> <b>also</b> <b>SMZ</b> <b>STB</b> <b>STP</b> <b>SZ</b> <b>V</b> <b>Z</b>	<p><u>S</u>et <u>T</u>ime for auto<u>Z</u>ero measurement</p> <p>STZ sets the time, in seconds, the Nephelometer spends measuring filtered air during a zero baseline measurement. See Figures 7-1 and 7-2.</p> <p><b>STZ</b><i>ttt</i>  where:  <i>ttt</i> = the range (1 - 9999 seconds)</p> <p><b>Note:</b> <i>The STZ value must be less than STP - 2(STB).</i></p> <p><b>Examples</b></p> <p>To set the zero baseline measurement to 300 seconds:</p> <pre>STZ300</pre> <p>To echo the zero baseline measurement:</p> <pre>STZ</pre>
--	---



<b>SV</b>	<p><u>S</u>et <u>V</u>oltage on the photomultiplier tube</p> <p>SV sets the voltage that drives the photomultiplier tube (PMT). This voltage has an effect on the PMT's gain and noise level.</p> <p><b>SVcbbb</b>  where:  c = color channel  B= blue  R = red  G = green  bbb = voltage level applied to the PMT (0 -1200)</p> <p><b>Examples</b>  To set the blue channel PMT to 850 volts:  SVB850</p> <p>To echo the red channel PMT voltage:  SVR</p>
-----------	---

<b>SX</b>  see also <b>SA</b> <b>SB</b> <b>RX</b>	<p><u>S</u>et <u>e</u>Xternal host analog value</p> <p>SX sets the external host analog value to a specified voltage. This value appears on the PROGRAMMABLE I/O hardware output (factory default) and can be configured with the SB command.</p> <p><b>SXvvv</b>  where:  vvv = output voltage in millivolts (0 - 5000 mv for Programmable I/O (BNC))</p> <p><b>Examples</b>  To set the external host analog value to 2.50 volts:  SX2500</p> <p>To set the external host analog value to 5.0 volts:  SX5000</p> <p><b>Note:</b> Auxiliary Programmable I/O will read maximum (5 volts).</p> <p>To echo the auxiliary PROGRAMMABLE I/O output:  SX</p>
---	--

<b>SZ</b>  <b>see</b> <b>also</b> <b>SMZ</b> <b>STB</b> <b>STP</b> <b>STZ</b> <b>V</b> <b>Z</b>	<u>S</u> et <u>Z</u> ero baseline (background)  SZ allows you to externally set zero baseline (background) values, normally measured during an autozero measurement. SZ also clears previous zero data taken in air-chop mode. In air-chop mode, the supplied values are first in the running average of zero baseline values.
--	--

## Action Commands

Action commands control mechanical components of the Nephelometer. If you enter an action command without a parameter, the mechanical state is echoed back.

<b>B</b>	<p><u>B</u>lower control</p> <p><b>Bxxx</b> where: xxx = value between 0 and 255, 0 turns the blower off; 255 represents full power.</p> <p><b>Example</b> To set the blower to one-half power: B128</p> <p>To echo the blower power: B</p>
<b>F</b>  <b>see also</b> <b>L</b> <b>SP</b>	<p><u>F</u>an control</p> <p><b>Fc</b> where: c = 1 or 0. 1 turns the cooling fan on; 0 turns the cooling fan off.</p> <p><b>Note:</b> <i>Turning off the fan also reduces the lamp power to zero.</i></p>
<b>H</b>	<p><u>H</u>heater control</p> <p><b>Hc</b> where: c = 1 or 0. 1 activates the heater control. The heater then operates to maintain the sample temperature at or above the inlet temperature. 0 turns the heater control off.</p>
<b>L</b>  <b>see also</b> <b>F</b> <b>SP</b>	<p><u>L</u>amp control</p> <p><b>Lc</b> where: c = 1 or 0 1 turns the lamp on to the power set by the SP command; 0 turns the lamp off.</p>

<b>PD</b>	<p><u>P</u>ower <u>D</u>own</p> <p>PD places the Nephelometer in a minimum power ("sleep") mode. When the PD command is given, this sequence occurs:</p> <ol style="list-style-type: none"> <li>1. Lamp power off.</li> <li>2. Heater power off (if on).</li> <li>3. Chopper and shutter motors stop.</li> <li>4. Valve switched to filtered air position (zero).</li> <li>5. After a 5 second delay, fan power off.</li> <li>6. Blower power off.</li> </ol> <p><b>PD</b></p>
-----------	--

<b>PU</b>	<p><u>P</u>ower <u>U</u>p</p> <p>PU is the only command accepted after a PD (Power Down) command. When the PU command is given, the hardware is restored to the operating state present before the PD command was given. When the PU command is given, this sequence occurs:</p> <ol style="list-style-type: none"> <li>1. Blower power on</li> <li>2. Fan power on.</li> <li>3. Lamp power on.</li> <li>4. Chopper and shutter motors on.</li> <li>5. Heater power on (if on before the PD command).</li> <li>6. Ten seconds after the PU command is given and after the heater has stabilized, the Nephelometer begins an autozero measurement and resumes operation in normal mode. If the Nephelometer was in manual mode before the PD command, the valve remains in the filtered air position.</li> </ol> <p><b>PU</b></p>
-----------	--

<b>V</b>  <b>see also SMZ STB STP STZ SZ Z</b>	<u>V</u> alve position  <b>Vc</b> where: $c = N$ for <u>N</u> ormal position and $Z$ for <u>Z</u> eroing position (filtered air) VN changes the position of the valve to the normal measurement position. VN clears the averaged measured scattering values, but zero scattering values are not adjusted. VZ overrides the autozero mode and the average measured scattering values are cleared.  <b>Note:</b> A "FAULT" message (instead of "ERROR") is returned after the V command if the valve position is not completed within 10 seconds.  <b>Example</b> To echo the current position of valve NORMAL, ZERO, or FAULT:  v
--	---

<b>Z</b>  <b>see also SMZ STB STP STZ SZ V</b>	<u>Z</u> ero command  Z causes the Nephelometer to perform a zero baseline measurement. In air-chop mode, the zero baseline measurement is averaged with the previous baseline data as defined in the SMZ command.  <b>Z</b>
--	--

## Read Commands (Polled)

Read commands are polled, which means the Nephelometer sends data in response to a specific request from the computer

**Note:** *Some of the commands directly affect or are affected by other commands. Refer to other commands where indicated.*

<b>RA</b>	<p><u>R</u>ead raw <u>A</u>nalog bit values</p> <p>RA allows the reading of pressure, sample temperature, inlet temperature, and relative humidity in analog bit value format.</p> <p><b>RA</b></p> <table><thead><tr><th><b>Example</b></th><th><b>Response</b></th></tr></thead><tbody><tr><td>RA</td><td>26539,13245,13456,980</td></tr></tbody></table>	<b>Example</b>	<b>Response</b>	RA	26539,13245,13456,980
<b>Example</b>	<b>Response</b>				
RA	26539,13245,13456,980				
<b>RB</b>	<p><u>R</u>ead <u>B</u>arometric pressure (mbar)</p> <p><b>RB</b></p> <table><thead><tr><th><b>Example</b></th><th><b>Response</b></th></tr></thead><tbody><tr><td>RB</td><td>1013.2</td></tr></tbody></table>	<b>Example</b>	<b>Response</b>	RB	1013.2
<b>Example</b>	<b>Response</b>				
RB	1013.2				
<b>RD</b>	<p><u>R</u>ead scatter <u>D</u>ata</p> <p>The response to the RD command is in the Data Record (D). See the "Unpolled Record Formats" section in this chapter.</p> <p><b>RD</b></p>				

<b>RF</b>	<u>Read status Flags</u>																																																								
	<p>RF returns a four-character hexadecimal value representing the state of the Nephelometer. The values for the sixteen flags are as follows:</p> <table border="1"> <thead> <tr> <th><u>Bit Position</u></th> <th><u>0 (Normal)</u></th> <th><u>1(Fault)</u></th> </tr> </thead> <tbody> <tr> <td>0000 0000 0000 0001</td> <td>Lamp at power</td> <td>Not within 10% of SP setting</td> </tr> <tr> <td>0000 0000 0000 0010</td> <td>Valve Ok</td> <td>Valve fault</td> </tr> <tr> <td>0000 0000 0000 0100</td> <td>Chopper Ok</td> <td>Chopper fault</td> </tr> <tr> <td>0000 0000 0000 1000</td> <td>Shutter Ok</td> <td>Shutter fault</td> </tr> <tr> <td>0000 0000 0001 0000</td> <td>Heater Ok or inactive</td> <td>Heater active but not stabilized</td> </tr> <tr> <td>0000 0000 0010 0000</td> <td>Pressure within range</td> <td>Pressure out of range</td> </tr> <tr> <td>0000 0000 0100 0000</td> <td>Sample Temp within range</td> <td>Sample Temp out of range</td> </tr> <tr> <td>0000 0000 1000 0000</td> <td>Inlet temp within range</td> <td>Inlet temp out of range</td> </tr> <tr> <td>0000 0001 0000 0000</td> <td>RH within range</td> <td>RH out of range</td> </tr> <tr> <td>0000 0010 0000 0000</td> <td>Unused</td> <td></td> </tr> <tr> <td>0000 0100 0000 0000</td> <td>Unused</td> <td></td> </tr> <tr> <td>0000 1000 0000 0000</td> <td>Unused</td> <td></td> </tr> <tr> <td>0001 0000 0000 0000</td> <td>Unused</td> <td></td> </tr> <tr> <td>0010 0000 0000 0000</td> <td>Unused</td> <td></td> </tr> <tr> <td>0100 0000 0000 0000</td> <td>Unused</td> <td></td> </tr> <tr> <td>1000 0000 0000 0000</td> <td>Unused</td> <td></td> </tr> </tbody> </table> <p><b>RF</b></p> <table border="1"> <thead> <tr> <th><b>Examples</b></th> <th><b>Responses</b></th> </tr> </thead> <tbody> <tr> <td>RF</td> <td>0002 indicates all flags are normal and the valve is in an unknown position.</td> </tr> <tr> <td>RF</td> <td>0080 indicates the measured inlet temperature is out of range, possibly due to improper calibration or sensor failure.</td> </tr> </tbody> </table>	<u>Bit Position</u>	<u>0 (Normal)</u>	<u>1(Fault)</u>	0000 0000 0000 0001	Lamp at power	Not within 10% of SP setting	0000 0000 0000 0010	Valve Ok	Valve fault	0000 0000 0000 0100	Chopper Ok	Chopper fault	0000 0000 0000 1000	Shutter Ok	Shutter fault	0000 0000 0001 0000	Heater Ok or inactive	Heater active but not stabilized	0000 0000 0010 0000	Pressure within range	Pressure out of range	0000 0000 0100 0000	Sample Temp within range	Sample Temp out of range	0000 0000 1000 0000	Inlet temp within range	Inlet temp out of range	0000 0001 0000 0000	RH within range	RH out of range	0000 0010 0000 0000	Unused		0000 0100 0000 0000	Unused		0000 1000 0000 0000	Unused		0001 0000 0000 0000	Unused		0010 0000 0000 0000	Unused		0100 0000 0000 0000	Unused		1000 0000 0000 0000	Unused		<b>Examples</b>	<b>Responses</b>	RF	0002 indicates all flags are normal and the valve is in an unknown position.	RF
<u>Bit Position</u>	<u>0 (Normal)</u>	<u>1(Fault)</u>																																																							
0000 0000 0000 0001	Lamp at power	Not within 10% of SP setting																																																							
0000 0000 0000 0010	Valve Ok	Valve fault																																																							
0000 0000 0000 0100	Chopper Ok	Chopper fault																																																							
0000 0000 0000 1000	Shutter Ok	Shutter fault																																																							
0000 0000 0001 0000	Heater Ok or inactive	Heater active but not stabilized																																																							
0000 0000 0010 0000	Pressure within range	Pressure out of range																																																							
0000 0000 0100 0000	Sample Temp within range	Sample Temp out of range																																																							
0000 0000 1000 0000	Inlet temp within range	Inlet temp out of range																																																							
0000 0001 0000 0000	RH within range	RH out of range																																																							
0000 0010 0000 0000	Unused																																																								
0000 0100 0000 0000	Unused																																																								
0000 1000 0000 0000	Unused																																																								
0001 0000 0000 0000	Unused																																																								
0010 0000 0000 0000	Unused																																																								
0100 0000 0000 0000	Unused																																																								
1000 0000 0000 0000	Unused																																																								
<b>Examples</b>	<b>Responses</b>																																																								
RF	0002 indicates all flags are normal and the valve is in an unknown position.																																																								
RF	0080 indicates the measured inlet temperature is out of range, possibly due to improper calibration or sensor failure.																																																								

<b>RI</b>	<u>Read Inlet temperature (K)</u>			
	<b>RI</b>			
	<table border="1"> <thead> <tr> <th><b>Example</b></th> <th><b>Response</b></th> </tr> </thead> <tbody> <tr> <td>RI</td> <td>295.4</td> </tr> </tbody> </table>	<b>Example</b>	<b>Response</b>	RI
<b>Example</b>	<b>Response</b>			
RI	295.4			

<b>RL</b>	<u>Read Lamp voltage and current (volts and amps)</u>			
	<b>RL</b>			
	<table border="1"> <thead> <tr> <th><b>Example</b></th> <th><b>Response</b></th> </tr> </thead> <tbody> <tr> <td>RL</td> <td>12.5,4.1 indicates the lamp is operating at 12.5 volts and 4.1 amps.</td> </tr> </tbody> </table>	<b>Example</b>	<b>Response</b>	RL
<b>Example</b>	<b>Response</b>			
RL	12.5,4.1 indicates the lamp is operating at 12.5 volts and 4.1 amps.			

<b>RN</b>	<u>Read Noise</u> (sensitivity) levels
	Read noise levels based on photon counting statistics. This command requires clean filtered air and varies with the lamp power (SP), averaging time (STA), and PMT voltages (SV).
	<b>RN</b>
	<b>Example                      Response</b> RN                                  ????.

<b>RO</b>	<u>Read accumulated On</u> time of Nephelometer (not including power down time)
	<b>RO</b>
	<b>Example                      Response</b> RO                                  3425,48 indicates the Nephelometer has been on for 3425 hours and 48 minutes.

<b>RP</b>	<u>Read all Photon</u> counts (for blue, green, and red)
	<b>RP</b>
	<b>Examples                      Response</b> RP                                  B, G, or R photon count records. See the Photon Count Records in the "Unpolled Record Formats" section of this chapter.

<b>RPG</b>	<u>Read Photon</u> counts for <u>Green</u> (only)
	<b>RPG</b>
	<b>Example                      Response</b> RPG                                  See the Photon Count Records in the "Unpolled Record Formats" section of this chapter.

<b>RR</b>	<u>Read Relative</u> Humidity (0 - 99.9%)
	<b>RR</b>
	<b>Example                      Response</b> RR                                  15.6



<b>RS</b>	Read <u>S</u> ample temperature (K)
	<b>RS</b>
	<b>Example            Response</b>
	RS                    298.0

<b>RT</b>	Read <u>T</u> ime and date
	<b>RT</b>
	<b>Example            Response</b>
	RT                    See the Time Record (T) in the "Unpolled Record Formats" section of this chapter.

<b>RU</b>	Read all enabled <u>U</u> npolled records
	<b>RU</b>
	<b>Example            Response</b>
	RU                    This command will send all unpolled records that have been enabled with the corresponding unpolled command.

<b>RV</b>	Read firmware <u>V</u> ersion
	<b>RV</b>
	<b>Example            Response</b>
	RV                    Version 1.3 July 4, 1994

<b>RX</b>  see also <b>SA</b> <b>SB</b> <b>SX</b>	Read au <u>X</u> iliary PROGRAMMABLE I/O input port (millivolts)
	<b>RX</b>
	<b>Example            Response</b>
	RX                    2543 (indicating an input voltage of 2.543 volts)

<b>RY</b>	Read Auxiliar <u>Y</u> data record
	<b>RY</b>
	<b>Example            Response</b>
	RY                    See the Status Record (Y) in the "Unpolled Record Formats" section of this chapter.

<b>RZ</b>	<u>Read Zero</u> background data record	
	<b>RZ</b>	
	<b>Example</b> RZ	<b>Response</b> See the Zero Background Data Record (Z) in the "Unpolled Record Formats" section of this chapter

## Unpolled Commands

Use unpolled commands to cause the Nephelometer to automatically output data records at specific intervals. In unpolled mode, 0 disables a record and 1 enables the record. During unpolled operation, records that have been enabled are sent at the end of each averaging time.

**Note:** *UE is the only command allowed during unpolled operation. All other commands receive a “no response.”*

<b>UB</b>	<u>U</u> npolled operation <u>B</u> egins The only command accepted after UB is UE. <b>UB</b>
<b>UD</b>	<u>U</u> npolled <u>D</u> ata record UD0 disables the record; UD1 enables a record to be sent at the end of each averaging time. See Data Record (D) in the "Unpolled Record Formats" section of this chapter.
<b>UE</b>	<u>U</u> npolled mode <u>E</u> nds UE is the only command allowed during unpolled operation. All other commands receive a “no response.” <b>UE</b>
<b>UP</b>	<u>U</u> npolled <u>P</u> hoton counts <b>UP<sub>n</sub></b> where: $n = 0, 1, \text{ or } 3$ 0 = disabled 1 = sends photon counts for green channel only 3 = sends all three photon count records (blue, green, and red) See Photon Count Records in the "Unpolled Record Formats" section of this chapter.
<b>UT</b>	<u>U</u> npolled <u>T</u> ime and date UT1 enables the T record sent at the end of each averaging time; UT0 disables the record. See Time Record (T) in the "Unpolled Record Formats" section of this chapter.

<b>UY</b>	<p><u>U</u>npolled <u>A</u>uxiliar<u>Y</u> status data record</p> <p>UY1 enables the record sent at the end of each averaging time; UY0 disables the record. See Auxiliary Status Record (Y) in the "Unpolled Record Formats" section of this chapter.</p>
-----------	--

<b>UZ</b>	<p><u>U</u>npolled <u>Z</u>ero background data record</p> <p>UZ1 enables sending of the zero background data record one time after the end of a zero background measurement. UZ0 disables sending of the zero background data record.</p> <p>See Zero Background Data Record (Z) in the "Unpolled Record Formats" section of this chapter.</p>
-----------	--

## Unpolled Record Formats

The following are examples of unpolled record formats. These records are shown with the data delimiter set to space (SD1). Delimiters other than a space, cause the record to be packed with no added white space.

### Time Record (T)

T YYYY MM DD HH NN SS

where:

T = Time record  
 YYYY = Year  
 MM = Month  
 DD = Day  
 HH = Hour  
 NN = Minute  
 SS = Second

## Photon Count Records (B, G, and R)

B CCCCCCCC SSSSSSSS DDDDDD RRRRR LLLLLLLL MMMMMMMM NNNNNN VVVVV PPPP.P TTT.T

where:

B = Photon count record for B , R or G (Blue, Red, or Green)  
CCCCCCCC = Photon counts from calibrator (total scatter)  
SSSSSSSS = Photon counts from measure (total scatter)  
DDDDDD = Photon counts from dark  
RRRRR = Revolutions of chopper for total scatter measurement  
LLLLLLL = Photon counts from calibrator(unused backscatter cycle)  
MMMMMMMM = Photon counts from measure (backscatter)  
NNNNNN = Photon counts from dark (backscatter cycle)  
VVVVV = Revolutions of chopper for backscatter measurement  
PPPP.P = Pressure in millibar  
TTT.T = Sample temperature in degrees K

**Notes:** *In total scatter only mode L, M and N are set to 0 and RR reflects revolutions of the calibrator.*

*It is possible that for averaging times in excess of 1000 seconds, the photon count data will overflow the space delimited fixed format. In this case the format is extended to include the required extra digits.*

## Data Record (D)

D wxyz TTTT +B.BBBE-B +G.GGGe-G +R.RRRRe-R +A.AAAe-A +B.BBBE-B +C.CCCe-C

where:

D = Total scatter data record  
w = Current mode. N if in normal measurement mode,  
Z if in Zero mode; B if in blanking mode  
x = T if in total scatter mode; B if in backscatter mode  
y = X Currently unused  
z = X Currently unused  
TTTT = Time remaining in current state  
+B.BBBE-B = Scattering coefficient in blue  
+G.GGGe-G = Scattering coefficient in green  
+R.RRRRe-R = Scattering coefficient in red  
+A.AAAe-A = Scattering coefficient in blue (backscatter)  
+B.BBBE-B = Scattering coefficient in green (backscatter)  
+C.CCCe-C = Scattering coefficient in red (backscatter)

**Note:** In normal mode the Rayleigh scattering signal is subtracted to give the scattering coefficient, whereas in zero mode it is not. In blanking mode the scattering coefficients retain their value from the previous mode.

## Auxiliary Status Record (Y)

Y 9999999 PPP.P TTT.T III.I RR.R VV.V AA.A BBBB FFFF

where:

Y = Auxiliary status record  
9999999 = Sensitivity based on green channel (proton frequency)  
PPPP.P = Barometric pressure (mbar)  
TTT.T = Sample temperature (degrees K)  
III.I = Inlet temperature (degrees K)  
RR.R = Relative humidity (%)  
VV.V = Lamp voltage  
AA.A = Lamp current  
BBBB = BNC input voltage (millivolts)  
FFFF = Status flags (hex)

## Zero Background Data Record (Z)

Z +B.BBBE-B +G.GGGe-G +R.RRRe-R +A.AAAe-A +B.BBBE-B +C.CCCe-C  
+D.DDDe-D +E.EEEe-E +F.FFFe-F

where:

Z = Zero background record  
+B.BBBE-B = Scattering value from last zero (blue)  
+G.GGGe-G = Scattering value from last zero (green)  
+R.RRRe-R = Scattering value from last zero (red)  
+A.AAAe-A = Scattering value from last zero (blue  
backscatter)  
+B.BBBE-B = Scattering value from last zero (green  
backscatter)  
+C.CCCe-C = Scattering value from last zero (red  
backscatter)  
+D.DDDe-D = Rayleigh scattering value from last zero (blue)  
+E.EEEe-E = Rayleigh scattering value from last zero  
(green)  
+F.FFFe-F = Rayleigh scattering value from last zero (red)

## T, B, G, R, D, Y and Z Records

(Shown in order of data transmitted)

```
T 1994 12 31 23 59 59
B 99999999 99999999 999999 10990 99999999 999999999 999999 10991 1013.7 299.7
G 999999999 999999999 9999999 10990 9999999 999999999 999999 10991 1013.7 299.7
R 999999999 999999999 9999999 10990 9999999 999999999 999999 10991 1013.7 299.7
D NTXX 9999 +1.111e-1 +2.222e-2 +3.333e-3 +4.444e-4 +5.555e-5 +6.666e-6
Y 9999999 222.2 333.3 444.4 55.5 66.6 77.7 8888 9999
Z +1.111e-1 +2.222e-2 +3.333e-3 +4.444e-4 +5.555e-5 +6.666e-6 +7.777e-7 +8.888e-8
+9.999e-9
```

## How to Input Commands and Troubleshoot the Results

Use the following information as a guide to inputting software commands and for troubleshooting possible problems.

### Input Guidelines

- Input all alpha characters as capital letters (SMZ, *not* smz).
- Separate parameters with commas, not spaces.
- If you are in a command string, use the <Backspace> key to back up and make changes. Do *not* use <arrow> keys.
- At the end of a command string, press <Enter> to complete the string.



## Troubleshooting Input

Use Table 6-3 as a troubleshooting guide.

**Table 6-3**  
Troubleshooting Software Commands

<b>Symptom</b>	<b>Possible Problem</b>	<b>Refer to</b>
"Error" message after pressing <Enter>	<ul style="list-style-type: none"> <li>• An invalid command; command does not exist.</li> <li>• An invalid parameter, which includes too many parameters or a parameter that is out-of-range.</li> <li>• Incorrect syntax</li> </ul>	<ul style="list-style-type: none"> <li>• Table 6-3??? in this section.</li> <li>• The command showing the range and an example.</li> <li>• "Input Guidelines" in this section.</li> </ul>
No response after pressing <Enter>	<ul style="list-style-type: none"> <li>• In unpolled mode</li> <li>• Serial cable</li> <li>• Incorrect COM port</li> <li>• Incorrect baud rate</li> <li>• RS232 chip on the Nephelometer</li> <li>• Nephelometer is locked up</li> <li>• In power down mode (PD command)</li> </ul>	<ul style="list-style-type: none"> <li>• Use the UE command to exit unpolled mode. Reenter the command if an "OK" is returned.</li> <li>• Check the cable and the cable connection. See Chapter 2, "Unpacking and Setting Up the Hardware."</li> <li>• Check the COM port specified in the software.</li> <li>• Software must be set at 9600 baud. Also check computer hardware.</li> <li>• Contact TSI. Refer to Chapter 9, "Contacting Customer Service."</li> <li>• Remove power from the Nephelometer, then apply power to the instrument. If the problem continues, contact TSI.</li> <li>• Enter the PU (power up) command to start the Nephelometer. Reenter the command if an "OK" is returned.</li> </ul>



## CHAPTER 7

# Theory of Operation

This chapter contains the theory of operation for the Model 3550/3560 Series Integrating Nephelometer with these main sections:

- History
- Theory of Integration
- Signal Processing
- Detection Optics
- Calibration
- Theory of Wavelength


---

## History

**Note:** Background information for this section has been taken from the paper: “Three-Wavelength Nephelometer Suitable for Aircraft Measurement of Background Aerosol Scattering Coefficient,” Bodhaine, Ahlquist, & Schnell, *Atmospheric Environment*, Vol 25A, No 10, pp 2268-2276, 1991.

The integrating Nephelometer has been used extensively for measurement of the aerosol light scattering/extinction coefficient ( $\sigma_{sp}$ ) and visual range in ground-based and airborne applications. The technique of nephelometry was begun by Beuttell and Brewer in 1949, refined by Crosby and Koerber in 1963, and improved by Charlson and Ahlquist in 1967. Multiwavelength nephelometers, predecessors of the TSI nephelometers, were developed beginning in 1969. The high-sensitivity integrating nephelometer (Figure 7-1) was standardized and commercialized by TSI in 1993 by incorporating design improvements and technology advances developed over more than two decades.

The integrating nephelometer is a high-sensitivity device capable of detecting the scattering properties of aerosol particles. The nephelometer detects by measuring the light scattered by the aerosol and then subtracting light scattered by the walls of the measurement chamber, light scattered by the gas, and electronic noise inherent in the detectors.



**Figure 7-1**  
Nephelometer Schematic

The three-color detection version of the TSI nephelometers detects scattered light intensity at three wavelengths. Normally the scattered light is integrated over an angular range of 7–170° from the forward direction. But with the addition of the backscatter shutter feature to the Nephelometer, this range can be adjusted to either 7–170° or 90–170° to give total scatter and backscatter signals. This provides additional useful information about the particle scattering behavior.

---


## Construction

The integrating nephelometer is constructed primarily of aluminum parts for purposes of rigidity, ruggedness, and weight. The main body of the nephelometer consists of 10 cm diameter thin walled aluminum tubing, which is approximately 90 cm long. Aperture plates are set along the axis of the tube, held and spaced by 8 cm diameter tubing. Receiving optics are located at one end of the main tube with a light trap located at the other end to provide a very dark reference against which to view the light scattered by particles and gas.

A 75 watt quartz-halogen lamp, with a built-in elliptical reflector, provides illumination for the aerosol. The reflector focuses the light onto one end of an optical pipe where the light is carried into the internal cavity of the instrument. The optical pipe is used to thermally isolate the lamp from the sensing volume.

The output end of the optical light pipe is an opal glass diffuser that acts as a cosine (Lambertian) light source. Within the measuring volume, the first aperture on the detection side of the instrument limits the light integration to angles greater than  $7^\circ$ , measured from the horizontal at the opal glass. On the other side, a shadow plate limits the light to angles less than  $170^\circ$ . The measurement volume is defined by the intersection of this light with a viewing volume cone defined by the second and fourth aperture plates on the detection side of the instrument. The fourth aperture plate incorporates a lens to collimate the light scattered by aerosol particles so that it can be split into separate wavelengths.

The nephelometer uses a reference chopper to calibrate scattered signals. The chopper makes a full rotation 23 times per second. The chopper consists of three separate areas labeled: signal, dark, and calibrate (Figure 7-2).



**Figure 7-2**  
Reference Chopper Zones

The signal section covers 180° and simply allows all light to pass through unaltered. The dark section covers 100° and is a very black background that blocks all light. This section provides a measurement of the photomultiplier tube (PMT) background noise. The third section, calibrate, covers 80°. The light source directly illuminates this section providing a measure of lamp stability over time. To reduce the lamp intensity to a level that will not saturate the photomultiplier tubes, the calibrate section incorporates a neutral density filter that blocks approximately 99.9% of the incident light.

To subtract the light scattered by the gas portion of the aerosol, a high-efficiency particulate air (HEPA) filter is switched periodically in line with the inlet. This allows compensation for changes in the background scattering of the nephelometer, and in gas composition that will affect Rayleigh scattering of air molecules with time. The HEPA filter is switched into line using an automated valve, activated by a gear motor and positioned, using an optical position sensor. When the HEPA filter is not in line with the inlet, a small amount of filtered air leaks through the light trap to keep the apertures and light trap free of particles. A smaller HEPA filter allows a small amount of clean air to leak into the sensor end of the chamber between the lens and second aperture. This keeps the lens clean and confines the aerosol light scatter to the measurement volume only.

---

## Theory of Integration

The property of interest in using the integrating nephelometer is usually the extinction of light over a short distance. Light extinction is caused by the properties of both light absorption and light scattering caused by air or other gas and by the particles within the gas. The total extinction coefficient,  $\sigma$ , is the sum of these properties. The effect is generally described by the Beer-Lambert law:

$$\frac{I}{I_0} = e^{(-\sigma_{ext} x)} \quad \text{[equation 7-1]}$$

where:

$$\sigma_{ext} = \sigma_{scat} + \sigma_{abs} = \ln(I/I_0)/x$$

$x$  = distance in meters (length of light path)  
 $I$  = intensity of light after distance  
 $I_0$  = intensity of incident light

$\sigma_{scat}$  and  $\sigma_{abs}$  are called the scattering and absorption coefficients and are described by:

$$\sigma_{scat} = \sigma_{rg} + \sigma_{sp} \quad \text{[equation 7-2]}$$

$$\sigma_{abs} = \sigma_{ag} + \sigma_{ap} \quad \text{[equation 7-3]}$$

where:

$\sigma_{rg}$  = term for Rayleigh scattering  
 $\sigma_{sp}$  = term for scattering by particles  
 $\sigma_{ag}$  = term for absorption by gases  
 $\sigma_{ap}$  = term for absorption by particles

The integrating nephelometer measures the  $\sigma_{scat}$  term and generates the  $\sigma_{sp}$  term by subtracting the  $\sigma_{rg}$  term from  $\sigma_{scat}$  using (equation 7-2).

Since it is often assumed that  $\sigma_{sp} \gg \sigma_{ap}$  or that  $\sigma_{sp} \gg \sigma_{ap}$  for areas other than those having high concentrations of soot particles, the measure of  $\sigma_{sp}$  often gives a good estimate of the aerosol contribution to  $\sigma_{ext}$ . Also,  $\sigma_{ag}$  is usually negligible.


The basic geometry of the TSI integrating nephelometer was first described by Beuttel and Brewer (1949). Since this geometry integrates the intensity of the scattered light over close to  $4\pi$  steradians, it gives a good approximation of the scattering component of extinction,  $\sigma_{sp}$ . Figure 7-3 shows a schematic of this geometry which results in the integration of the angular scattering function,  $\beta(\phi)$ , such that the sensor detects the flux  $B$  due to the light source with intensity  $I_0$  located at a distance  $y$  from the detector axis:

$$B = I_0/y \int_{\phi_1}^{\phi_2} \beta(\phi) \sin \phi d\phi \quad \text{[equation 7-4]}$$

Since the nephelometer has a Lambert (cosine  $\phi$ ) source, and if we make the assumption that  $\phi_1 \approx 0$  and that  $\phi_2 \approx \pi$ , then:

$$B \approx (I_0/y) * (\sigma_{scat}/2\pi) \quad \text{[equation 7-5]}$$

These equations are derived formally by Middleton (1958) and by Butcher and Charlson (1972).



**Figure 7-3**  
Nephelometer Geometry



---

## Theory of Wavelength Dependence

The theory of wavelength dependence of light scattering is simple to define in the broad sense, and has substantial detail in application. The simple explanation is that smaller particles scatter differently than larger particles which results in different amounts of scattering at different wavelengths.

Light scattered by particles less than 0.05  $\mu\text{m}$  is described in relatively simple terms by Rayleigh's theory. Particles larger than 10  $\mu\text{m}$  must be described by geometric optics. Particles between these sizes can be analyzed using the much more complicated Mie scattering theory. The scattering of these particles is complicated because the size of particle and the wavelength of light have the same order of magnitude.

However, if we consider a polydisperse aerosol, which is usually the case for atmospheric aerosols, the wavelength dependence in the Mie region is dictated entirely by the size distribution and can be analyzed in a simpler fashion than the full blown Mie theory.

All materials have a describing characteristic called the index of refraction ( $m$ ). This is the ratio of the speed of light in a vacuum ( $c$ ) to the speed of light in a particular material ( $V_p$ ).

$$m = \frac{c}{V_p} \quad \text{[equation 7-6]}$$

$c \cong 3 \times 10^{10}$  cm/s, is the speed of light in a vacuum.

The scattering of particles is governed by the ratio of the particle size ( $d$ ) to the wavelength of light ( $\lambda$ ). This dimensionless number is called the size parameter  $a$ .

$$a = (\pi * d) / \lambda \quad \text{[equation 7-7]}$$

As can be seen from the equation, for  $a \ll 1$  particles are in the Rayleigh region and for  $a \cong 1$  they are in the Mie region.

Rayleigh's equation is as follows:

$$I_{\theta} = \frac{I_0 \pi^4 d^6 (m^2 - 1)^2}{8 R^2 \lambda^4 (m^2 + 2)^2} (1 + \cos^2 \theta) \quad d < 0.05 \mu m$$

[equation 7-8]

Where  $m$  is the index of refraction and  $R$  is the distance of the particle, at angle  $\theta$ , from the light.

It is already evident even for the case of Rayleigh scattering that particle size relative to wavelength is a key concern.

A simple description of the way that particle size distribution controls the wavelength dependence of  $\sigma_{sp}$  for  $a \cong 1$  can be stated, starting with the empirical observation that  $\sigma_{sp}$  for atmospheric aerosols is often close to a power-law function of wavelength:

$$\sigma_{sp} = C \lambda^{-\dot{a}}$$

[equation 7-9]

where:

$C$  = a quantity describing the amount of scattering  
 $\dot{a}$  = Angstrom exponent

The number distribution of atmospheric aerosol often can be crudely described as a power-law function of size:

$$dN/d\log r = c r^{-\beta}$$

[equation 7-10]

where:

$N$  = number concentration  
 $r$  = particle radius  
 $c$  = related to total concentration  
 $\beta$  = slope of the number distribution on a log-log plot.

If this is approximately true, which is usually the case for  $a \cong 1$ , then it can be shown that:

$$\dot{a} = \beta - 2$$

[equation 7-11]

Thus, what controls  $\dot{a}$  is the slope of the number distribution; in other words the relative amounts of fine and coarse particles. Typically values of  $\dot{a}$  for gases are very nearly 4, for urban aerosols are on the order of 2, for rural haze 1-2, and close to 0 for very coarse aerosols like wind blown dust or sea salt. For large values of  $\dot{a}$ , there will be a large difference in scattering at different wavelengths and for small values of  $\dot{a}$ , for example near 0, scattering at different wavelengths are nearly equal.

Because the Model 3550/3560 Series Nephelometers can calculate the logarithm of  $\sigma_{sp}$  as an analog output, a simple numerical exercise yields the quantity  $\dot{a}$ :

since

$$\begin{aligned}\sigma_{sp} &= C\lambda^{-\dot{a}} \\ \log \sigma_{sp} &= \log C - \dot{a} \log \lambda\end{aligned}\quad \text{[equation 7-12]}$$

so taking the derivative ( $d\log C = 0$ )

$$d\log \sigma_{sp} = -\dot{a} d\log \lambda \quad \text{[equation 7-13]}$$

Since the Nephelometer measures at different specific wavelengths, we can approximate the derivative with a finite difference:

$$\dot{a} \cong \frac{\Delta \log \sigma_{sp}}{\Delta \log \lambda} = \frac{\log \sigma_{sp}(\lambda_2) - \log \sigma_{sp}(\lambda_1)}{\log \lambda_2 - \log \lambda_1} \quad \text{[equation 7-14]}$$

Why measure  $\dot{a}$ ?

One of the main reasons to measure  $\sigma_{sp}$  is to quantify and understand how it governs the transfer of light through the air. Examples are the partial transmission of sunlight to the earth's surface, transmission of lasers or other light beams horizontally or the appearance of haze in front of distant vistas. Because  $\sigma_{sp}$  is a function of wavelength,  $\lambda$ , different amounts of blue versus red light are transmitted. Usually, atmospheric aerosols scatter more blue light than green or red light, that is  $\sigma_{sp}$  for blue light is greater than  $\sigma_{sp}$  for green light etc. However since  $\dot{a}$  is variable and visible light spans quite a range of wavelengths (roughly 400-700 nm) it is necessary to know  $\dot{a}$  in order to quantify the transmission of all the wavelengths.

Measurement of  $\dot{a}$  also provides a continuous record of variations in the size distribution via estimation of variations of  $\beta$ .

In conclusion for  $a \ll 1$ , e.g., for  $\lambda = 0.6 \mu\text{m}$ ,  $d \ll 0.1 \mu\text{m}$ , the Rayleigh theory is typically valid. The intensity of the Rayleigh scatter is proportional to the sixth power of the particle diameter and inversely proportional to the fourth power of the incident wavelength.

For  $\alpha \cong 1$ , where the particle diameter and the wavelength are of the same magnitude, say  $0.1\text{--}1 \mu\text{m}$ , atmospheric aerosol scattering can often be described by a simple power-law function, which is governed solely by the slope of the number distribution of the aerosol on a log-log plot.

---

## Signal Processing

This section describes the microprocessors and equations used in the Nephelometer firmware to calculate various scattering coefficients from the raw scatter data and calibration data.

### Description

All electronics are located inside the nephelometer. Two microprocessors are used to control all functions. A Motorola 68HC711D3 microprocessor is used to control the brushless DC motors used to spin the reference chopper and optional backscatter shutter. A Motorola 68HC16Z1 microprocessor is used to control all other signal processing and input/output operations. All setup and configuration is done through the serial interface to an external computer with all parameters stored in the Nephelometer's battery-backed-up RAM.

### Equations/Calculations

This section contains firmware and calibration equations for the Nephelometer.

## Firmware Calculations

The Nephelometer calculations begin with the raw signals received from the photomultiplier tubes (PMT) for each color, for each section of the reference chopper and for each rotation of the backscatter shutter. Therefore, there are eighteen values. In addition, the number of cycles of the chopper are given for each color. This is the data that is sent from the nephelometer as raw count data. Examples are shown in Table 7-1.

**Table 7-1**  
Examples of Raw Count Data

Color	Total Scatter				Backscatter			
	Calibrate	Signal	Dark	Cycles	Calibrate	Signal	Dark	Cycles
Blue	523939	12691	28	693	413847	6350	16	693
Green	1022163	12185	52	693	807927	6146	27	693
Red	514975	5271	1038	693	401071	3835	1021	693

These values are the raw photon counts for each portion of the reference chopper for each color. Since the angular width of each portion of the chopper is different, the values are normalized for gate width and number of cycles during the measurement. This converts the photon counts to cycles/second (Hz) using the equation:

$$C_s = (360 * C * S) / (G * N) \quad \text{[equation 7-15]}$$

where:  $C_s$  = scaled count rate (Hz)

$C$  = raw photon count rate

$S$  = speed of reference chopper (22.994 RPS)

$G$  = gate width in degrees

40° for calibrate section

140° for signal section

60° for dark section

$N$  = number of revs in measurement

**Notes:** The gate widths defined in equation 7-15 do not add up to 360° as one might imagine. This is because there is a 40° period between each of the sections to blank the overlap.

The calibrate signal taken during the backscatter cycle is not entirely used in the measurement cycle because the chopper is not illuminated during the backscatter cycle.

Table 7-2 shows the comparisons between gate width and physical dimensions of the three sections.

**Table 7-2**  
Comparisons Between Gate Width and Physical Dimensions

<b>Shutter Area</b>	<b>Gate Width [°]</b>	<b>Physical Width [°]</b>	<b>Blank Width [°]</b>
Calibrate	40	80	40
Signal	140	180	40
Dark	60	100	40
Total	240	360	120

**Note:** Blank width is 20° at the beginning and 20° at the end of each gate.

Next, each of the scaled count rates is corrected for dead time (due to photomultiplier pulse width) by the approximation:

$$F = C_s * (C_s * K_1 + 1) \quad \text{[equation 7-16]}$$

where:  $K_1$  = constant based on pulse width (picoseconds/counts)  
 $C_s$  = scaled count rate (Hz)  
 $F$  = corrected count rate (Hz)

When the numbers from the Table 7-1 are run through equation 7-15 and equation 7-16, the results are shown in Table 7-3.

**Table 7-3**  
Counts Normalize to Photon Frequency (Hz)

<b>Color</b>	<b>Total Scatter</b>			<b>Backscatter</b>		
	<b>Calibrate</b>	<b>Signal</b>	<b>Dark</b>	<b>Calibrate</b>	<b>Signal</b>	<b>Dark</b>
Blue	156950	1083	6	123890	542	3
Green	307105	1040	10	242430	524	5
Red	154257	450	207	120056	327	203

Next, the nephelometer calculates averages of the data using a simple boxcar running average with up to 300 elements. Each of the eighteen values in Table 7-4 are averaged, depending on the averaging time selected.

**Table 7-4**  
Boxcar Averages

<b>Selected Boxcar Average</b>	<b>Seconds per Boxcar</b>
1 to 300	1
302 to 600	2
603 to 900	3
906 to 1800	6
1810 to 3000	10
3020 to 6000	20
6030 to 9000	30
9060 to 18000	60

Next, the nephelometer uses the scaled and corrected count rates to calculate the raw scattering signals for each of the three colors. Periodically, the aerosol sampled by the nephelometer is filtered to determine the background scattering portion of the total scattered light. At all times, the nephelometer automatically monitors the temperature and pressure to allow correction of these environmental parameters. This is important since the Rayleigh scattering by gas molecules is dependent on temperature and pressure. Also, the user periodically performs a manual calibration with a span gas such as CO<sub>2</sub>.

The nephelometer calculates two raw scattering signals: one for the aerosol and the other for filtered air:

$$B_f = K_2 * (S_f - D_f) / (C_f - D_f) \quad \text{[equation 7-17]}$$

where:  $B_f$  = filtered air raw scatter signal  
 $K_2$  = constant determined by span gas  
 $C_f$  = calibrate (Hz) scaled count rate for filtered air  
 $S_f$  = signal (Hz) scaled count rate for filtered air  
 $D_f$  = dark (Hz) scaled count rate including the dark measured during total scatter and backscatter

$$B_a = K_2 * (S_a - D_a) / (C_a - D_a) \quad \text{[equation 7-18]}$$

where:  $B_a$  = aerosol raw scatter signal  
 $K_2$  = constant determined by span gas  
 $C_a$  = calibrate (Hz) scaled count rate for aerosol  
 $S_a$  = signal (Hz) scaled count rate for aerosol  
 $D_a$  = dark (Hz) scaled count rate for aerosol

The Rayleigh scattering is calculated for the filtered air measurement and for the clean air using the equations including temperature and pressure:

$$R_f = (K_3 * P_f * T_s) / (T_f * P_s) \quad \text{[equation 7-19]}$$

where:  $R_f$  = filtered air Rayleigh scatter  
 $K_3$  = Rayleigh scatter of air for given color  
 $P_f$  = pressure of filtered air  
 $T_s$  = standard temperature (273.2 K)  
 $T_f$  = temperature of filtered air  
 $P_s$  = standard pressure (1013.3 mbar)

$$R_a = K_3 * (P_a / T_a) * (T_s / P_s) \quad \text{[equation 7-20]}$$

where:  $R_a$  = aerosol Rayleigh scatter  
 $K_3$  = Rayleigh scatter of air for given color  
 $P_a$  = pressure of aerosol  
 $T_a$  = temperature of aerosol  
 $P_s$  = standard pressure (1013.3 mbar)  
 $T_s$  = standard temperature (273.2 K)

These parameters allow us to calculate a variety of interesting numbers for each of the three colors:

$$W = B_f - R_f \quad \text{[equation 7-21]}$$

where:  $W$  = filtered air wall scatter  
 $B_f$  = filtered air raw scatter signal  
 $R_f$  = filtered air Rayleigh scatter

$$B_s = B_a - W \quad \text{[equation 7-22]}$$

where:  $B_s$  = total scatter (particles + Rayleigh)  
 $B_a$  = aerosol raw scatter signal  
 $W$  = filtered air wall scatter

$$\sigma_{sp} = B_s - R_a \quad \text{[equation 7-23]}$$



where:  $\sigma_{sp}$  = aerosol scatter

$B_s$  = total scatter (particles + Rayleigh)

$R_a$  = aerosol Rayleigh scatter

Equations 7-21, 7-22, and 7-23 can also be calculated for backscatter only by using the backscatter signal from Table 7-3 and multiplying all Rayleigh scatter by  $K_4$ .

**Note:**  $W$  and  $R_a$  will be different for backscatter than for total scatter.

### Calibration Calculations

The previous equations use the four calibration factors:  $K_1$ ,  $K_2$ ,  $K_3$ , and  $K_4$ . These constants are obtained by calibration or measurement of the nephelometer.

The constant,  $K_1$ , is a one time calibration factor based on the measured pulse width of the photomultiplier tubes for each color. This parameter is set at the factory and used to correct for the dead time when a pulse is registered on the photomultiplier tube, preventing a second pulse from being registered.

$K_2$  is calculated during calibration using a equation similar to equation 7-8 and a table of wavelength-dependent, theoretical Rayleigh-scattering coefficients (from Table 2, Bodhaine, et. al.)

The constant  $K_3$  represents the known scattering value for air at a standard temperature and pressure condition for the wavelength in question.

$K_4$  is a constant near 0.5 representing the ratio of Rayleigh backscatter to Rayleigh total scatter.

Calculation of  $K_2$  using total scatter data (equations 7-22 and 7-23) can be expanded to:

$$\sigma_{sp} = B_l - W_t - R_l \quad \text{[equation 7-24]}$$

For low span gas

$$\sigma_{sp_l} = B_l - W_t - R_l = K_2 \frac{S_l - D_l}{C_l - D_l} - K_3 \frac{P_l T_s}{T_l P_s} - W_t = 0 \quad \text{[equation 7-25]}$$

For high span gas

$$\sigma_{sp_h} = B_h - W_t - R_h = K_2 \frac{S_h - D_h}{C_h - D_h} - K_{3_h} \frac{P_h T_s}{T_h P_s} - W_t = 0$$

[equation 7-26]

where:

- $B_l$  = low span gas raw scatter signal
- $B_h$  = high span gas raw scatter signal
- $R_l$  = low span gas Rayleigh scatter
- $R_h$  = high span gas Rayleigh scatter
- $W_t, W_b$  = wall scatter (total or backscatter)
- $C_l$  = low span gas calibration average scaled count rate
- $D_l$  = low span gas dark average scaled count rate
- $S_l$  = low span gas signal average scaled count rate
- $S_{lb}$  = low span gas backscatter signal average scaled count rate
- $P_l$  = low span gas calibration pressure
- $T_l$  = low span gas calibration temperature
- $K_{3_l}$  = low span gas scattering coefficients at standard pressure and temperature
- $C_h$  = high span gas calibration average scaled count rate
- $D_h$  = high span gas dark average scaled count rate
- $S_h$  = high span gas signal average scaled count rate
- $S_{hb}$  = high span gas backscatter signal average scaled count rate
- $P_h$  = high span gas calibration pressure
- $T_h$  = high span gas calibration temperature
- $K_{3_h}$  = high span gas scattering coefficients at standard pressure and temperature
- $P_s$  = standard pressure (273.2K)
- $T_s$  = standard temperature (1013.3 mbar)

Equations 7-25 and 7-26 are equal to zero because there are no particles present, and therefore, no scattering by particle ( $\sigma_{sp} = 0$ ). Scattering is only from the wallscatter and the Rayleigh scatter of the gas.

We have two equations, two unknowns  $K_2$  and  $W_t$  (wall). Subtract equation 7-25 from 7-26 to eliminate  $W_t$ .

To get:

$$K_2 \left[ \frac{S_h - D_h}{C_h - D_h} - \frac{S_l - D_l}{C_l - D_l} \right] - \frac{T_s}{P_s} \left[ K_{3h} \frac{P_h}{T_h} - K_{3l} \frac{P_l}{T_l} \right] \quad \text{[equation 7-27]}$$

Solve for  $K_2$ :

$$K_2 = \frac{\frac{T_s}{P_s} \left[ K_{3h} \frac{P_h}{T_h} - K_{3l} \frac{P_l}{T_l} \right]}{\frac{S_h - D_h}{C_h - D_h} - \frac{S_l - D_l}{C_l - D_l}} \quad \text{[equation 7-28]}$$

Generally, the low span gas in filtered air and the high span gas is  $\text{CO}_2$  or  $\text{CCl}_2\text{F}_2$ .

$K_4$  can be calculated using the backscatter data and  $K_2$  calculated from above.

$$\sigma_{sp_{lb}} = K_2 \left( \frac{S_{lb} - D_l}{C_l - D_l} \right) - K_4 K_{3l} \frac{P_l}{T_l} \frac{T_s}{P_s} - W_b = 0 \quad \text{[equation 7-29]}$$

$$\sigma_{sp_{hb}} = K_2 \left( \frac{S_{hb} - D_h}{C_h - D_h} \right) - K_4 K_{3h} \frac{P_h}{T_h} \frac{T_s}{P_s} - W_b = 0 \quad \text{[equation 7-30]}$$

Again, subtract equation 7-29 from 7-30 to eliminate  $W_b$  (wall).

Solve for  $K_4$ :

$$K_4 = \frac{K_2 \left[ \frac{S_{hb} - D_h}{C_h - D_h} - \frac{S_{lb} - D_l}{C_l - D_l} \right]}{\frac{T_s}{P_s} \left[ \frac{P_h}{T_h} K_{3h} - \frac{P_l}{T_l} K_{3l} \right]} \quad \text{[equation 7-31]}$$

---

## Detection Optics

The detection optics depend on the model of the nephelometer. The single wavelength unit (TSI Models 3551/61) consists of a broadband-coated 400 mm focal-length lens used to collimate the diverging light defined by apertures in the body of the instrument. This light is then passed through a 40 nm bandpass filter centered at 550 nm before entering a photomultiplier tube (green channel).

The three wavelength instrument (TSI Models 3553/63) contains the same lens as the single wavelength device. From the lens, the scattered light is separated by dichroic filters to three bandpass filters and separate photomultiplier tubes. The first wavelength is reflected by a color splitter that passes 500-800 nm light and reflects 400-500 nm light, through a 40 nm bandpass filter, centered at 450 nm, into a photomultiplier (blue channel). Light that passes through the first color splitter is split again by a filter that passes 500-600 nm light and reflects 600-800 nm light, through a 40 nm bandpass filter centered at 700 nm to a PMT (red channel). The light that passes through both color splitters (500-600 nm) continues to a 40 nm bandpass filter, centered at 550 nm and a PMT (green channel).

---

## Calibration

The Nephelometer is a photon counting instrument. Particles and gas molecules scatter photons that are detected and counted using highly sensitive photomultiplier tubes. These photon counts are converted into counting frequencies and correlated to light scattering coefficients using calibration constants. These constants are determined by filling the nephelometer's interior measurement volume with two span gases, one at a time, that have largely different scattering coefficients and measuring a few key parameters. By measuring the photon frequency, temperature, and pressure for each span gas measurement, and knowing the scattering coefficient for each of the span gases at a known temperature and pressure (e.g., STP), and using two simple equations, the calibration constants are easily calculated. For more detail on the calculations see "Calibration Calculations" earlier in this chapter.

# CHAPTER 8

## Performing Maintenance

Use the procedures in this chapter to maintain the Model 3550/3560 Series Integrating Nephelometer. Most maintenance procedures do not have specific time guidelines, but should be performed based on software or hardware failures, degradation in overall performance, or as a result of the special environments in which the Nephelometer may be used.

Table 8-1 gives an overview of the procedures in this chapter. The organization of the table reflects the order in which the procedures are presented.

**Table 8-1**  
Maintenance Overview

<b>Maintenance Procedure</b>	<b>Perform:</b>
<b>With the top and bottom covers on:</b>	
Calibrate the Nephelometer	Periodically, or if the reference chopper is dirty or scratched.
Clean or replace the fan filters	If a visual check shows the filters are clogged.
<b>With the top cover off:</b>	
Replace the main microprocessor EPROM	As part of a firmware update.
Replace the motor control microprocessor	As part of a firmware update.
Replace the lamp	If there is a status failure in the software, or periodically based on the power rating and the hours of use.
Replace the aerosol filters	Periodically, or if there is a significant drop in pressure.
Clean the light pipe lens	As needed, depending on the cleanliness of the environment in which the Nephelometer is operated.
Check for leaks	If the instrument is disassembled and is to be used in a low pressure environment.
<b>With the top and bottom covers off:</b>	
Clean the reference chopper	If the calibrate signal changes dramatically while using the lamp at a constant power.
Clean the light pipe outlet and the backscatter shutter	As needed, depending on the cleanliness of the environment in which the Nephelometer is operated.
Clean the flocked paper	If the background signal rises over time.
<b>With the photomultiplier tube (PMT) cover off:</b>	
Replace PMTs	If a failure occurs.
Check, clean, or replace bandpass filters	Periodically, or if using in a high-humidity environment.
Clean the lens	If there is an overall reduction in the signal of all three wavelengths over time.

---

## Removing Nephelometer Covers: Top, Bottom, and PMT

This section give instructions for removing Nephelometer top, bottom or PMT covers, necessary for most maintenance or hardware troubleshooting procedures (Figure 8-1).

**Note:** *The steps needed to remove covers are also included, when appropriate, in each maintenance procedure.*



**Figure 8-1**  
Nephelometer Showing Top, Bottom, and PMT Covers

### Removing the Top Cover

To remove the top cover of the Nephelometer:


1. Remove power from the Nephelometer and from any instruments connected to the Nephelometer.
2. Loosen the four screws fastening the Nephelometer cover and remove the cover (Figure 8-2).



### Caution

To avoid damage to the electronic circuitry, use electrostatic discharge (ESD) precautions when removing the top or bottom Nephelometer cover:

- Use only a table top with a grounded conducting surface.
- Wear a grounded, static-discharging wrist strap.



**Figure 8-2**

Top Cover Screws, Bottom Cover Screws, Fan Cover Screws, and I/O Panel Screws

## Removing the Bottom Cover

To remove the bottom cover of the Nephelometer:

- 1.** Remove power from the Nephelometer and from any instruments connected to the Nephelometer.
- 2.** Remove the top cover using the procedure in this section.



### Caution

To avoid damage to the electronic circuitry, use electrostatic discharge (ESD) precautions when removing the top or bottom Nephelometer cover:

- Use only a table top with a grounded conducting surface.
- Wear a grounded, static-discharging wrist strap.

3. The bottom cover of the Nephelometer is pliable and fits very tightly around the bottom of the instrument. Remove the bottom cover in one of two ways:
- The more difficult method is to pry off the bottom cover without removing the I/O panel and the fan panel. Unscrew the four bottom cover screws (Figure 8-2). *Carefully* remove the bottom cover, working alternately at both ends. Pry the edges of the cover apart, and ease the cover away from the body of the Nephelometer.
  - The easier method takes more time. After removing the I/O panel, the fan panel and the end plate, slide the bottom cover off the body of the instrument. Using Figure 8-2 as a reference:
    - a. Unscrew the three screws securing the I/O panel to the backbone, loosen the panel, and let the panel hang over the side of the instrument.
    - b. Unscrew the two screws securing the fan panel to the backbone, loosen the panel, and let the panel hang over the side of the instrument.
    - c. Unscrew the three screws securing the end plate to the inlet section of the instrument and remove the end plate.
    - d. Unscrew the four bottom cover screws, and *carefully* slide the cover off the end of the instrument.
- Note:** *Attach the bottom cover of the Nephelometer reversing the steps in step 3 of this procedure. Remember that the cover is pliable and fits tightly around the bottom of the instrument.*



## Removing the PMT Cover

To remove the PMT cover:


1. Remove power from the Nephelometer and from any instruments connected to the Nephelometer.
2. Unscrew the four screws attaching the PMT cover to the sensor and remove the cover (Figure 8-3).



### Caution

To avoid damage to the electronic circuitry, use electrostatic discharge (ESD) precautions when removing the Nephelometer cover:

- Use only a table top with a grounded conducting surface.
- Wear a grounded, static-discharging wrist strap.



**Figure 8-3**  
PMT Cover

---

## Calibrating the Nephelometer

Calibrate the Nephelometer before an intensive experiment, calibrate periodically to verify no drift has occurred, calibrate if the reference chopper is dirty or scratched, or if you clean the chopper as part of periodic maintenance.

Perform calibration using TSI software. The "Calibration" option under the Main Menu screen allows you to easily calibrate the Nephelometer using two span gases, comparing the results between air (low span) and CO<sub>2</sub> (high span).

Refer to the "Performing Calibration" section in Chapter 4. Hardware setup instructions are included as part of the calibration process.


---

## Cleaning Fan Filters

The Nephelometer circuitry is cooled by a fan, which has two fan filters located on the outside of the Nephelometer (Figure 8-4). Running the fan with clogged filters can significantly reduce the life of the lamp.

If power is applied to the Nephelometer and the lamp is working properly, but you cannot see light (at any angle) coming through the filters, the foam inserts are clogged and should be cleaned or replaced.

**Note:** You can clean fan filters while power is applied to the Nephelometer.



**Figure 8-4**  
Fan Covers and Foam Inserts

### Tools and Parts

You need a vacuum or a compressed air source, and if you replace the foam inserts, TSI P/N 1602071. One set of foam inserts is included in the accessory kit.

1. Snap off the fan covers and remove the open-cell foam inserts from the covers (Figure 8-4).
2. To clean the foam inserts, use a vacuum to remove accumulated dirt or a compressed air source to blow away the dirt.

**Note:** Cleaning solvents will dissolve the foam inserts.

3. Replace the foam inserts in the fan covers and snap the fan covers back in place.

---

# Replacing the Main Microprocessor EPROM

The Nephelometer sensor contains a 68HC16Z1 main microprocessor that controls signal processing and input/output operations. The main microprocessor has separate battery-backed-up RAM-storage and EPROM chip.

Over time, TSI may update the firmware to enhance the performance of the Nephelometer. Updating the main microprocessor firmware requires replacing a separate EPROM chip on the digital PC board.

**Note:** *Replacing an electronic chip should only be performed by someone who is technically qualified and who is familiar with the operation of the Nephelometer.*

## **Tools and Parts**

You need a Phillips-head screwdriver, a replacement chip, a small flat screwdriver, a static-discharging wriststrap, and an antistatic mat.

To replace the main microprocessor EPROM:

1. Apply power to the Nephelometer and start the Nephelometer software (Chapter 4).
2. Select "Configuration Menu" from the Main Menu.
3. Select "Read Config. Data from Nephelometer" from the Configuration Menu screen.
4. After data is read from the Nephelometer, select "Save Config. Data to Disk File" from the Configuration Menu screen.
5. After data is saved, select "Print Config. Data" from the Configuration Menu screen to make a hard copy of the data record.
6. Confirm that data was written and saved to the NEPHCNFG.DAT file and exit the software.
7. Remove power from the Nephelometer and from any instruments connected to the Nephelometer.


8. Loosen the four screws fastening the Nephelometer top cover and remove the cover (Figure 8-5).



### Caution

To avoid damage to the electronic circuitry, use electrostatic discharge (ESD) precautions when removing the Nephelometer cover:

- Use only a table top with a grounded conducting surface.
- Wear a grounded, static-discharging wrist strap.



**Figure 8-5**  
Top Cover Screws

9. Locate the main microprocessor EPROM on the digital circuit board (Figure 8-6). Remove the main microprocessor EPROM by inserting the tip of a small flat screwdriver under the short side of the chip and *gently* twisting.
10. Remove the new chip from the static protective wrapper.
11. Align the new chip in the socket, making sure notch in the microprocessor EPROM matches the position in Figure 8-6. Apply pressure to the center of the chip until it snaps down into the socket. Press *firmly* to make sure it is seated properly and **check to see that there are no bent or unseated pins.**

- 12.** Apply power to the Nephelometer and observe the LEDs: POWER lights when power is applied, VALVE lights as the valve turns, and STATUS lights when the Nephelometer is ready for operation.

**Note:** *If the LEDs do not function, remove power from the Nephelometer and make sure the chip is seated correctly.*

- 13.** Remove power from the Nephelometer, and replace the top cover, attaching the cover with the four screws.
- 14.** Apply power to the Nephelometer and start the Nephelometer software.
- 15.** Select "Configuration Menu" from the Main Menu.
- 16.** Select "Read Config. Data from Disk File" from the Configuration Menu screen.
- 17.** After data is read from the Nephelometer, select "Write Config Data to Nephelometer" from the Configuration Menu screen.
- 18.** After data is written, press <Esc> to go to the Main Menu and choose "Data Collection" to verify that the Nephelometer is operating properly.

---

# Replacing the Motor Control Microprocessor

The Nephelometer sensor contains an HC711D3 microprocessor that controls the motor functions of the calibrator and backscatter shutters. Over time, TSI may update the firmware to enhance the performance of the Nephelometer. Updating the firmware may require replacing the motor microprocessor chip on the digital PC board.

**Note:** *Replacing an electronic chip should only be performed by someone who is technically qualified and who is familiar with the operation of the Nephelometer.*

## Tools and Parts

You need a Phillips-head screwdriver, a replacement chip, a 44-pin chip extractor (provided with the replacement chip), a static-discharging wriststrap, and an antistatic mat.

To replace the microprocessor:

1. Remove power from the Nephelometer and from any instruments connected to the Nephelometer.
2. Loosen the four screws fastening the Nephelometer top cover and remove the cover (Figure 8-5).




## Caution

To avoid damage to the electronic circuitry, use electrostatic discharge (ESD) precautions when removing the Nephelometer cover:

- Use only a table top with a grounded conducting surface.
- Wear a grounded, static-discharging wrist strap.

3. Locate the motor microprocessor chip on the digital circuit board (Figure 8-6). Use the 44-pin chip extractor to remove the chip.
4. Remove the new chip from the static protective wrapper.

- Align the new chip in the socket, making sure the clipped corner of the motor microprocessor matches the position in Figure 8-6. Apply pressure to the center of the chip until it snaps down into the socket. Press *firmly* to make sure it is seated properly and check to see that there are no bent or unseated pins.



**Figure 8-6**  
Microprocessor and EPROM Locations on the Digital Circuit Board

- Apply power to the Nephelometer and observe the LEDs: POWER lights when power is applied, VALVE lights as the valve turns, and STATUS lights when the Nephelometer is ready for operation.

**Note:** *If the LEDs do not function, remove power from the Nephelometer and make sure the chip is seated correctly.*

- Remove power from the Nephelometer, and replace the top cover, securing the cover with the four screws.




## Replacing the Lamp

The lamp is a projector-type halogen bulb with a built-in elliptical dichroic mirror. The lamp is rated at 75 watts with a recommended maximum of 12 volts. The lamp has an estimated life of 2000 hours of continuous operation at the maximum power (75 watts).

### Tools and Parts

You need a Phillips-head screwdriver and a replacement lamp (P/N 2201111). A replacement lamp is included in the accessory kit.

1. Remove power from the Nephelometer and from any instruments connected to the Nephelometer.
2. Loosen the four screws fastening the Nephelometer top cover and remove the cover (Figure 8-7).




**Figure 8-7**  
Top Cover Screws

3. Unscrew the two screws attaching the lamp shield. Remove the shield, supporting the shield so that it does not hang by the sensor wires (Figure 8-8).



### Caution

Lamp and shield at high temperatures, which can cause burns.  
To avoid personal injury, disconnect power to the Nephelometer and allow the halogen lamp and the lamp shield to cool before handling.



**Figure 8-8**  
Lamp Shield

- 4.** Allow the lamp to cool and push up the lamp lever to force the lamp from the socket.
- 5.** Plug in a new lamp and push down the lamp lever to secure the lamp.
- 6.** Replace the lamp shield and attach the shield with the two screws.
- 7.** Replace the Nephelometer top cover and attach with the four screws.
- 8.** After applying power to the Nephelometer, check the operation of the lamp either by observing the light inside the sensor or by selecting “Data Collection” from the Main Menu of the TSI software.

---

## Replacing Aerosol Filters

The Nephelometer has two aerosol filters:

- ❑ A large white HEPA filter that can be switched into the inlet to filter all sample air coming into the sensor.
- ❑ A small blue DQ filter that is used to purge a small flow of clean air into the light pipe apertures. This airflow keeps the lens clean.

There is no software indication or error message if either filter is not functioning. However, if there is a significant increase in pressure drop through the filter from the time you set up Nephelometer to the present, check the white HEPA filter and fittings to see if the inlet flow is blocked.


Replace the white HEPA filter every six months of operation and replace the blue DQ filter every 12 months. If the Nephelometer is sampling very dirty aerosols, shorten these times. If the Nephelometer is sampling very clean aerosols, extend the times.

### **Tools and Parts**

You need a Phillips-head screwdriver, adjustable wrench, small needle-nose pliers, RTV silicone sealer, and replacement filters, white HEPA: TSI P/N 1602051, and blue DQ: TSI P/N 1602080. Replacement filters are included in the accessory kit.

To replace the white HEPA filter or the blue DQ filter:


- 1.** Remove power from the Nephelometer and from any instruments connected to the Nephelometer.
- 2.** Loosen the four screws fastening the top cover of the Nephelometer and remove the cover (Figure 8-9).



**Figure 8-9**  
Top Cover Screws


**3. To replace the white HEPA filter:**

- a.** Loosen the end plate by removing the two bottom cover screws (Figure 8-10).
- b.** Unscrew the three screws fastening the end plate and remove the end plate.




**Figure 8-10**  
Bottom Cover Screws Securing End Plate

- c.** Use an adjustable wrench to remove the nylon nut connecting the elbow fitting to the filter outlet (Figure 8-11).
- d.** Pull back on the elbow to separate the elbow fitting from the filter.
- e.** Unscrew the filter from the brass fitting on the inlet block (Figure 8-11).
- f.** Coat the threads on the brass fitting on the inlet block with RTV silicone sealer.
- g.** Making sure the arrow on the new filter is pointing away from the inlet block, screw the new filter into the brass fitting on the inlet block (Figure 8-11).
- h.** Remove the nylon fitting from the old filter, coat the threads of the fitting with RTV silicone sealer, and screw the nylon fitting into the outlet of the new filter.
- i.** Use the nylon nut to attach the elbow fitting to the outlet of the filter.
- j.** Reattach the end plate using the three screws.
- k.** Replace the two bottom cover screws.



**Figure 8-11**  
White HEPA Filter

- 4.** To replace the blue DQ filter:
  - a.** Remove the  $\frac{1}{8}$  in. ID silicone tubing from both ends of the blue filter (Figure 8-12).
  - b.** Making sure the arrow of the new filter is pointing away from the inlet block, install the new DQ filter, reattaching the silicone tubing.



**Figure 8-12**  
Blue DQ Filter

- 5.** Replace the top cover of the Nephelometer and attach with the four screws.

---


## Cleaning the Light Pipe Lens

The light pipe is a solid glass rod that transfers light from the lamp to the measurement volume, providing a thermal break between the lamp and the measurement volume. The light pipe lens, which receives a high volume of cooling air, may require frequent cleaning.

### Tools


You need a Phillips-head screwdriver, soft cloth, cotton swabs, isopropyl alcohol, and O-ring grease.

1. Remove power from the Nephelometer and from any instruments connected to the Nephelometer.
2. Loosen the four screws fastening the top cover of the Nephelometer and remove the cover (Figure 8-13).



**Figure 8-13**  
Top Cover Screws

3. Unscrew the two screws attaching the lamp shield (Figure 8-14). Remove the shield, supporting the shield so that it does not hang by the sensor wires.



**Figure 8-14**  
Lamp Shield

4. Allow the lamp to cool and push up the lamp lever to force the lamp from the socket.




### **C a u t i o n**

Lamp and shield at high temperatures, which can cause burns.  
To avoid personal injury, disconnect power to the Nephelometer and allow the halogen lamp and the lamp housing to cool before handling.



5. Locate the light pipe lens below the lamp socket. Use the cotton swab and alcohol to clean the lens (Figure 8-15).



**Figure 8-15**  
Light Pipe Lens

6. Plug in the lamp and push down the lamp lever to secure the lamp.
7. Replace the lamp shield securing it with the two screws.
8. Replace the top cover of the Nephelometer and attach with the four screws.
9. After applying power to the Nephelometer, check the operation of the lamp either by observing the light inside the sensor or by typing the RF command on the computer and checking the status flags (see "Read Commands" in Chapter 6).

---

## Checking for Leaks

This section describes procedures for checking the vacuum integrity of the Nephelometer.

### **Tools needed**

A vacuum pump (capable of 18 in. Hg vacuum), manometer or vacuum gauge, leak check solution (bubble solution) pump capable of 5 psi pressure, 1" NPT plugs.

To check if the instrument has a leak:

1. Remove the blower from the instrument and plug the lower connection port.
2. Connect a vacuum pump and vacuum gauge to the inlet.
3. Draw the pressure down to about 18 in. Hg below atmospheric pressure.
4. Wait for about ten minutes for any pressure/temperature changes with the instrument.
5. Record the reading on the vacuum gauge.
6. Wait five minutes and record again.
7. If the gauge has dropped by more than 0.2 in. Hg in five minutes, there is probably a leak.

To check for leaks:

1. Make sure that the vacuum system itself does not leak.
2. Remove the top and bottom covers from the instrument.
3. Connect the pump in place of the vacuum pump and pressurize the instrument to **no more than 5 psi**.
4. Drip bubble solution over all silicone sealed joints or O-ring joints until the leak is found.
5. Repair/replace seals or O-ring as needed.
6. Retest with vacuum as described above.
7. Reassemble covers.

---

## Cleaning the Reference Chopper

The reference chopper contains three areas: one area provides a signal from light scattered by an aerosol, one area provides an indication of the lamp power, and one area provides a measure of the PMT dark current.


Make sure the chopper is operating properly and clean the calibrate portion of the chopper if, over time, the calibrate signal from the Nephelometer rises (from light scattered by dirt) or falls (from light blocked by dirt) significantly using the same lamp power.

**Note:** *Recalibrate the Nephelometer after you clean the reference chopper. Refer to Chapter 4, "Using Nephelometer Software (DOS)," for the calibration procedure using TSI software commands.*

### **Tools and Parts**

You need a Phillips-head screwdriver, isopropyl alcohol, and cotton swabs.

- 1.** Remove power from the Nephelometer and from any instruments connected to the Nephelometer.
- 2.** Loosen the four screws fastening the top cover of the Nephelometer and remove the cover (Figure 8-16).




**Figure 8-16**  
Top Cover Screws, Bottom Cover Screws, Fan Cover Screws, and I/O Panel Screws

**3.** The bottom cover of the Nephelometer is pliable and fits very tightly around the bottom of the instrument. Remove the bottom cover in one of two ways:

- ❑ The more difficult method is to pry off the bottom cover without removing the I/O panel and the fan panel. Unscrew the four bottom cover screws (Figure 8-16). *Carefully* remove the bottom cover, working alternately at both ends. Pry the edges of the cover apart, and ease the cover away from the body of the Nephelometer.
- ❑ The easier method takes more time. After removing the I/O panel, the fan panel and the end plate, slide the bottom cover off the body of the instrument. Using Figure 8-16 as a reference:
  - a.** Unscrew the three screws securing the I/O panel to the backbone, loosen the panel, and let the panel hang over the side of the instrument.
  - b.** Unscrew the two screws securing the fan panel to the backbone, loosen the panel, and let the panel hang over the side of the instrument.






**Figure 8-17a**  
PMT End Plate Screws

- 8.** Using one hand to hold the PMT box, use the other hand to slide the aperture assembly from the outlet section of the Nephelometer.


**Note:** *Be careful not to let the inner tube drag inside the outer tube and scrape the paint.*

- 9.** Set the aperture assembly on a flat surface, PMT end down (Figure 8-18).




**Figure 8-18**  
Aperture Assembly

- 10.** Using an adjustable wrench, remove the two nuts and washers attaching the plate to the aperture assembly (Figure 8-18).
- 11.** The aperture plate fits snugly. As you remove the plate, disconnect the motor connector (J13) from the circuit board.
- 12.** To access the reference chopper:
  - a.** Set the plate assembly, circuit board side up, on a flat surface (Figure 8-19).
  - b.** Remove the three screws and washers attaching the circuit board and plate 1 to plate 2 (Figure 8-19).
  - c.** Separate plate 1 and circuit board from plate 2 and turn plate 1 so that the shutter is facing you (Figure 8-20).



**Figure 8-19**  
Aperture Plates and Circuit Board



**Figure 8-20**  
Reference Chopper



13. Use a cotton swab and alcohol to clean *only* the frosted or reflective aluminum surfaces of the reference chopper.

**Note:** *Although TSI recommends isopropyl alcohol to clean the chopper, alcohol will remove black paint from components.*

14. **Optional**

If you need better access to the frosted surface of the reference chopper, loosen the screw in the counterweight and remove the chopper from plate 1.

**Note:** *If you remove the reference chopper from plate 1, reattach the chopper by placing two sheets of paper as a spacer between the chopper and plate 1 so the hub of the chopper is very close to plate 1 without touching. Then tighten the setscrew and remove the paper. Make sure the chopper turns freely.*

15. Use the three screws and washers to reattach the circuit board and plate 1 to plate 2. Using a cotton swab, turn the chopper by the hub or plate to make sure that the chopper does not come in contact with any plate surface.
16. Reconnect the motor connector (J13) to the circuit board and seat the plate assembly in the aperture assembly.  
**Note:** *Make sure the motor connector cable does not block an aperture inside the aperture assembly.*
17. Reattach the two nuts and washers securing the aperture plate to the aperture assembly.
18. Reinsert the aperture assembly into the outlet end of the Nephelometer being careful not to bend the flocked paper insert in the middle of the Nephelometer sensor.
19. Reattach the PMT end plate with the three screws.
20. Reattach the analog and digital circuit boards to the end plate using the two screws (Figure 8-17).
21. Attach the bottom cover of the Nephelometer reversing the steps in step 3 of this procedure. Remember that the cover is pliable and fits tightly around the bottom of the instrument.

- 22.** Replace the top cover of the Nephelometer and attach with the four screws.
- 23.** Recalibrate the Nephelometer after you clean the reference chopper. Refer to Chapter 4, "Using Nephelometer Software (DOS)," for the calibration procedure using TSI software commands.

---

## Cleaning the Light Pipe Outlet and the Backscatter Shutter

The light pipe is a solid glass rod that transfers light from the lamp to the measurement volume, providing a thermal break between the lamp and the measurement volume. The backscatter shutter, easily accessed while you are cleaning the outlet to the light pipe, collects dust when it is spinning in backscatter mode.


### **Tools**

You need a Phillips-head screwdriver, soft cloth, cotton swabs, isopropyl alcohol, and O-ring grease.

1. If you are going to clean the outlet of the light pipe and the backscatter shutter, park the shutter in "Sleep" mode using the PD software command (refer to the "Action Commands" in Chapter 6).


**Note:** *The PD command positions the backscatter shutter so that you can remove the lamp base assembly from the Nephelometer without damaging the shutter.*

2. Remove power from the Nephelometer and from any instruments connected to the Nephelometer.
3. Loosen the four screws fastening the top cover of the Nephelometer and remove the cover (Figure 8-21).



**Figure 8-21**  
 Top Cover Screws, Bottom Cover Screws, Fan Cover Screws, and I/O Panel Screws


- 4.** Unscrew the two screws attaching the lamp shield (Figure 8-22). Remove the shield, supporting the shield so that it does not hang by the sensor wires.



**Figure 8-22**  
 Lamp Shield


- 5.** The bottom cover of the Nephelometer is pliable and fits very tightly around the bottom of the instrument. Remove the bottom cover in one of two ways:
- The more difficult method is to pry off the bottom cover without removing the I/O panel and the fan panel. Unscrew the four bottom cover screws (Figure 8-21). *Carefully* remove the bottom cover, working alternately at both ends. Pry the edges of the cover apart, and ease the cover away from the body of the Nephelometer.
  - The easier method takes more time. After removing the I/O panel, the fan panel and the end plate, slide the bottom cover off the body of the instrument. Using Figure 8-21 as a reference:
    - a.** Unscrew the three screws securing the I/O panel to the backbone, loosen the panel, and let the panel hang over the side of the instrument.
    - b.** Unscrew the two screws securing the fan panel to the backbone, loosen the panel, and let the panel hang over the side of the instrument.
    - c.** Unscrew the three screws securing the end plate to the inlet section of the instrument and remove the end plate.
    - d.** Unscrew the four bottom cover screws, and *carefully* slide the cover off the end of the instrument.

- 6.** Using Figure 8-23 as a reference:
- a.** Remove the connector for the lamp on the analog board.
  - b.** Remove the ribbon cable connector from the inlet block.
  - c.** Remove the motor controller ribbon cable connector from the lamp block.
  - d.** Disconnect these sensor connectors (as marked): 1 (lamp temperature), 2 (valve position), 3 (valve motor), 4 (heater), and 6 (fan).
  - e.** Disconnect the blue DQ filter silicone tubing at the elbow connector on the backbone and at the elbow connector on the inlet.




**Figure 8-23**  
Disconnect Connectors and Tubing

- 7.** To separate the Nephelometer inlet section from the outlet section, loosen the three screws attaching the two sections to each other (Figure 8-24) and *carefully* pull the sections apart (Figure 8-25). The flocked paper insert will be visible.




**Figure 8-24**  
Nephelometer Sections Together



**Figure 8-25**  
Nephelometer Sections Apart

- 8.** To clean the light pipe outlet:
  - a.** Unscrew the four screws, one on each corner, on the aluminum lamp base. Two of the screws are recessed (Figure 8-26).
  - b.** Check to make sure the backscatter shutter is located underneath the light pipe. If not, rotate the shutter with your fingers.




**Figure 8-26**  
Lamp Base



- c. Using one hand to support the backscatter shutter from below, use the other hand to *carefully* remove the lamp base assembly, angling the assembly as you pull it out of the sensor body (Figure 8-27).

**Note:** *Be careful not to catch the shadow lightstop plate on the backbone.*




**Figure 8-27**  
Removing the Lamp Base Assembly From the Sensor

- d. Locate the light pipe outlet on the underside of the lamp base assembly (Figure 8-28). Use a cotton swab and alcohol to clean the outlet, making sure the swab does *not* touch black components.


**Note:** Although TSI recommends isopropyl alcohol to clean the outlet, alcohol removes black paint from components.

- e. Use a soft dry cloth to clean the backscatter shutter. If necessary, use a soft cloth and a mild detergent solution followed by a soft dry cloth.



**Figure 8-28**  
Light Pipe Outlet and Backscatter Shutter

- f. Grease the O-ring on the backbone of the sensor (Figure 8-29) and replace the O-ring.
- g. Rotate the backscatter shutter so it is positioned under the light pipe and *gently* reinsert the lamp base assembly, at an angle, into the sensor.
- h. Tighten the four screws, one at each corner of the aluminum lamp base.



**Figure 8-29**  
O-ring on the Sensor Backbone

- 9.** Put a light coating of grease on the O-ring between the Nephelometer inlet and outlet sections and replace the O-ring in the groove (Figure 8-25).
- 10.** Use the three screws to attach the Nephelometer outlet to the inlet.

- 11.** Using Figure 8-23 as a reference:
  - a.** Attach the connector for the lamp to the analog board.
  - b.** Attach the ribbon cable connector to the inlet block.
  - c.** Attach the motor controller ribbon cable connector to the lamp block.
  - d.** Connect these sensor connectors (as marked): 1 (lamp temperature), 2 (valve position), 3 (valve motor), 4 (heater), and 6 (fan).
  - e.** Connect the blue DQ filter silicone tubing at the elbow connector on the backbone and at the elbow connector on the inlet.
- 12.** Attach the bottom cover of the Nephelometer reversing the steps in step 5 of this procedure. Remember that the cover is pliable and fits tightly around the bottom of the instrument.
- 13.** Replace the top cover of the Nephelometer and attach with the four screws.

---


## Cleaning the Flocked Paper

Black flocked paper inside the main tube has a direct effect on the amount of wallscatter light. If the paper collects dust, more light is scattered, affecting counts from the wall.

### Tools and Parts

You need a Phillips-head screwdriver, a vacuum or a compressed air source, and O-ring grease.


1. Remove power from the Nephelometer and from any instruments connected to the Nephelometer.
2. Loosen the four screws fastening the top cover of the Nephelometer and remove the cover (Figure 8-30).



**Figure 8-30**


Top Cover Screws, Bottom Cover Screws, Fan Cover Screws, and I/O Panel Screws

- 3.** The bottom cover of the Nephelometer is pliable and fits very tightly around the bottom of the instrument. Remove the bottom cover in one of two ways:
  - The more difficult method is to pry off the bottom cover without removing the I/O panel and the fan panel. Unscrew the four bottom cover screws (Figure 8-30). *Carefully* remove the bottom cover, working alternately at both ends. Pry the edges of the cover apart, and ease the cover away from the body of the Nephelometer.
  - The easier method takes more time. After removing the I/O panel, the fan panel and the end plate, slide the bottom cover off the body of the instrument. Using Figure 8-30 as a reference:
    - a.** Unscrew the three screws securing the I/O panel to the backbone, loosen the panel, and let the panel hang over the side of the instrument.
    - b.** Unscrew the two screws securing the fan panel to the backbone, loosen the panel, and let the panel hang over the side of the instrument.
    - c.** Unscrew the three screws securing the end plate to the inlet section of the instrument and remove the end plate.
    - d.** Unscrew the four bottom cover screws, and *carefully* slide the cover off the end of the instrument.
- 4.** If the Nephelometer is resting on its feet, support the underside of the instrument.
- 5.** Using Figure 8-31 as a reference:
  - a.** Remove the connector for the lamp on the analog board.
  - b.** Remove the ribbon cable connector from the inlet block.
  - c.** Remove the motor controller ribbon cable connector from the lamp block.
  - d.** Disconnect these sensor connectors (as marked): 1 (lamp temperature), 2 (valve position), 3 (valve motor), 4 (heater), and 6 (fan).
  - e.** Disconnect the blue DQ filter silicone tubing at the elbow connector on the backbone and at the elbow connector on the inlet.




**Figure 8-31**  
Disconnect Connectors and Tubing

- 6.** To separate the Nephelometer inlet section from the outlet section, loosen the three screws attaching the two sections to each other (Figure 8-32) and *carefully* pull the sections apart (Figure 8-33). The flocked paper will be visible.




**Figure 8-32**  
Nephelometer Sections Together



**Figure 8-33**  
Nephelometer Sections Apart



7. Pull the flocked paper from the outlet section (Figure 8-34).



**Figure 8-34**  
Flocked Paper Outside the Nephelometer

8. To clean the flocked paper, use a vacuum to take up or a compressed air source to blow away the dirt.

**Note:** *If the paper has bald patches, it should be replaced. Contact TSI for assistance.*

9. Put a light coating of grease on the O-ring between the Nephelometer inlet and outlet sections and replace the O-ring in the groove (Figure 8-33).

10. *Carefully* reinstall the flocked paper, fitting the paper around the PMT aperture plate in the outlet section of the Nephelometer.

11. Use the three screws to attach the Nephelometer outlet section to the inlet section.

12. Using Figure 8-31 as a reference:

a. Attach the connector for the lamp to the analog board.

b. Attach the ribbon cable connector to the inlet block.

c. Attach the motor controller ribbon cable connector to the lamp block.

- d.** Connect these sensor connectors (as marked): 1 (lamp temperature), 2 (valve position), 3 (valve motor), 4 (heater), and 6 (fan).
  - e.** Connect the blue DQ filter silicone tubing at the elbow connector on the backbone and at the elbow connector on the inlet.
- 13.** Attach the bottom cover of the Nephelometer reversing the steps in step 3 of this procedure. Remember that the cover is pliable and fits tightly around the bottom of the instrument.
  - 14.** Replace the top cover of the Nephelometer and attach with the four screws.

---


# Replacing a Temperature or Humidity Sensor

The Nephelometer has two built-in temperature sensors and a humidity sensor. The sensors are mounted on PC boards to make service easy. Although the sensors are rugged, they may have to be replaced sometime during the life of the instrument. The first temperature sensor is located at the inlet of the instrument. The second temperature sensor and the humidity sensor are located near the outlet, on the same PC board.

### Tools and Parts

You need a Phillips-head screwdriver, a replacement sensor mounted on its PC board, and O-ring grease.


1. Remove power from the Nephelometer and from any instruments connected to the Nephelometer.
2. Loosen the four screws fastening the top cover of the Nephelometer and remove the cover (Figure 8-35).



**Figure 8-35**

Top Cover Screws, Bottom Cover Screws, Fan Cover Screws, and I/O Panel Screws

- 3.** Disconnect the ribbon cable from the appropriate sensor (Figure 8-36).
- 4.** Remove the four screws from the appropriate PC board and remove the sensor assembly.



**Figure 8-36**  
Humidity and Temperature Sensor Locations

- 5.** Lightly grease the sealing O-ring if necessary.
- 6.** Replace the sensor assembly with the replacement sensor.
- 7.** Replace the four screws.
- 8.** Replace the cover and run the Nephelometer software to verify that the sensor is operating correctly.
- 9.** Calibrate the sensor, if needed (refer to SC command in Chapter 6).

---

# Replacing a Photomultiplier Tube


As a photomultiplier tube (PMT) ages, it becomes less sensitive. You may find that you can compensate for the gain by increasing the voltage. At this time, there is no data on life expectancy of a PMT or how the tube may be affected by different environments. If there is a substantial increase in dark current (see the RP command under "Read Commands" in Chapter 6), replace the PMT.

**Notes:** *Replacing a PMT should only be attempted by someone who is technically qualified and who is familiar with the operation of the Nephelometer.*

## Tools and Parts

You need a Phillips-head screwdriver and a replacement PMT:

- Blue PMT R1527P TSI P/N 3009004
- Green PMT R1527P TSI P/N 3009004
- Red PMT R2949 TSI P/N 3009005



## WARNING

High voltages that can shock or burn.

The box that houses the photomultiplier tubes (PMT) contains a high voltage source. Components within this box should only be serviced by a qualified technician. Disconnect all power to the Nephelometer and contact TSI personnel before disassembling the PMT box.




## Caution

Light sensitive components.

To avoid damage to the photomultiplier tubes (PMT), remove power from the Nephelometer and open the PMT box in subdued lighting. Exposing the photomultiplier tubes to direct light may affect tube performance.

1. Remove power from the Nephelometer and from any instruments connected to the Nephelometer.
2. Unscrew the four screws attaching the PMT cover to the sensor and remove the cover (Figure 8-37)




**Figure 8-37**  
PMT Cover

- 3.** Each PMT fits inside a magnetic shield (Figure 8-38). To remove the shield and PMT, support the sensor frame with one hand and use your other hand to *gently* rock the shield from side to side as you pull it from its socket.
- 5.** To remove the PMT from the shield, loosen the three setscrews at the base of the shield (Figure 8-38) and withdraw the tube.
- 6.** *Carefully* insert a new PMT into the shield. Position the tube so that the light-sensing portion (criss-crossed wires) of the tube faces the slotted opening on the shield. *Gently* tighten the three setscrews so that the end of the tube is  $\frac{1}{8}$  in. from the end of the shield.

**Notes:** *If the tube is not inserted far enough, the setscrews will touch glass and may shatter the tube when they are tightened.*

*Tighten the setscrews only enough to support the tube. If you overtighten the screws you may crack the plastic base.*



**Figure 8-38**  
Red, Blue and Green PMTs with Setscrews

- 7.** Plug the tube into the socket, noting the alignment of the center pin of the tube into its mating socket.
- 8.** Replace the PMT cover and attach with the four screws.

---


## Checking, Cleaning, and Replacing Bandpass Filters

Three bandpass filters are centered at 450 nm (blue), 550 nm (green), and 700 nm (red) wavelengths. Although located within a protective can, the filters may need periodic cleaning if they accumulate dust or dirt. In addition, exposure to high humidity can attack the coatings on the filters causing them to cloud.

### Tools and Parts

You need a Phillips-head screwdriver, an S-shaped Phillips-head screw driver, and a bandpass filter:

- Blue filter P/N 2502486
- Green filter P/N 2502485
- Red filter P/N 2502487



### WARNING

High voltages that can shock or burn.

The box that houses the photomultiplier tubes (PMT) contains a high voltage source. Components within this box should only be serviced by a qualified technician. Disconnect all power to the Nephelometer and contact TSI personnel before disassembling the PMT box.




### Caution

Light sensitive components.

To avoid damage to the photomultiplier tubes (PMT), remove power from the Nephelometer and open the PMT box in subdued lighting. Exposing the photomultiplier tubes to direct light may affect tube performance.


1. Remove power from the Nephelometer and from any instruments connected to the Nephelometer.
2. Remove the PMT cover by unscrewing the four screws attaching the cover to the sensor (Figure 8-39).






**Figure 8-39**  
PMT Cover

- 3.** To check the bandpass filters, you must remove the PMTs with their magnetic shields (Figure 8-40). To remove a shield and tube, support the sensor frame with one hand and use your other hand to *gently* rock the shield from side to side as you pull it from its socket.
- 4.** Visually inspect the filter but *avoid* touching the filter with your fingers.
  - If there is dust, use a lens cloth to remove the dust. If the filter looks clean and transparent, it does *not* need to be replaced. Go to step 7 of this procedure.
  - If the filter is cloudy, remove the filter for a closer inspection. Go to step 5 of this procedure.



**Figure 8-40**  
Red, Blue and Green PMTs with Setscrews

- 5.** Use the S-shaped Phillips-head screwdriver to unscrew the two screws holding the filter in place, and remove the filter (Figure 8-41).
- 6.** Hold the filter up to the light. If the filter is *not* transparent, or if there is clouding around the edges, replace the filter. Attach the new filter in the holder using the two screws and washers.
- 7.** Plug the tube into the socket, noting the alignment of the center pin of the tube into its mating socket.
- 8.** Replace the PMT cover and attach with the four screws.



**Figure 8-41**  
Bandpass Filter

---


## Cleaning the Lens

The lens collimates light from apertures that define the viewing volume in the body of the instrument. The light is collimated before it is split into separate colors to ensure that each photomultiplier tube sees the same area of light scatter from the viewing volume.

**Note:** *Under normal circumstances, you do not have to periodically check or clean the lens. Check the lens if there is an overall reduction in the signal in all three wavelengths over time.*

### Tools and Parts

You need a Phillips-head screwdriver, flashlight, lens cloth, and cotton swabs.



### WARNING

High voltages that can shock or burn.

The box that houses the photomultiplier tubes (PMT) contains a high voltage source. Components within this box should only be serviced by a qualified technician. Disconnect all power to the Nephelometer and contact TSI personnel before disassembling the PMT box.




### Caution

Light sensitive components.

To avoid damage to the photomultiplier tubes (PMT), remove power from the Nephelometer and open the PMT box in subdued lighting. Exposing the photomultiplier tubes to direct light may affect tube performance.

1. Remove power from the Nephelometer and from any instruments connected to the Nephelometer.
2. Remove the PMT cover by unscrewing the four screws attaching the cover to the sensor (Figure 8-42).



**Figure 8-42**  
PMT Cover


- 3.** To access the lens, you must remove the red PMT tube with its magnetic shield (Figure 8-43). To remove the tube, support the sensor frame with one hand and use your other hand to *gently* rock the shield from side to side as you pull the shield and tube from its socket.

**Note:** *If you need further access to the lens, you can remove the screw attaching the lens holder.*

- 4.** Use a flashlight to inspect the lens for signs of dust.
- 5.** If the lens is dusty, *carefully* clean the lens with a cotton swab you have covered with lens paper.

**Note:** *Although TSI normally recommends isopropyl alcohol to clean lenses, alcohol or acetone removes black paint from components.*

- 6.** Replace the tube into its socket, noting the alignment of the center pin of the tube with its mating socket.
- 7.** Replace the PMT cover and attach with the four screws.



**Figure 8-43**  
Red, Blue and Green PMTs with Setscrews

## CHAPTER 9

# Contacting Customer Service

This chapter gives directions for contacting people at TSI Incorporated for technical information and directions for returning the Model 3550/3560 Series Integrating Nephelometer for service.

---

## Technical Contacts at TSI

- ❑ If you have any difficulty installing the Nephelometer, or if you have technical or application questions about this instrument, contact an applications engineer at TSI Incorporated, 1-800-678-2708 or (651) 490-2833.
- ❑ If the Nephelometer fails, or if you are returning it for service, contact TSI Particle Instrument Division, Customer Service, at 1-800-874-3893 (USA) or (651)-490-3893.

---

## Returning the Nephelometer for Service

Call Customer Service for specific return instructions. Customer Service will need this information when you call:

- ❑ The instrument model number
- ❑ The instrument serial number
- ❑ A purchase order number (unless under warranty)
- ❑ A billing address
- ❑ A shipping address.

Use the original packing material to return the instrument to TSI. If you no longer have the original packing material, seal off the inlet and fan on the sensor to prevent debris from entering the instrument and ensure that the indicator lights and the connectors on the instrument back panel are protected.





# APPENDIX A

# Specifications

Table A-1 lists the specifications for the Nephelometer sensor and Table A-2 lists the power supply. These specifications are subject to change.

**Table A-1**  
Specifications for the Nephelometer Sensor

<b>Wavelengths</b>	
Models 3551.....	550 nm (green)
Models 3563.....	450 nm (blue), 550 nm (green), and 700 nm (red/infrared)
<b>Bandwidth</b> .....	50 nm (all wavelengths)
<b>Sensitivity at 30-sec averaging time (aerosol scattering coefficient, <math>\sigma_{sp}</math>)</b>	
Blue and green wavelengths .....	$1.0 \times 10^{-7} \text{ m}^{-1}$
Red/Infrared wavelength .....	$3.0 \times 10^{-7} \text{ m}^{-1}$
<b>Averaging time</b> .....	1 to 4096 sec (selectable)
<b>Drift</b> .....	$<2.0 \times 10^{-7} \text{ m}^{-1}$ at 30-sec averaging time for up to one hour after filtered-air reference measurement for green wavelength
<b>Optical background signal</b>	
Blue and green wavelengths .....	$<5.0 \times 10^{-5} \text{ m}^{-1}$
Red/Infrared wavelength .....	$<1.0 \times 10^{-5} \text{ m}^{-1}$
<b>Angular integration</b> .....	7 to 170°
<b>Backscatter shutter</b> .....	Changes angular integration from <u>7 to 170°</u> to <u>90 to 170°</u> (Models 3561 and 3563 only)
<b>Reference shutter</b> .....	Allows measurement of light intensity of a reference object illuminated by main lamp or of photodetector's dark signal
<b>Filtered air reference chopper</b> .....	High-efficiency particle filter switches into sample air stream automatically on host computer demand or at intervals selected by user
<b>Response time</b> .....	<10 sec
<b>Recommended flowrate</b> .....	20 to 200 L/min
<b>Inlet/Outlet Dimensions</b> .....	25 mm diameter
<b>Particle transport efficiency</b> .....	>95% of unit-density particles from 0.05 to 5 $\mu\text{m}$ in diameter

*(continued)*

**Table A-1**Specifications for the Nephelometer Sensor (*continued*)

<b>Temperature and pressure sensors...</b>	Built-in sensors allow corrections for changes in the Rayleigh-scattering coefficient of air within sample volume
<b>Humidity sensor.....</b>	Measures relative humidity of sample from 5 to 95% $\pm$ 5%
<b>Time and date .....</b>	Provided by internal, real-time clock with battery backup
<b>Vacuum integrity .....</b>	<10 mm Hg/hr at a negative pressure of 700 mm Hg (not including optional blower)
<b>Dimensions.....</b>	1100 mm $\times$ 300 mm $\times$ 250 mm (43 in. $\times$ 12 in. $\times$ 10 in.)
<b>Weight.....</b>	<18 kg (<40 lb)
<b>Power requirements .....</b>	Operates on 24.0 $\pm$ 4.0 VDC at <5.0 A (125 W maximum), supplied by power supply (included)
<b>Environmental conditions .....</b>	Indoor use Altitude up to 2000 m (6500 ft) Ambient temperature 5–40°C Ambient humidity 0–90% RH noncondensing Overvoltage category II Pollution degree II

**Table A-2**

Specifications for the Power Supply

<b>Output .....</b>	24.0 $\pm$ 4.0 VDC at <5.0 A
<b>Dimensions.....</b>	300 mm $\times$ 150 mm $\times$ 100 mm (12 in. $\times$ 8 in. $\times$ 5 in.)
<b>Weight.....</b>	<5 kg (<11 lb)
<b>Power requirements .....</b>	85–260 VAC, 50–60 Hz at <100 W
<b>Fuse (not replaceable by user) .....</b> <i>(Internal—<b>not</b> accessible by operator)</i>	~F 7A FB/250V

TSI and the TSI logo are registered trademarks of TSI Incorporated.

APPENDIX B

# **Bibliography**

Prepared for  
TSI Incorporated  
by  
Robert J. Charlson, Ph. D.

## **Instrument Principles, Design and Calibration**

- 1.** Beuttell, R.G. and A. W. Brewer, (1949), "Instruments for the Measurement of the Visual Range," *Journal of Scientific Instruments and of Physics in Industry*, 26, 357-359.
- 2.** Ahlquist, N. C. and R. J. Charlson, (1967), "A New Instrument for Evaluating the Visual Quality of Air," *APCA Journal*, 17, 467-469.
- 3.** Ensor, D. S. and A. P. Waggoner, (1970), "Angular Truncation Error in the Integrating Nephelometer," *Atmospheric Environment*, 4, 481-487.
- 4.** Hasan, H. and C. W. Lewis, (1983), "Integrating Nephelometer Response Corrections for Bimodal Size Distributions," *Aerosol Science and Technology*, 2, 443-453.
- 5.** Ruby, M.G. and A. P. Waggoner, (1981), "Intercomparison of Integrating Nephelometer Measurements," *Environmental Science and Technology*, 15, 109-113.
- 6.** Bodhaine, B. A., (1979), "Measurement of the Rayleigh Scattering Properties of Some Gases with a Nephelometer," *Applied Optics*, 18, 121-125.
- 7.** Heintzenberg, J. and L. Bäcklin, (1982), "A High Sensitivity Integrating Nephelometer for Airborne Air Pollution Studies," *Atmospheric Environment*, 17, 433-436.
- 8.** Charlson, R. J., (1980), "The Integrating Nephelometer," *Atmospheric Technology*, Number 12, 10-14.
- 9.** Bodhaine, B. A., N. C. Ahlquist, and R. C. Schnell, (1991), "Three-Wavelength Nephelometer Suitable for Aircraft Measurement of Background Aerosol Scattering Coefficient," *Atmospheric Environment*, 25A, 2267-2276.
- 10.** Heintzenberg, J., (1978), "The Angular Calibration of the Total Scatter/Backscatter Nephelometer, Consequences and Applications," *Staub-Reinhalt Luft*, 38, 62-63.

- 11.** Heintzenberg, J. and P. S. Bhardwaja, (1976), "On the Accuracy of the Backward Hemispheric Integrating Nephelometer," *Journal of Applied Meteorology*, 15, 1092-1096.
- 12.** Middleton, W. E. K., (1952), in *Vision Through the Atmosphere*, pp. 203-206. The University of Toronto Press, Toronto, Canada.
- 13.** Rae, J. B. and J. A. Garland, (1970), "A Stabilized Integrating Nephelometer for Visibility Studies," *Atmospheric Environment*, 4, 219-223.
- 14.** Ruppertsberg, G. H., (1967), "Principles et Procédés de Mesure Automatique de la Visibilité," *International Association of Lighthouse Authorities Bulletin*, No. 31, January.
- 15.** Heintzenberg, J. and H. Quenzel, (1973), "On the Effect of the Loss of Large Particles on the Determination of Scattering Coefficients with Integrating Nephelometers," *Atmospheric Environment*, 7, 503-507.
- 16.** Heintzenberg, J. and H. Quenzel, (1973), "Calculations on the Determination of the Scattering Coefficient of Turbid Air with Integrating Nephelometer," *Atmospheric Environment*, 7, 509-519.
- 17.** Nyeki, S. A. P., I. Colbeck, and R. M. Harrison, (1992), "A Portable Aerosol Sampler to Measure Real-Time Atmospheric Aerosol Mass Loadings," *Journal of Aerosol Science*, 23, S687-S690.
- 18.** Charlson, R. J. and N. C. Ahlquist, (1971), *Integrating Nephelometer*, U.S. Patent No. 3,563,661.
- 19.** Charlson, R. J. and N. C. Ahlquist, (1972), *Method and Apparatus for Making an In-Situ Determination of the Chemical Properties of Atmospheric Aerosols*, U.S. Patent No. 3,700,333.
- 20.** Ahlquist, N. C., A. P. Waggoner, and R. J. Charlson, (1976), *Photon-Counting Integrating Nephelometer*, U.S. Patent No. 3,953,127.

## **Applications of the Integrating Nephelometer to Studies of Atmospheric Visibility**

- 1.** Charlson, R. J., H. Horvath and R. F. Pueschel, (1967) "The Direct Measurement of Atmospheric Light Scattering Coefficient for Studies of Visibility and Pollution," *Atmospheric Environment*, 1, 469-478.
- 2.** Horvath, H. and K. E. Noll, (1969), "The Relationship Between Atmospheric Light Scattering Coefficient and Visibility," *Atmospheric Environment*, 3, 543-550.
- 3.** Waggoner, A. P., R. E. Weiss, N. C. Ahlquist, D. S. Covert, S. Will, and R. J. Charlson, (1981), "Optical Characteristics of Atmospheric Aerosols," *Atmospheric Environment*, 15, 1891-1909.
- 4.** Dzubay, T. G., R. K. Stevens, C. W. Lewis, D. H. Hern, W. J. Courtney, J. W. Tesch, and M. A. Mason, (1982), "Visibility and Aerosol Composition in Houston, Texas," *Environmental Science and Technology*, 16, 514-525.
- 5.** Trijonis, J. C. (Principal author), (1990), "Visibility: Existing and Historical Conditions - Causes and Effects," *National Acid Precipitation Assessment Program Report 24*, Vol. III, Superintendent of Documents, Government Printing Office, Washington, D.C. 20402-9325.
- 6.** Duce, R. A. et al., Committee on Haze in National Parks and Wilderness Areas, (1993), "Protecting Visibility in National Parks and Wilderness Areas," National Research Council, National Academy Press, Washington, D.C., pp. 252-255; 343-351.

## **Light Scattering Measurements Related to Aerosol Mass Concentration**

- 1.** Charlson, R. J., N. C. Ahlquist and H. Horvath, (1968), "On the Generality of Correlation of Atmospheric Aerosol Mass Concentration and Light Scatter," *Atmospheric Environment*, 2, 455-464.

2. Waggoner, A. P., and R. E. Weiss, (1980), "Comparison of Fine Particle Mass Concentration and Light Scattering Extinction in Ambient Aerosol," *Atmospheric Environment*, 14, 623-626.
3. Charlson, R. J., (1969), "Atmospheric Visibility Related to Aerosol Mass Concentration: A Review," *Environmental Science and Technology*, 3, 913-918.
4. Horvath, H. and R. J. Charlson, (1969), "The Direct Optical Measurement of Atmospheric Air Pollution," *American Industrial Hygiene Association Journal*, 30, 500-509.
5. Bennett, B. G., J. G. Kretschmar, G. G. Akland, and H. W. de Koning, (1985), "Urban Air Pollution Worldwide," *Environmental Science and Technology*, 19, 298-304.
6. White, W. H., E. S. Macias, R. C. Nininger, and D. Schorran, (1993), "Size-Resolved Measurements of Light Scattering by Ambient Particles in the Southwestern U.S.," *Atmospheric Environment*, \_\_, \_\_-\_\_. See also paper 5 under "Applications of the Integrating Nephelometer to Studies of Atmospheric Visibility," and references therein.

### **Vertical Column and Horizontal Profile Studies with the Integrating Nephelometer**

1. Ahlquist, N. C. and R. J. Charlson, (1968), "Measurement of the Vertical and Horizontal Profile of Aerosol Concentration in Urban Air with the Integrating Nephelometer," *Environmental Science and Technology*, 2, 363-366.
2. Charlson, R. J., N. C. Ahlquist, H. Selvidge, and P. B. MacCready, Jr., (1969), "Monitoring of Atmospheric Aerosol Parameters with the Integrating Nephelometer," *APCA Journal*, 19, 937-942.
3. Blumenthal, D. L., W. H. White and T. B. Smith, (1978), "Anatomy of a Los Angeles Smog Episode: Pollutant Transport in the Daytime Sea Breeze Regime," *Atmospheric Environment*, 12, 893-907.

4. White, W. H., J. A. Anderson, D. L. Blumenthal, R. B. Husar, N. V. Gillani, J. D. Husar, and W. E. Wilson, Jr., (1976), "Formation and Transport of Secondary Air Pollutants: Ozone and Aerosols in the St. Louis Urban Plume," *Science*, *194*, 187-189.
5. Bodhaine, B. A., N. C. Ahlquist, and R. C. Schnell, (1991), "Three-Wavelength Nephelometer Suitable for Aircraft Measurement of Background Aerosol Scattering Coefficient," *Atmospheric Environment*, *25A*, 2267-2276.

### **Apportioning of Light Scattering Coefficient to Chemical Species and to Sources**

1. Waggoner, A. P., A. J. Vanderpol, R. J. Charlson, S. Larsen, L. Granat, and C. Trägårdh, (1976), "Sulphate-Light Scattering Ratio as an Index of the Role of Sulphur in Tropospheric Optics," *Nature*, *261*, 120-122.
2. Hasan, H. and T. G. Dzubay, (1983), "Apportioning Light Extinction Coefficients to Chemical Species in Atmospheric Aerosol," *Atmospheric Environment*, *17*, 1573-1581.
3. Macias, E. S., J. O. Zwicker and W. H. White, (1981), "Regional Haze Case Studies in the Southwestern U.S.—II. Source Contributions," *Atmospheric Environment*, *15*, 1987-1997.
4. Appel, B. R., Y. Tokiwa, J. Hsu, E. L. Kothny, and E. Hahn, (1985), "Visibility as Related to Atmospheric Aerosol Constituents," *Atmospheric Environment*, *19*, 1525-1534.
5. Dzubay, T. G., R. K. Stevens, C. W. Lewis, D. H. Hern, W. J. Courtney, J. W. Tesch, and M. A. Mason, (1982), "Visibility and Aerosol Composition in Houston, Texas," *Environmental Science and Technology*, *16*, 514-525.
6. Sloane, C. S., (1988), "Forecasting Visibility Impairment: A Test of Regression Estimates," *Atmospheric Environment*, *22*, 2033-2045.
7. White, W. H., (1986), "On the Theoretical and Empirical Basis for Apportioning Extinction by Aerosols: A Critical Review," *Atmospheric Environment*, *20*, 1659-1672.



## Use of the Integrating Nephelometer in Studies of Physical and Chemical Properties of Aerosols

1. Charlson, R.J., A. H. Vanderpol, D. S. Covert, A. P. Waggoner, and N. C. Ahlquist, (1974), "Sulfuric Acid-Ammonium Sulfate Aerosol: Optical Detection in the St. Louis Region," *Science*, *184*, 156-158.
2. Vanderpol, A. H., F. D. Carsey, D. S. Covert, R. J. Charlson and A. P. Waggoner, (1975), "Aerosol Chemical Parameters and Air Mass Character in the St. Louis Region," *Science*, *190*, 570.
3. Weiss, R. E., A. P. Waggoner, R. J. Charlson, and N. C. Ahlquist, (1977), "Sulfate Aerosol: Its Geographical Extent in the Midwestern and Southern United States," *Science*, *195*, 979-981.
4. Charlson, R. J., D. S. Covert, T. V. Larson, and A. P. Waggoner, (1978), "Chemical Properties of Tropospheric Sulfur Aerosols," *Atmospheric Environment*, *12*, 39-53.
5. Sloane, C. S., M. J. Rood, and C. F. Rogers, (1991), "Measurements of Aerosol Particle Size: Improved Precision by Simultaneous Use of Optical Particle Counter and Nephelometer," *Aerosol Science and Technology*, *14*, 289-301.
6. Weiss, R. E., T. V. Larson, and A. P. Waggoner, (1982), "In-Situ Rapid Response Measurement of  $\text{H}_2\text{SO}_4/(\text{NH}_4)_2\text{SO}_4$  Aerosols in Rural Virginia," *Environmental Science and Technology*, *16*, 525-532.
7. Bhardwaja, P. S., J. Heintzenberg, and R. J. Charlson, (1976), "Refractive Index of Atmospheric Aerosols," *Proceedings of the Symposium on Radiation in the Atmosphere*, pp. 91-93. IAMAP, Garmisch.
8. Kapustin, V. N., G. V. Rozenberg, N. C. Ahlquist, D. S. Covert, A. P. Waggoner, and R. J. Charlson, (1980), "Characterization of Nonspherical Atmospheric Aerosol Particles with Electrooptical Nephelometry," *Applied Optics*, 1345-1348.

9. Rood, M. J., D. S. Covert, and T. V. Larson, (1987), "Hygroscopic Properties of Atmospheric Aerosol in Riverside, California," *Tellus*, 39B, 383-397.
10. Rood, M. J., D. S. Covert, and T. V. Larson, (1987), "Temperature and Humidity Controlled Nephelometry: Improvements and Calibration," *Aerosol Science and Technology*, 7, 57-65.
11. Covert, D. S., (1988), "North Pacific Marine Background Aerosol: Average Ammonium to Sulfate Molar Ratio Equals 1," *Journal of Geophysical Research*, 93, 8455-8458.
12. Rood, M. J., M. A. Shaw, T. V. Larson, and D. S. Covert, (1989), "Ubiquitous Nature of Ambient Metastable Aerosol," *Nature*, 337, 537-539.

### **Observation of the Effects of Humidity on Light Scattering Coefficient; Measurements on Fogs and Clouds**

1. Charlson, R. J., D. S. Covert, T. V. Larson, (1984), "Observation of the Effect on Humidity on Light Scattering by Aerosols," in *Hygroscopic Aerosols*, Lothar H. Ruhnke and Adarsh Deepak (Eds.), A. Deepak Pub. Co., 35-44.
2. Ogren, J. A., J. Heintzenberg and R. J. Charlson, (1985), "In-Situ Sampling of Clouds with a Droplet to Aerosol Converter," *Geophysical Research Letters*, 12, 121-124.
3. Ogren, J. A., J. Heintzenberg, A. Zuber, K. J. Noone, and R. J. Charlson, (1989), "Measurements of the Size-Dependence of Solute Concentrations on Cloud Droplets," *Tellus*, 41B, 24-31.

### **Observation of Aerosols in Remote Locations with the Integrating Nephelometer**

1. Porch, W. M., R. J. Charlson, and L. F. Radke, (1970), "Atmospheric Aerosol: Does a Background Level Exist?" *Science*, 70, 315-317.

2. Bodhaine, B. A. and J. C. Bortniak, (1981), "Four Wavelength Nephelometer Measurements at South Pole," *Geophysical Research Letters*, 8, 539-542.
3. Bodhaine, B. A., J. M. Harris and G. A. Herbert, (1981), "Aerosol Light Scattering and Condensation Nuclei Measurements at Barrow, Alaska," *Atmospheric Environment*, 15, 1375-1389.
4. Heintzenberg, J., (1980), "Particle Size Distribution and Optical Properties of Arctic Haze," *Tellus*, 32, 251-260.
5. Clarke, A. D. and R. J. Charlson, (1985), "Radiative Properties of the Background Aerosol: Absorption Component of Extinction," *Science*, 229, 263-265.
6. Bodhaine, B. A., (1989), "Barrow Surface Aerosol: 1976-1986," *Atmospheric Environment*, 23, 2357-2369.
7. Bodhaine, B. A., (1989), "Aerosol Monitoring of Mauna Loa Observatory," *Proceedings of the International Conference on Global and Regional Environmental Atmospheric Chemistry*, 3-10 May 1989, Beijing, China.
8. Bodhaine, B. A., (1983), "Aerosol Measurements at Four Background Sites," *Journal of Geophysical Research*, 88, 10,753-10,768.

## The Use of Integrating Nephelometer Data in Radiation and Climate Studies

1. Russell, P. B., J. M. Livingston, and E. E. Uthe, (1979), "Aerosol-Induced Albedo Change: Measurement and Modelling of an Incident," *Journal of Atmospheric Science*, 36, 1587-1608.
2. Rosen, J. M., B. A. Bodhaine, J. F. Boatman, J. J. DeLuise, M. J. Post, Y. Kim, R. C. Schnell, P. J. Sheridan, and D. M. Garvey, (1992), "Measured and Calculated Optical Property Profiles in the Mixed Layer and Free Troposphere," *Journal of Geophysical Research*, 97, 12,837-12,850.
3. Charlson, R. J., J. Langner, H. Rodhe, C. B. Leovy, and S. G. Warren, (1991), "Perturbation of the Northern Hemispheric

Radiative Balance by Backscattering from Anthropogenic Sulfate Aerosols," *Tellus*, 43AB, 152-163.


4. Charlson, R. J., S. E. Schwartz, J. M. Hales, R. D. Cess, J. A. Coakley, Jr., J. E. Hansen, and D. J. Hofmann, (1992), "Climate Forcing by Anthropogenic Aerosols," *Science*, 255, 423-430.
5. Penner, J. E., R. J. Charlson, J. M. Hales, N. Laulainen, R. Liefer, T. Novakov, J. Ogren, L. F. Radke, S. E. Schwartz, and L. Travis, (1993), "Quantifying and Minimizing Uncertainty of Climate Forcing by Anthropogenic Aerosols," U.S. Department of Energy Report DOE/NBB-0092T.

# APPENDIX C


## Schematics

This appendix contains the following drawing and schematics:

- ❑ Nephelometer sensor connectors
- ❑ 9060583, revision B, 1 sheet
- ❑ 9060584, revision B, 1 sheet
- ❑ 9060585, revision B, 4 sheets
- ❑ 9060586, revision A, 1 sheet
- ❑ 9060587, revision B, 1 sheet
- ❑ 9060601, revision D, 3 sheets







REVISIONS				
REV	ECO #	DESCRIPTION	DATE	INIT.
A	93-3085	MINOR COMP. CHGS. COLORS.	1/6/94	GV RH
B	94-3017	COMP CHGS TO MATCH UP PWM	5/9/94	RH FO

CONN. J12 GOES TO CONN. J12 OF PREAMP PCB

TO PHT PCB CONN. J4

TO PHT PCB CONN. J5

TO PHT PCB CONN. J7

TO PHT PCB CONN. J6

NOTE: C11-C31, C38-C58, AND C66-C86 MUST BE 200V

PCB 2605955 REV B


DATE: 4-27-93  
SCALE: NONE

**TSI** TSI INCORPORATED  
ST. PAUL, MN.

SCHEMATIC DIAGRAM -  
355X/6X NEPHELOMETER  
HVPS PCB

DRAWN BY: R. HOLM  
CHK: F. QUANT

SIZE: D  
DRAWING NO.: 9060583  
SHEET 1 OF 1



REVISIONS				
LTR	ECO #	DESCRIPTION	DATE	INIT
A	93-3081	REASSIGNED COLORS	1/6/94	GV RH
B	94-3018	CUT GAIN ON U1,2,3	5/10/94	RH FD

FROM DIGITAL PCB

J11 4 BLU-  
 J11 5 BLU+  
 J11 9 GND

J12 3 GRN\_HV\_RET  
 J12 2 GRN\_HV\_CON  
 J12 7 BLU\_HV\_RET  
 J12 6 BLU\_HV\_CON  
 J12 5 RED\_HV\_RET  
 J12 4 RED\_HV\_CON

CONN. J12 TO J12 ON HVPS PCB

THE RED, GRN, BLU HV RET. ARE TIED TO GROUND FROM THE DIGITAL PCB

J11 1 GRN-  
 J11 2 GRN+  
 J11 6 GND

TP1  
 TP2  
 J12 8 +15V  
 J12 1 GND  
 J12 3 GNDL  
 J12 9 +5V  
 J12 10 GND

CONN. J12 TO J12 ON HVPS PCB

J11 7 RED-  
 J11 8 RED+  
 J11 19 GND  
 J11 20 GND


TP3  
 TP4  
 TP5  
 TP6

PCB 2605956 REV B


DATE	4-27-93	TSI INCORPORATED ST. PAUL, MN.
SCALE	NONE	
SCHEMATIC DIAGRAM - 355X/6X NEPHELOMETER PRE-AMP PCB		
DRAWN BY	R. HOLM	SIZE <b>9060584</b>
CHK	F. QUANT	
SHEET 1 OF 1		REV B




REVISIONS				
REV	ECO #	DESCRIPTION	DATE	APPROVED
A	93-3083	UPDATE;CHG NAME FR. MAIN	1/10/94	GV FQ
B	96-3020	ADDED TEST PINS CS & SS	1/9/97	TK




### SDI-12



### RS-232



PCB 2605957 REV B	
DATE 7-19-93	SCALE NONE
<b>TSI</b> TSI INCORPORATED ST. PAUL, MN.	
<b>SCHEMATIC DIAGRAM - 355X/6X NEPHELOMETER DIGITAL PC BOARD</b>	
DRAWN BY T.KERRICK	REV D
CHK R.HOLM 1258	DRAWING NO. 9060585
	SHEET 1 OF 4




PCB 2605957 REV B	
DATE 7-19-93	TSI TSI INCORPORATED ST. PAUL, MN.
SCALE NONE	
SCHEMATIC DIAGRAM - 355X/6X NEPHELOMETER DIGITAL PCB	
DRAWN BY T.KERRICK	SIZE D
CHK R.HOLM 1258	DRWING NO. 9060585
	SHEET 2 OF 4
	REV B


A B C D E F G H J K L M N P Q R

9  
8  
7  
6  
5  
4  
3  
2  
1  
0


### MOTOR CONTROLLER



### CHOPPER





### SHUTTER




PDB 2605957 REV B			
DATE	7-19-93	TSI TSI INCORPORATED ST. PAUL, MN.	
SCALE	NONE		
SCHEMATIC DIAGRAM - 355X/6X NEPHELOMETER DIGITAL PCB			
DESIGNED BY	T. KERRICK	SIZE	D
CHECKED BY	R. HOLM 1258	DRAWING NO.	9060585
		SHEET	3 OF 4
		REV	B

A B C D E F G H J K L M N P Q R


# DIGITAL-ANALOG




## PMT BOX




## CHOPPER




## SPARES



## DEBUGGER






## SHUTTER



PCB 2605957 REV B			
DATE	7-19-93	TSI INCORPORATED ST. PAUL, MN.	
SCALE	NONE		
SCHEMATIC DIAGRAM - 355X/6X NEPHELOMETER DIGITAL PCB			
DRWN BY	T. KERRICK	SIZE	D
CHK	R. HOLM 1258	DRAWING NO.	9060585
		SHEET	4 OF 4
		REV	B



REVISIONS				
LTR	ECO #	DESCRIPTION	DATE	INIT
A	93-3075	CORRECTED NET LABELING	1/6/94	GV RH
B	94-3016	UPDATED TO MATCH PCB REV	5/10/94	RH FO




PCB 2605999 REV B

DATE	1-27-93	 TSI INCORPORATED ST. PAUL, MN.
DATE	NONE	
SCHEMATIC DIAGRAM - PMT SOCKET PCB		
DESIGN BY	R. HOLM	D 9060587 SHEET 1 OF 1
CHKD BY	F. QUANT	
REV	B	


REVISIONS

REV	ECO #	DESCRIPTION	DATE	INIT.
D	97-3124	INCREASED TRACE WIDTHS ON LAMP DRIVE CIRCUIT CHANGED CR2 & R34 VALUES SEE ECO	11/14/97	TK
D	3-22641480-RH	U12(MP7610BS) D/A IC OBSOLETE REMOVE RELATED CIRCUIT	3/14/2001	PC


+ -15 VOLT SUPPLY




+5 VOLT SUPPLY



LAMP DRIVE



8 CHAN ANALOG OUTPUT




BLOWER


FAN

VALVE

HEATER



BNC\_I/O




NO LONGER INSTALLED AS OF 3/2001

BOM 1035507	PCB 2605966D	ASY 9081779
DATE 8-5-93	SCALE NONE	TSI TSI INCORPORATED ST. PAUL, MN.
SCHEMATIC DIAGRAM - 355X/6X NEPHELOMETER PWRSPLY PCB		
DRAWN BY T.KERRICK	CHEK R.HOLM 1258	REV D
SHEET 1 OF 3	DRAWING NO. 9060601	REV D



REVISIONS

REV	ECO #	DESCRIPTION	DATE	INIT.
D	97-3124	INCREASED TRACE WIDTHS ON LAMP DRIVE CIRCUIT CHANGED CR2 & R34 VALUES SEE ECO	11/14/97	TK
D	3-22641480-RH	U12(MP7610BS) D/A IC OBSOLETE REMOVE RELATED CIRCUIT	3/14/2001	PC





PRESSURE

LAMP TEMP



TEMP\_OUT

VALVE POSITION



TEMP\_IN

BOM 1035507 PCB 2605966D ASY 9081779

DATE 8-5-93  
SCALE NONE  
TSI TSI INCORPORATED ST. PAUL, MN.

SCHEMATIC DIAGRAM -  
355X/6X NEPHELOMETER  
PWRSPLY PCB


DRAWN BY T.KERRICK  
CHK R.HOLM 1258  
SHEET 2 OF 3  
DRAWING NO. 9060601  
REV D




REVISIONS

REV	ECO #	DESCRIPTION	DATE	INIT.
D	97-3124	INCREASED TRACE WIDTHS ON LAMP DRIVE CIRCUIT CHANGED CR2 & R34 VALUES SEE ECO	11/14/97	TK
D	3-22641480-RH	U12(MP7610BS) D/A IC OBSOLETE REMOVE RELATED CIRCUIT	3/14/2001	PC

PWR SPLY-MAIN




BACK PANEL




NO LONGER USED  
AS OF 3/2001


LAMP




POWER




VALVE MOTOR




LAMP\_TEMP



TEMP IN/OUT/RH



HEAT/BLOWER/FAN



BOM 1035507 PCB 2605966D ASY 9081779

DATE	8-5-93	 TSI INCORPORATED ST. PAUL, MN.
SCALE	NONE	

SCHEMATIC DIAGRAM -  
355X/6X NEPHELOMETER  
PWRSPLY PCB

DRAWN BY	T.KERRICK	STEP	D	DRAWING NO.	9060601	REV	D
CHEK	R.HOLM 1258	SHEET	3	OF	3		



## APPENDIX D

# Windows Program Files

The program files required for the Windows version of the Nephelometer software are described below.

**Note:** *In addition to these files, your directory will also contain data files and possibly other files.*

The .vbx and .dll files are support files required for program operation. The .vbx files are:

cmdialog.vbx	csopt.vbx
cschk.vbx	csspin.vbx
cscmd.vbx	cstext.vbx
cscombo.vbx	graph.vbx
csdialog.vbx	mscomm.vbx
csform.vbx	ssdata2.vbx
csgroup.vbx	threed.vbx

The .dll files are:

gswdll.dll  
qpro200.dll  
vbrun300.dll

The .exe files are required for program execution. The .exe files are:

gsw.exe	Required for graphical display.
nephmain.exe	The main execution file for the program.

The .dat files are data files. The only two data files provided with the program are:

nephcnfg.dat	The current Nephelometer configuration information.
rayscat.dat	
sample.dat	The sample log data file.



# Reader's Comments

Please help us improve our manuals by completing and returning this questionnaire to the address listed in the "About This Manual" section. Feel free to attach a separate sheet of comments.

**Manual Title** \_\_\_\_\_ **P/N** \_\_\_\_\_ **Rev.** \_\_\_\_\_

1. Was the manual easy to understand and use?

Yes       No

Please identify any problem area(s) \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

2. Was there any incorrect or missing information? (please explain) \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

3. Please rate the manual according to the following features:

	<b>Good</b>	<b>Adequate</b>	<b>Poor</b>
Readability	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Accuracy	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Completeness (is everything there?)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Organization (finding what you need)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Quality and number of illustrations	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Quality and number of examples	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Comments: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

4. Which part(s) of this manual did you find most helpful? \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

5. Rate your level of experience with the product:

Beginning       Intermediate       Expert

6. Please provide us with the following information:

Name \_\_\_\_\_ Address \_\_\_\_\_

Title \_\_\_\_\_

Company \_\_\_\_\_